



## PuraLev® Life Science Pump Series




### PuraLev® 600SU (Single-Use)

3.1 bar	(45 psi)
75 liters/min	(20 gallons/min)

Low Shear Design - High Cell Viability



**Figure 1:** Schematic of the main elements of the MagLev centrifugal pump.



**Figure 2:** Single-use pump head with motor

## INTRODUCTION

Levitronix® has developed a revolutionary pump that has no bearings to wear out or seals to break. Based on the principles of magnetic levitation, the pump's impeller is suspended, contact-free, inside a sealed casing and is driven by the magnetic field of the motor (Figure 1). The impeller and casing are both fabricated from biocompatible (FDA, USP-VI, BSE/TSE and Animal free) gamma sterilizable polypropylene (PP) and together they make up the disposable pump head. A simple and intuitive exchange of the single use pump head is achieved with a bayonet socket type mounting procedure (see Figure 4). Flow rate or pressure are precisely controlled by electronically regulating the rotor speed, which eliminates any pulsation. With the lack of mechanical bearings plus the self-contained pump head design, the risk of contamination is drastically reduced. The absence of narrow gaps between the impeller and pump casing, plus the low-shear pump design allows the gentle pumping of sensitive liquids. The pump casing is fabricated with Triclamp fittings and can be easily inserted and removed with an intuitive bayonet socket.

## SYSTEM BENEFITS

- Low shear-forces
- Reduced risk of contamination due to the self-contained design with magnetic bearings
- No particle generation
- No over-pressure situations (compared to roller pumps)
- No narrow gaps between the impeller and pump casing where bacteria could be entrapped
- Pump head is gamma sterilizable
- Biocompatibility of wet materials: FDA, USP-VI, Animal/BSE/TSE free
- Bayonet socket design for easy and intuitive exchange of disposable pump head (see Figure 4)
- Small size
- Dry running capability
- Proven technology in the medical (disposable blood pumps) and semiconductor (high-purity pumps) industries
- High flow capability with compact design
- Pulsation free

## APPLICATIONS

- Pumping of shear-sensitive liquids and cells
- Bioprocessing (for example perfusion)
- Recirculation and transfer applications in bioreactors
- Filtration

## STAND-ALONE SYSTEM CONFIGURATION

The stand-alone configuration of the *PuraLev® 600SU* pump system consists of a controller with an integrated user panel allowing the operator to set the speed manually (see *Figure 5*). The speed is automatically stored in the internal EEPROM of the controller. As an option, the speed can also be set with an analogue signal (see specification for *Position 3a* in *Table 2*).







## EXTENDED SYSTEM CONFIGURATION


The extended version of the *PuraLev® 600SU* pump system (*Figure 6*) consists of a controller with an extended PLC interface. The PLC interface allows the speed to be set via an external signal, facilitating precise closed-loop flow or pressure control when either a flow or pressure sensor is integrated into the system (see specification of *Position 3b* in *Table 2*). A computer can be connected via a USB interface to allow communication with *Levitronix® Service Software*. Hence parameterization, firmware updates and failure analysis are possible.

## ATEX / IECEx SYSTEM CONFIGURATION


An *ATEX / IECEx* certified motor together with the pump head allows installation of motor and pump head within an *ATEX Zone 2* area (see *Figure 7*). The *ATEX / IECEx* motor (*Pos. 2b* in *Table 2*) comes with special connectors and relevant extension cables (*Pos. 5a* and *5b* in *Table 3*). An *Ex* conform solution is needed for the motor cables to leave the *ATEX* area. One option is an *ATEX* certified cable sealing system as listed in *Table 4* (see *Pos. 9*) and shown in *Figure 11*.

- *ATEX / IECEx* certified for *Category 3G* and *3D* (*Zone 2* for Gas and *Zone 22* for Dust).
- *ATEX / IECEx* marking of motor with pump head:



 II 3G Ex ec h mc IIC T4 Gc  


 II 3D Ex h tc IIIC T105°C Dc



**Figure 3:** Pressure/flow curves (DCP-600.2 pump head)



**Figure 4:** Intuitive 3-step pump head mounting procedure with bayonet type socket (PHS-600.1)

# SYSTEM CONFIGURATIONS




Figure 5: System configuration for standalone operation (Speed setting with integrated user panel)




Figure 6: Extended operation (flow or pressure control) with extended controller






Figure 7: System Configuration for ATEX / IECEx applications

# DIMENSIONS OF MAIN COMPONENTS




Controller  
LPC-600.1-01

4a




Controller  
LPC-600.2-01



4b




**Figure 8: Dimensions of controllers**  
Note 1: Non-tolerated dimensions are for reference only.



**45° Mounting Position**



**Vertical Mounting Position**



**Figure 9: Dimensions of motor with single-use pump head**  
Note 1: Non-tolerated dimensions are for reference only.



# ORDER INFORMATION

System Name	Article #	Pump Head Socket	Motor	Controller	Note
PLD-600SU.1	100-90654	PHS-600.1	LPM-600.5	LPC-600.1-01	Adaptor/Extension (0.5 - 10m) cables according to Table 3 have to be ordered as separate article with specified length. ATEX Cable Sealing System can be ordered according to Table 4. Certifications: CE, IECEx CB scheme, ETL (NRTL), ATEX and IECEx.
PLD-600SU.2	100-90655		LPM-600.5	LPC-600.2-01	
PLD-600SU.4	100-90657		LPM-600.4	LPC-600.1-01	
PLD-600SU.5	100-90658		LPM-600.4	LPC-600.2-01	

**Table 1:** Standard system configurations with motor, pump head socket and controller

Pos.	Component	Article Name	Article #	Characteristics	Value / Feature
1a	Single-Use (SU) Pump Heads	DCP-600.2	100-90784	Impeller / Pump Housing In-/Outlet Fittings	Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) Triclamp 1" (ASME BPE 2009)
1b		DCP-600.2-G25 (Gamma irradiated with Dosage ≥ 25 kGy)	100-91079	Max. Flow Max. Diff.-Pressure Max. Viscosity	75 liters/min / 20 gallons/min 3.1 bar / 45 psi 50 cP
1c	SU Pump Head with Sterile Fittings	DCP-600.2-SF1-G25	100-91294	Wet Pump Volume/Surface Max. Liquid Temp. Sterilization Methods	113 ml / 370 cm <sup>2</sup> 60°C / 140°F Gamma radiation up to 40 kGy.
2a	Motor (ATEX / IECEx)	LPM-600.5	100-10039	Pump Type (A) / Adaptor (C) Sterile Fittings (B) Fitting Compatibility Applied Gamma Dosage	DCP-600.2 / Triclamp reducer in Polypropylene AseptiQuik® X from CPC® with part # AOX33024 <sup>1</sup> Various types including autoclavable versions available at CPC® ≥ 25 kGy
2b	Motor (ATEX / IECEx)	LPM-600.4	100-10038	Housing Cable / Connectors	Epoxy (anti-corrosive) coated Aluminum, waterproofed (IP67) 2x 3m cables with PVC jacket / 2x circular (M23, IP-67)
3	Pump Head Socket	PHS-600.1	100-90696	ATEX / IECEx Marking	CE Ex ec h mc IIC T4 Gc / CE Ex ec h tc IIC T105°C Dc
4a	Standalone Controller (User Panel)	LPC-600.1-01	100-30039	Housing Cable / Connectors	ETFE (chemical resistant) coated Aluminum, waterproofed (IP67) 2x 3m cables with FEP jacket / 2x circular (M23, IP-67)
4b	Extended Controller (PLC and USB)	LPC-600.2-01	100-30040	Mounting Type Material Assembly Screws	Bayonet type with locking pin Anodized Aluminum 4 pcs M6x16 mm (Stainless Steel)
				Voltage / Power Housing Rating Interfaces	48V DC / 600 W IP20 Panel to set speed (automatic storage on internal EEPROM)
				Standard Firmware	D9.25
				Interfaces	PLC with - 1x analog input ("Speed") - 1x digital input ("Enable") - 1x digital output ("Status")
				Standard Firmware	D9.48

**Table 2:** Specification of standard components (Note 1: CPC® and AseptiQuik® are registered marks of the Colder Product Company)

Pos.	Component	Article Name	Article #	Characteristics	Value / Feature
5a	Extension Adaptor Cable for Sensors	MCAS-600.2-05 (0.5m)	190-10226	Jacket Material Connectors	PVC Circular Wallmountable, Metallic (IP-67) to D-SUB
		MCAS-600.2-30 (3m)	190-10238		
		MCAS-600.2-50 (5m)	190-10127		
		MCAS-600.2-70 (7m)	190-10105		
		MCAS-600.2-100 (10m)	190-10239		
5b	Extension Adaptor Cable for Power	MCAP-600.2-05 (0.5m)	190-10227	Jacket Material Connectors	PVC Circular Wallmountable, Metallic (IP-67) to COMBICON
		MCAP-600.2-30 (3m)	190-10240		
		MCAP-600.2-50 (5m)	190-10126		
		MCAP-600.2-70 (7m)	190-10106		
		MCAP-600.2-100 (10m)	190-10241		

**Table 3:** Specification of adaptor/extension cables

Pos.	Component	Article Name	Article #	Characteristics	Value / Feature
6a 6b	Air Cooling Module	ACM-600.2 ACM-600.3 (for ATEX)	190-10140 190-10410	Housing / Connection Port Air Pressure	a) PP (+ 40% Talkum) b) PP EL-S (conductive black PP) / NPT 1/4" -1 - 3 bar (14 - 43 psi)
7a	Fan Cooling Module	FCM-600.1	190-10401	Housing / Cable Spec. Supply Spec. / IP Rating	PP (+ 20% Talkum) white / PP jacket, 3m, circular sealed M12 connector (PP). 24 VDC, 3.4 W / IP-65 (fan is IP68 rated).
7b	Fan Cool. Module Cable	FCC-1.1-50 (5 m) FCC-1.1-100 (10 m)	190-10407 190-10408	Specification	PP cable jacket with circular M12 connector (PP) to open wires
8 (A-F)	ATEX Cable Sealing System	ACS-A.1 (Roxtec)	100-90292	Sleeve (a) and Gasket (b) Frame (c) 2x Cable Module (d)	Stainless Steel and EPDM Roxylon (EPDM rubber) Roxylon (EPDM rubber) Note: Lubricant (e) and measurement plates (f) are included.
9	AC/DC Power Supply	TSP 600-148-M (M = Modified Levitronix design from Traco)	100-40013 (Traco ID Number: T1068-01D)	Voltage / Power Output Voltage Input Certification or Standards	48 VDC / 600 W 85 - 265 VAC (automatic detection) CB, UL, CSA, Semi F47

**Table 4:** Specification of accessories

# COMPONENTS



Figure 10: Pump system with standard components




Figure 11: Accessories

Levitronix® is the world-wide leader in magnetically levitated bearingless motor technology. Levitronix® was the first company to introduce bearingless motor technology to the Semiconductor, Medical and Lifescience markets. The company is ISO 9001 certified. Production and quality control facilities are located in Switzerland. In addition, Levitronix® is committed to bring other highly innovative products like the LEVIFLOW® flowmeter series to the market.



#### Headquarter and European Contact

Levitronix GmbH  
Bändliweg 30  
CH-8048 Zurich  
Switzerland

Phone: +41 44 974 4000  
E-Mail: [salesEurope@levitronix.com](mailto:salesEurope@levitronix.com)

#### US Contact

Levitronix Technologies Inc.  
10 Speen Street, Suite 102  
Framingham, Massachusetts 01701  
USA

Phone: +1 508 861 3800  
E-Mail: [salesUS@levitronix.com](mailto:salesUS@levitronix.com)

#### Japan Contact

Levitronix Japan K.K.  
Wing Eight 5floor, 4-16-4  
Asakusabashi, Taito-ku  
Tokyo, 111-0053 Japan

Phone: +81 3 5823 4193  
E-Mail: [salesJapan@levitronix.com](mailto:salesJapan@levitronix.com)

#### Taiwan Contact

Levitronix Taiwan  
5F, No. 251, Dong Sec. 1, Guangming  
6th Rd., Chu Pei City,  
Hsin-Chu 302, Taiwan, R.O.C.

Phone: +886 3 657 6209  
E-Mail: [salesAsia@levitronix.com](mailto:salesAsia@levitronix.com)

*This document and its content are the property of Levitronix® and shall not be reproduced, distributed, disclosed or used for manufacturing or sale of Levitronix® products without the expressed written consent of Levitronix®.*