

PuraLev® Life Science Integrated Flow Controller Series



PuraLev® iF30SU (Single-Use)

Pump Pressure / Flow: 1.0 bar / 7.7 l/min

Single-Use Flow Sensor LFS-03SU: 0.8 l/min

Single-Use Flow Sensor LFS-06SU: 8 l/min

Low Shear Design - High Cell Viability




Figure 1: Schematic of the main elements of the MagLev centrifugal pump




Figure 2: Integrated MagLev flow control system with ultrasonic flow sensors

INTRODUCTION

With the *PuraLev® iF30SU* flow control system *Levitronix®* combines its unique magnetic levitation pump technology with its ultrasonic flow measurement technology. The result is a highly integrated precise flow controller with an integrated pressure source. The centrifugal pump, as a pressure source, has no bearings to wear out or seals to break down and fail. The pump impeller is suspended, contact-free, inside a casing and is driven by the magnetic field of the motor (*Figure 1*).

The pump controller, motor and flow converter are integrated into the driver housing. This reduces cabling and setup effort significantly. Fluid flow rate is precisely controlled by electronically regulating the impeller speed without pulsation. The pump head can be easily inserted and removed with an intuitive bayonet socket.

With the lack of mechanical bearings plus the self-contained pump head design, the risk of contamination is drastically reduced. The absence of narrow gaps between the impeller and pump casing, plus the low-shear pump design allows the gentle pumping of sensitive liquids.

SYSTEM BENEFITS

- High precision, dynamics and turndown ratio.
- No dependency on external pressure source.
- Extremely low particle generation due to the absence of mechanically contacting parts.
- Reduced risk of contamination due to the self-contained design with magnetic bearings and ultrasonic technology.
- Very gentle to sensitive fluids due to low-shear design.
- No narrow gaps and fissures where particles or microorganisms could be entrapped.
- Smooth, continuous flow without pressure pulsation.
- No over-pressure situations (compared to roller pumps).
- Biocompatible and gamma sterilizable single-use parts.
- Dry running capability.
- High flow capability with compact design.
- Very low integration costs as no external controller is needed for flow control.
- Proven pump and ultrasonic flow measurement technology.

APPLICATIONS

- Pumping of shear-sensitive liquids and cells.
- Bioprocessing (for example perfusion).
- Recirculation and transfer applications in bioreactors.
- Filtration.

FLOW CONTROL CONCEPT

Figure 4 illustrates the flow control concept. Flow control, pressure generation and flow measurement is done with one unit. This allows realization of sophisticated flow control algorithms and optimizations to various situations.

There is a linear relationship between flow and speed (Figure 3). The speed is precisely controlled with a high resolution over a wide pump speed range. This allows a flow control with high resolution and high turndown ratio compared to non-linear flow control with valve type flow control concepts. Additionally, the highly dynamic speed controller allows fast flow step responses.

As the speed is monitored and the pressure cannot increase uncontrolled at a given speed, there is no need to protect the hydraulic circuit against over-pressure situations as for example for roller pumps in tube clogging situations.

The versatility of *Levitronix*® flow control systems goes far beyond the capabilities of simple flow controllers. In addition to the flow control function, the *Levitronix*® control firmware comes with several condition monitoring features to monitor the integrity of the fluid circuit. *Levitronix*® flow control systems can generate alarms for preventive filter exchange, no-flow conditions or line clogging. Dynamic Condition Trending (DCT) enables failure prediction and scheduling of preventive maintenance

SYSTEM CONFIGURATION – “STAND-ALONE”

Figure 6 and Figure 10 illustrate a “Plug and Play” stand-alone system with integrated user panel and buttons to set the flow manually. The driver also contains a PLC interface for remote flow control by analog and digital signals. Various accessories are available like a desktop power supply with relevant power cable and signal cables to connect to the PLC.

SYSTEM CONFIGURATION – “EASYCONNECT”

The “EasyConnect” models (see Figure 7 and Figure 12) with according cable accessories are designed to realize various interface configurations with minimal setup effort.

Two Fieldbus connectors (IN and OUT) allow to chain multiple flow controllers. Therefore, blending configurations as shown in Figure 8 can be realized.

The PLC interface allows remote control by analog/digital signals but also connections of other signal and sensors.

The Fieldbus interface allows remote control over a PC, a User Panel or other devices with Modbus protocol.

SYSTEM CONFIGURATION – “OEM”

The “OEM” models are designed for a compact integration with one integrated hybrid connector containing all available interface signals (see Figure 9 and Figure 14). Basically, all configurations of the “EasyConnect” models are possible allowing the users with integration capabilities to adapt the cable to their needs.




Figure 3: Pressure/flow curves for aqueous liquids (similar to water)




Figure 4: Simplified block schematics of flow control with PuraLev® IF30SU




Figure 5: Flow control system models

SYSTEM CONFIGURATIONS




Figure 6: Standard "Stand-Alone" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)




Figure 7: Standard "EasyConnect" system configuration with main accessories
(See section "Order Information" for details to numbered components and other options)

SYSTEM CONFIGURATIONS




Figure 8: Flowcontrol array configuration with "EasyConnect" models for blending applications
(See section "Order Information" for details to numbered components and other options)






Figure 9: Standard "OEM" system configuration
(See section "Order Information" for details to numbered components and other options)

MODEL DESCRIPTION – STAND-ALONE




Interface	PIN Name	Description	Standard Designation	Hardware Specification
Power Supply	P+	+ 24 VDC	Supply	Voltage: 24 VDC Power: 35 W
	P-	Power Input Ground / Earth		
	NC	Not connected.	--	--
PLC 6	Ain	Analog Input (Current Input)	Reference (Set) Flow	Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation)
	Ain_GND	Analog In. GND	--	Reference for Ain
	Dout	Digital Output 1	Status	Open drain, max. 24V, 100mA Reference ground is GND
	GND	Analog Ground	--	Reference for Dout
	Din1	Digital Input 1	Enable (Reset)	Galvanic separation with optocoupler 2.2 k Ω input resistance, 5-24V for active input
	Din_COM	Com. Digi. Input	--	Reference for digital input.
Flow Sensor	6 Pins	Flow Sensor	--	Compatible to LEVIFLOW [®] flow sensors of the LFS-SU-F1 series.
Display and Buttons	--	Display	Flow and Status Display	--
	--	Up/Down	Setting Flow	--
	--	On/Off	Enable/Disable	--

Figure 10: Interface specifications of standard “Stand-Alone” model



Pump Head Inlet and Outlet Fitting:
Barb 3/8" for tubing with typical ID = 1/4"



Note: Pump Head Socket can be mounted and rotated in steps of 45° to change pump outlet orientation.






Figure 11: Basic dimensions and description of standard “Stand-Alone” model

MODEL DESCRIPTION – EASYCONNECT




Connector	PIN Name	Description	Standard Designation	Hardware Specification
Power Supply	P+	+ 24 VDC	Supply	Voltage: 24 VDC
	P-	Ground / Earth		Power: 35 W
	NC	Not connected.		--
PLC 12	Dout1	Digital Output 1	Status Pump	Open drain, max. 24V, 100mA
	Dout2	Digital Output 2	Status Flow Sen	Reference ground is GND
	Din1	Digital Input 1	Enable (Reset)	Galvanic separation with optocoupler
	Din2	Digital Input 2	Zero Adjust	2.2 kΩ input resistance, 5-24V for active input
	Din_COM	Com. Digi. Input	--	Reference for digital input.
	Ain1	Analog Input 1 (Current Input)	Reference Value (Set Flow/Speed)	Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation)
	Ain2	Analog Input 2 (Voltage Input)	Free Configurable	Analog voltage input: 0 – 10V (7.9 kOhm, no galvanic isolation)
	Ain_GND	Analog In. GND	--	Reference for Ain1 and Ain2
	Aout1	Analog Output (Voltage Output)	Actual Flow	0 – 10V (no galvanic isolation) GND is reference
	GND	Analog Ground	--	Reference for Aout1, Dout1, Dout2 and Pout
	Pout	Output +24VDC	Supply Output	For supply of external devices (e.g. sensors). (Current 200mA together with Pout o Fieldbus OUT)
	NC	Not connected.	--	--
	GND	Ground	--	Reference for Pout.
	Pout	Output +24VDC	Supply Output	For supply of external devices (user panels) (Current 200mA together with Pout of PLC 12)
Fieldbus OUT	RS485+	RS485 +	Fieldbus	Modbus RTU protocol
	RS485-	RS485 -		
	Internal	Internal Bus	Do not connect	Internal bus needed to connect pumps for serial pumping.
	Internal	Internal Bus	Do not connect	
Fieldbus IN	GND	Ground	--	Connected to PLC12 GND and reference for supply
	NC	Not connected.	--	--
	RS485+	RS485 +	Fieldbus	Modbus RTU protocol
	RS485-	RS485 -		
	Internal	Internal Bus	Do not connect	Internal bus needed to connect pumps for serial pumping.
	Internal	Internal Bus	Do not connect	
Flow Sensor	6 Pins	Flow Sensor	--	Compatible to LEVIFLOW® flow sensors of the LFS-SU-F1 series.

Figure 12: Interface specifications of standard "EasyConnect" model



Pump Head Inlet and Outlet Fitting:
Barb 3/8" for tubing with typical ID = 1/4"



Note: Pump Head Socket can be mounted and rotated in steps of 45° to change pump outlet orientation.





Figure 13: Basic dimensions and description of standard "EasyConnect" model





Connector	Designation	Description	Standard Designation	Hardware Specification
Hybrid Driver	P+	+ 24 VDC	Supply	Voltage: 24 VDC
	P-	Power Input Ground / Earth		P- to be connected to earth
	Ain1	Analog Input 1 (Current Input)	Reference Value (Set Flow/Speed)	Analog current input: 4 – 20 mA (450 Ohm shunt input, no galvanic isolation)
	Ain2	Analog Input 2 (Voltage Input)	Free Configurable	Analog voltage input: 0 – 10V (7.9 kOhm, no galvanic isolation)
	Ain_GND	Analog Input Ground	--	Reference for Ain1 and Ain2
	Din1	Digital Input 1	Enable (Reset)	Galvanic separation with optocoupler
	Din2	Digital Input 2	Zero Adjust	2.2 kΩ input resistance, 5-24V for active input
	Din_COM	Common Digital Input	--	--
	Aout1	Analog Output (Voltage Output)	Actual Speed	0 – 10V (no galvanic isolation) GND is reference
	Dout1	Digital Output 1	Status Pump	Open drain, max. 24V, 100mA
	Dout2	Digital Output 2	Status Flow	Reference ground is GND
	GND	Analog Ground	--	Reference for Aout1, Dout1 and Dout2
	RS485+	RS485 +	Fieldbus	Modbus RTU protocol
	RS485-	RS485 -		
	Shield	Shielding	Shielding	To be connected to earth (see wire No. 2, P-)
Flow Sensor	6 Pins	Flow Sensor	--	Compatible to LEVIFLOW® flow sensors of the LFS-SU-F1 series.

Figure 14: Interface specifications of standard "OEM" model

Note 1: Power supply wire cross-section is 1.5 mm² and for signal wires 0.14 mm²

Note 2: For more detailed description of interfaces consult user manual



Note: Pump head socket can be mounted and rotated in steps of 45° to change pump outlet orientation.

Figure 15: Basic dimensions and description of standard "OEM" model

FLOW SENSOR SPECIFICATIONS




Figure 16: Dimensions for LFS-03SU-Z-F1 and LFS-06SU-Z-F1 flow sensors

Flow Controller Type	PuraLev® IF30SU with LFS-03SU-F1	PuraLev® IF30SU with LFS-03SU-F1-SC1	PuraLev® IF30SU with LFS-06SU-F1	PuraLev® IF30SU with LFS-06SU-F1-SC1
Characteristics				
Flow Range [lpm]	0 – 0.8	0 – 0.8	0 – 8 (7.7 for pump)	0 – 8 (7.7 for pump)
Accuracy of Reading (at 20°C fluid temperature) Note: Repeatability < Accuracy/2	> 35 ml/min: ±1% < 35 ml/min: Figure 17	> 6 ml/min: ±1% < 6 ml/min: ±0.06 ml/min	> 1.7 l/min: ±1% < 1.7 l/min: ±17 ml/min	> 0.075 l/min: ±1% < 0.075 l/min: ±0.75 ml/min
Response Time: Step from 10 – 90% of full scale.	< 1s ¹	< 1s ¹	< 1s ¹	< 1s ¹
Fluid Temperature / Ambient Temperature	Normal range: 10 – 60 °C (50 – 140 °F)		/ 0 – 40 °C (32 – 104 °F)	

Table 1: Specifications of flow controller systems PuraLev® IF30SU with LFS-SU-F1 single-use flow sensors

1: Values for to the specific hydraulic circuit optimized flow control parameters. Standard settings are tuned for general flow control stability and may be higher.




Figure 17: Accuracy and repeatability for single-use sensor LFS-03SU-F1

ORDER INFORMATION

System Name	Article #	Pump Head Socket	Flow Control Driver	Note
PLD-IF30SU.1	100-91223	PHS-i30.1	IFD-30.3-02	OEM – Driver with one hybrid connector, pump head socket
PLD-IF30SU.2	100-91224	PHS-i30.1	IFD-30.4-02 (MBP-i30.1 included)	EasyConnect - Driver with interface connectors, pump head socket.
PLD-IF30SU.3	100-91225	PHS-i30.1	IFD-30.5-02 (MBP-i30.1 included)	Stand-Alone - Driver with integrated user panel, pump head socket.

Table 2: Standard system configurations

Pos.	Component	Article Name	Article #	Characteristics	Value / Feature
1a 1b	Disposable Pump Head	DCP-30.2 (Barb) ² DCP-30.1 (Triclamp) ²	100-90968 100-90959	Material Impeller and Pump Housing Housing Sealing In-/Outlet Fittings Max. Flow Max. Diff.-Pressure Max. Viscosity Max. Liquid Temp. Wet Pump Volume/Surface Sterilization Methods	Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) Infrared welding Barb 3/8" or Triclamp 3/8" for tubing with typical ID = 1/4" 7.7 liters/min / 2.0 gallons/min 1.0 bar / 14.5 psi 40 cP 60 °C / 140 °F 7.7 ml / 55.9 cm ³ Gamma radiation up to 40kGy
2a	Integrated Flow Control Driver ("OEM Model")	IFD-30.3-02	100-10130	Voltage, Power Housing Interfaces Standard Firmware	24 VDC ±10%, 35 W Epoxy coated Aluminum, PP for bottom lid, IP65 ¹ PLC, RS485 with Modbus RTU protocol (see Figure 14 for details) and flow sensor K2.48
2b	Integrated Flow Control ("EasyConnect" Model)	IFD-30.4-02 (MBP-i30.1 included)	100-10131	Housing Interfaces Standard Firmware	Epoxy coated Aluminum, PP for bottom lid, IP65 2x Fieldbus RS485 with Modbus RTU protocol, PLC, power supply and flow sensor K2.48
2c	Integrated Flow Control ("Stand-Alone" Model)	IFD-30.5-02 (MBP-i30.1 included)	100-10132	Housing Interfaces Standard Firmware	Epoxy coated Aluminum, PP for bottom lid, IP65 User panel with 3 user buttons, PLC, power supply and flow sensor K2.48
3	Pump Head Socket	PHS-i30.1	100-90947	Mounting Type Material Assembly Screws	Bayonet type with locking pin Anodized Aluminum 4 pcs M3 x 6 mm (Stainless Steel, INOX A4)

Table 3: Specification of standard main components

Note 1: Designed and tested for IP67. Note 2: Gamma irradiated versions and versions with sterile fittings available.

Pos.	Component	Article Name	Article #	1% Accuracy Range	Fitting	Wet Material	Note
4a 4b	LEVIFLOW® Single-Use	LFS-03SU-Z-F1 LFS-03SU-Z-F1-SC1 ¹	100-30423 100-30470	35 – 800 ml/min 6 – 800 ml/min	Triclamp 3/8"	Polypropylene (FDA, USP Class VI, BSE/TSE/Animal free) Gamma stable for up to 40 kGy.	See Levitronix® technical brochure of LFS-SU single-use sensor series for more detailed specifications and for other configurations.
5a	Flow Sensors	LFS-06SU-Z-F1	100-30425	1.7 – 8 l/min			
5b		LFS-06SU-Z-F1-SC1	100-30445	0.075 – 8 l/min			

Table 4: Specification of LEVIFLOW® single-use high-precision (1% accuracy of reading) flow sensors compatible with IFD-30 drivers

Note 1: Extended calibration for wider 1% accuracy range. Note 2: All flow sensors available with gamma irradiation (see LEVIFLOW® product literature for more details).

Pos.	Component	Article Name	Article #	Characteristics	Value / Feature
6	Mounting Base Plate	MBP-i30.1	190-10313	Material / Mounting Screws	PP+GF / 2 pieces, stainless steel FEP coated, M3 x 10 mm
7a	AC/DC Power Supply	TPC 055-124 (Traco)	100-40014	Voltage Output / Input Basic Dimensions Certification or Standards	24 VDC with 55 W / 85 – 264 VAC, 47-63 Hz 45 x 90 x 96.5 mm (mountable on DIN rail 35 mm) UL, CSA, CB, Semi F47
7b	Desktop AC/DC Power Supply	AC/DC Power Supply VEC50US24 HR30	100-40015	Voltage Output / Input Basic Dimensions / Cable Spec. Safety Approvals	24VDC, 50W / 90 – 264 VAC, 47-63 Hz 116 x 52 x 31 mm / Cable length 1.2m IEC60950-1, EN60950-1, UL/cUL60950-1
7c	AC Mains Cables (for Desktop power supply 7b)	AMC-1.1 (2m) AMC-1.2 (2.5m) AMC-1.3 (2.5m) AMC-1.4 (2.5m) AMC-1.5 (2.5m)	190-10331 190-10332 190-10333 190-10334 190-10335	Approvals and Country Approvals and Country Approvals and Country Approvals and Country Approvals and Country	UL, cUL, US, Canada CB, Germany, Denmark, Norway, Finland, Belgium, Netherland, Sweden, Austria PSE, Japan Switzerland CE, United Kingdom
8	USB to RS485 Adaptor-TR Isolated	YN-485I-TR	100-30392	Structure/Design Purpose	USB connector (A) with termination resistor and cable (2m) with connector pair (B and C) for external RS485 wire connection. Magnetically isolated. Cable length is 2m. Included is a USB space saver cable (D). Communication over fieldbus of driver with PC
9	IPS Cable Power 3 Wires	ICP-1.1-50 (5 m)	190-10342	Cable Material / Wires Connection In / Connection Out Main Purpose	PVC jacket / 3x 0.5 mm ² (only 2 wires used, 1 is cut) Open wires / Circular Hirose type to driver Connection of power supply to "Stand-Alone" and "EasyConnect" drivers
10a	IPS Cable Signal 6 Wires	ICS-1.1-01 (0.1 m) ICS-1.1-10 (1 m) ICS-1.1-30 (3 m)	190-10343 190-10344 190-10345	Cable Material / Wires Connection In / Connection Out Main Purpose	PVC jacket / 6x 0.08 mm ² and shielding Circular Hirose type / Circular Hirose type Fieldbus connection to "EasyConnect" drivers and flow sensor connection.
10b 10c	IPS Cable Signal 6 Wires	ICS-1.2-50 (5 m) ICS-1.3-50 (5 m)	190-10346 190-10389	Cable Material / Wires Connection In / Connection Out Main Purpose ICS-1.2 / ICS-1.3	PVC jacket / 6x 0.08 mm ² and shielding Connector with screw type plug for open wire connection / Circular Hirose type Fieldbus connection to "EasyConnect" driver / To PLC of "Stand-Alone" driver.
11	IPS Cable Signal 12 Wires	ICS-2.1-50 (5 m)	190-10347	Cable Material / Wires Connection In / Connection Out Main Purpose	PVC jacket / 12x 0.14 mm ² and shielding Connector with screw type plug for open wire connection / Circular Hirose type General connection to PLC of "EasyConnect" drivers.
12	IPS Cable Hybrid 15 Wires	ICH-1.1-30 (3 m) ICH-1.1-50 (5 m)	190-10386 190-10341	Cable Material / Wires Connection In / Connection Out Main Purpose	PVC jacket / 2x 1.5 mm ² for supply wired, 13 x 0.14 mm ² for signal and shielding wire Open wires / Circular hybrid connector for driver connection General connection integrated driver connector of to "OEM" driver models.
13	Fieldbus Termination Connector	FTC-1.1	190-10348	Materials Main Purpose	PPS for connector housing and FPM for sealing. Termination of fieldbus.
14	Mounting Kit	LMK-1.2	100-91478	Material / Structure Main Purpose	Anodized Aluminium / Screw closing concept For mounting of LFS-03SU and LFS-06SU flow sensors.
15	User Panel	LUI-B.1-01	100-30448	Interface / Housing Rating Standard Firmware	RS485 / IP65 A3.00

Table 5: Specification accessories




Figure 18: Flow control system with standard main components



Figure 19: General standard accessories




Figure 20: Standard cables

Levitronix® is the world-wide leader in magnetically levitated bearingless motor technology. *Levitronix®* was the first company to introduce bearingless motor technology to the Semiconductor, Medical and Life Science markets. The company is ISO 9001 certified. Production and quality control facilities are located in Switzerland. In addition, *Levitronix®* is committed to bring other highly innovative products like the *LEVIFLOW®* flowmeter series to the market.



Headquarter and European Contact

Levitronix GmbH
Bandliweg 30
CH-8048 Zurich
Switzerland

Phone: +41 44 974 4000
E-Mail: salesEurope@levitronix.com

US Contact

Levitronix Technologies Inc.
10 Speen Street, Suite 102
Framingham, Massachusetts 01701
USA

Phone: +1 508 861 3800
E-Mail: salesUS@levitronix.com

Japan Contact

Levitronix Japan K.K.
Wing Eight 5floor, 4-16-4
Asakusabashi, Taito-ku
Tokyo, 111-0053 Japan

Phone: +81 3 5823 4193
E-Mail: salesJapan@levitronix.com

Taiwan Contact

Levitronix Taiwan
5F, No. 251, Dong Sec. 1,
Guangming 6th Rd., Chu Pei City,
Hsin-Chu 302, Taiwan, R.O.C.

Phone: +886 3 657 6209
E-Mail: salesAsia@levitronix.com

This document and its content is the property of Levitronix® and shall not be reproduced, distributed, disclosed or used for manufacturing or sale of Levitronix® products without the expressed written consent of Levitronix®.