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- Alias CB, Lopez MCG, Fernandez FGA, Sevilla IMF, Sanchez JLG & Molina-Grima E. (2004). Influence of power supply in the feasibility of phaeodactylum tricornutum cultures. *Biotechnol Bioeng* 87: 23-733.
- Bludszuweit C (1995). Three-dimensional numerical prediction of stress loading Of blood particles in a centrifugal pump.
- Artif Organs, 19: 590—596 Boccafoschi F, Bosetti M, Sandra PM, Leigheb M & Cannas M (2010). Effects of mechanical stress on cell adhesion: A possible mechanism for morphological changes. *Cell Adh Migr*. 4: 19-25.
- Chisti, Y. (1993). Animal cell culture in stirred bioreactors: Observations on scale-up. *Bioprocess Eng*. 9: 191-196.
- Chisti, Y. (2000). Animal-cell damage in sparged bioreactors. *Trends Biotechnol*. 18: 420-432.
- Chua LP, Song G, Lim TM & Zhou T (2006). Numerical analysis of the inner flow field of a biocentrifugal blood pump. *Artif Organs* 30: 467-477.
- Eibl R & Eibl D (2007). Disposable bioreactors for inoculum production and protein expression. In: *Methods in Biotechnolow*, Vol. 24: *Animal Cell Biotechnology*. Methods and Protocols, 2nd Ed., edited by: R. Pörtner, Humana Press, Inc., Totowa, NJ.
- Giersiepen M, Wurzinger LJ, Opitz R & Reul H (1990). Estimation of shear stress-related blood damage in heart valve prostheses — in Vitro comparison of 25 aortic valves. *Int J Artificial Organs*. 13: 300-306.
- Hahn J (2010). Hemolysis study — pump versus bearingless centrifugal pump. Test review, provided by Levitronix.
- Jaouen P, Vandanjon L & Querneneur F (1991). The shear of microalgal cell suspensions (*Tetraselmis suecica*) in tangential flow filtration systems: the role of pumps. *Bioresour Technol*. 68: 149-154.
- Kamaraju H, Wetzel K & Kelly WI (2010). Modeling shear-induced CHO cell damage in a rotary positive displacement pump. *Biocatalysts and Bioreactor Design*. Published online, DOI: 10.1002/btpr.479.
- Kretzmer G (2000). Influence Of stress on adherent cells. In: *Advances in Biochemical Engineering / Biotechnology* (ed. T. Scheper), Vol. 67: *Influence Of stress on cell growth and product formation*, Springer: Berlin, Heidelberg.
- Langer G & Deppe A (2000). Zum Verständnis der hydrodynamischen Beanspruchung von partikeln in turbulenten Rührerströmungen. *Chem Ing Tech*. 72: 31-41.

- Mardikar SH & Niranjan K (2000). Observations on the shear damage to different animal cells in a concentric viscometer. *Biotechnol Bioeng.* 68: 697-704.
- McQueen A & Bailey J (1989). Influence Of serum level, cell line, flow type and Viscosity on flow- induced of suspended animal cells. *Biotechnol Lett.* 11: 531-536.
- Mollet M, Ma N, Zhou Y, Brodkey R, Taticek R & Chalmers J. (2004). Bioprocess equipment: characterization of energy dissipation rate and its potential to damage cells. *Biotechnol Prog.* 20: 1437-1448.
- Ries C, Riesen N & Eibl R (2010). Usage of the Nucleocounter@ NC-300V for advanced cell analyses in insect cell derived processes. Poster presented at BioTech2010 Conference, ZHAW Wädenswil.
- Stintzing A, Pilz RD, Hempel DC & Krull R (2008). Mechanische Beanspruchungen in Mehrphasenreaktoren. *Chem Ing Tech.* 80: 1837-1842.
- Thomas CR (1993). Review: shear effects on cells in bioreactors. *Process Sol-Liq Suspensions.* 158-191.
- Transonic systems (2010). H-XL sterile tubing flow sensors for the HTI 1 O(R) bypass flow meter. Website, URL: <http://www.transonic.com/flowsensors.shtml>, last request: 09.12.2010.
- Wollny S (2010). Experimentelle und numerische Untersuchungen zur Partikelbeanspruchung in gerührten (Bio-) Reaktoren. (Dissertation) Technische Universität Berlin.
- Vandanjon N, Rossignol N, Jaouen P, Robert JM & Quemeneur F (1999). Effects of shear on two microalgae species. Contribution of pumps and valves in flow filtration systems. *Biotechnol Bioeng.* 63: 1-9.
- Yeleswarapu KK, Antaki JF, Kameneva MV & Rajagopal KR (1995). A mathematical model for shear-induced hemolysis. *Artif Organs,* 19: 576—582.
- Zhang J, Gellman B, Koert A, Dasse KA, Gilbert RJ, Griffith BP & Wu ZJ (2006). Computational and experimental evaluation of the fluid dynamics and hemocompatibility of the CentriMag blood pump. *Artif Organs.* 30: 168-177.

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NucleoCounter® NC-3000™

Vitality assay: Analysis of the level of cellular thiols using the NucleoCounter® NC-3000™ system

Product description

The NucleoCounter® NC-3000™ system enables the user to perform automated cell counting and analyses of a broad range of eukaryotic cells.

Application

This protocol for the NucleoCounter® NC-3000™ system enables the user to detect changes in the intracellular level of (reduced) thiols. Such changes may occur in apoptotic cells or cells undergoing other pathological processes. As the intracellular reducing power available to the cell is an indicator of overall health status, this assay provides a very easy and fast way to evaluate cell vitality.

Introduction

This application note describes a method for investigating apoptosis and cell health by determining the level of free thiols such as reduced glutathione. The tripeptide glutathione exists in two forms, a reduced state (GSH) and in an oxidized state; glutathione disulfide (GSSG). In the reduced state the thiol group of cysteine is able to donate a reducing equivalent ($H^+ + e^-$) to unstable molecules such as free radicals. In donating an electron, glutathione itself becomes reactive, but readily reacts with another reactive glutathione to form GSSG. GSH can be regenerated from GSSG by the enzyme glutathione reductase. (See **Figure 1**.)


Figure 1. Glutathione redox cycle. The flavin adenine dinucleotide (FAD)-dependent enzyme, glutathione reductase reduces GSSG to GSH. When the cell encounters free radicals, the antioxidant GSH reduces the free radicals, thereby itself becoming oxidized to GSSG.

GSH is the most abundant low molecular weight thiol in animal cells; thus, its oxidation status largely determines the thiol-disulfide status of the cell by thiol-disulfide interchange reactions. Moreover, GSH is involved in many cellular processes including quenching of free radicals, drug detoxification, cell signaling, and cell proliferation. Alterations in the concentration of intracellular GSH have been demonstrated as a common feature of many diseases including AIDS, neurodegenerative diseases, and cancer.

A decrease in cellular GSH concentration is an early hallmark in the progression of cell death in response to different apoptotic stimuli. Studies have shown a correlation between cellular GSH depletion and the progression of apoptosis. The decrease in GSH level in connection to apoptosis seems to be attributed to two mechanisms; A) direct GSH oxidation promoted by radicals and B) export of GSH through an ATP-dependent plasma membrane transport system which is triggered by the initiation of apoptosis. When GSH is depleted, the cytosol is shifted from a reducing to an oxidizing environment, which may lead to a further depletion of GSH.

This assay provides a very easy method to quantify the amount of free thiols at the single cell level. The stain VitaBright-48 (which is a component of **Solution 5**) immediately reacts with thiols forming a fluorescent product. By quantifying the fluorescence it is possible to determine the level of cellular thiols, and thus determine cell health.

Principle

In this application note, a method for measuring the cellular level of thiols is described. The cells to be investigated are mixed with **Solution 5**. The solution contains three different reagents: a stain staining all nucleated cells, a stain staining dead cells only and VitaBright-48™ which stains viable cells in an intensity-dependent manner depending on their level of thiols. A high fluorescence intensity of a particular cell indicates that the cell has a high level of thiols such as GSH.

The stained cells are immediately loaded into a NC-Slide: either the 2-chamber NC-Slide A2™ or the 8-chamber NC-Slide A8™. Samples are analyzed using the NC-3000 system. A fluorescence intensity histogram showing the distribution of thiol levels in all cells are

displayed on the PC screen. By comparing histograms of treated cells to controls the fraction of cells with low vitality (e.g. apoptotic or stressed cells) can be determined.

Procedures

If the cell line to be investigated is adherent or semi-adherent, then start by getting all cells into suspension using the preferred method of your laboratory (e.g. trypsin/EDTA treatment). Although the NucleoCounter® NC-3000™ is able to count aggregated cells, the accuracy is higher for single cell suspensions.

Materials needed

- Cells to be stained ^{1,2}
- **Solution 5** (VB-48-PI-AO)
- **NC-Slide A2™** or **NC-Slide A8™**

¹ Cells provided by the user.

² An untreated control should be included. Preferable, use logarithmically proliferating cells as control.

1. Pipette a representative cell sample from the cell suspension into a microcentrifuge tube. Add one volume of **Solution 5** into 20 volumes of the cell suspension. E.g., if the volume of the cell suspension is 190 µL then add 10 µL **Solution 5**. Mix by pipetting.
2. Engage NucleoCounter® NC-3000™ by starting the accompanying software.
3. Depending on the number of samples a 2-chamber slide (**NC-Slide A2™**) or an 8-chamber slide (**NC-Slide A8™**) can be used.
 - a. **NC-Slide A2™**: Load 30 µl of each of the cell suspensions into the chambers of the NC-Slide. Place the loaded NC-Slide on the tray of the NucleoCounter® NC-3000™ and select "**Vitality protocol**" and sample unit **NC-Slide A2™** and press RUN.
 - b. **NC-Slide A8™**: Load 8 µl of each of the cell suspensions into the chambers of the NC-Slide. Place the loaded NC-Slide on the tray of the NucleoCounter® NC-3000™ and select "**Vitality protocol**" and sample unit **NC-Slide A8™** and press RUN.

Cellular fluorescence is quantified by NC-3000 and the VitaBright-48 intensity is displayed in a histogram. Cells with a low level of thiols (e.g. apoptotic cells) also have a low intensity score, and are thus found in the lower end of the histogram. By placing the marker in the histogram, the user can divide the cell population into two (or more) subpopulations; e.g. divide the cell population into a fraction of cells with a low thiol level and another fraction with high thiol level, and if desired, an intermediate thiol level cell subpopulation. Typically the histogram shows two distinct populations, and in this case the marker is placed in between the two peaks (See example in **Figure 2**). The markers can be stored and retrieved in a new histogram. Thus, the position for the markers can be determined using an untreated control, and then retrieved for the treated cell sample. As nonviable cells always have a very low level of thiols, only the fluorescence intensity of viable cells is represented in the histogram.


Figure 2. Fluorescence intensity histogram for untreated MCF7 cells (upper left panel) and nocodazole treated MCF7 cells (upper right panel). Nocodazole treatment causes a decrease in the level of thiols in a subpopulation of cells (cells with low fluorescence intensity). The marker is placed in between the two populations in the first histogram (untreated MCF7 cells). Storing and retrieving the stored marker in the second histogram (nocodazole treated MCF7 cells) makes quantitative comparison possible. Below is fluorescence intensity histogram for healthy Jurkat cells (lower left panel) and starved Jurkat cells (lower right panel). Starvation causes a decrease in the thiol level of the cell population seen as a decrease in fluorescence intensity.

Note

To assure reliable results, it is recommended that the total cell concentration of the cell suspension should be in the range of $5 \cdot 10^4$ cells/mL to $5 \cdot 10^6$ cells/mL. If the concentration of cells is below $5 \cdot 10^4$ cells/mL then the cell concentration may be increased by centrifugation followed by resuspension of the pellet using growth media or PBS. The resuspended cell sample is then treated as described above.

If the total cell concentration is above $5 \cdot 10^6$ cells/mL, the cell suspension can be diluted with growth media or PBS to achieve the desired concentration. The diluted cell sample is then treated as described in the procedure.