


# How can a homogeneous H<sub>2</sub>O<sub>2</sub> Distribution enable time- and load-independent Decontamination in a Material Airlock

Comparative Analysis of Static vs. Active Distribution with  
Puro™ MagLev

# Room Layout of Filling Line



A-Zone

B-Zone

D-Zone


# Material Transfer Chambers (MTCs)

## The Gateway to Grade A/B

- Material Transfer Chambers are critical for the safe transfer of tools, materials, and equipment between cleanliness zones.
- Prevents cross-contamination.
- “Sterilizes” surfaces via VHP (vaporized  $H_2O_2$ ).
- Validated decontamination process ensures efficient surface decontamination of all validated loading patterns.



# Load: Empty Cart



View on right side of cart




View on front side of cart



View from D-Zone

# Load: Production (D)



View on right side of cart



View on left side of cart




View from D-Zone

# Material Transfer Chambers (MTCs)

## The crucial Point – System Design

- Load-independent decontamination process depends on homogeneous  $\text{H}_2\text{O}_2$  distribution, generated by turbulences.
- Distribution depends on design.
  - If given, the load is flexible.
    - Amount of items of validated loading patterns can be reduced
    - High flexibility
    - Adaptive to production need
  - If not given, the load is fixed and not variable.
    - Validated loading patterns are allowed to use, but
    - No flexibility in reducing amount of loading items.
    - Challenging for production, since loading can not reflect availability of items and need for production process.



# How to improve $\text{H}_2\text{O}_2$ Distribution

## Active Turbulence via Puro

### Why Puro?

- does not require lubricants or sealings
- no mechanical bearing and thus no friction
- no mechanical contact between fan rotor and motor, therefore, no wear and no particle generation
- Maintenance-free
- non-absorbing high-performance plastic
- no other material in contact with  $\text{H}_2\text{O}_2$
- can be operated at extremely high speeds



# How to evaluate H<sub>2</sub>O<sub>2</sub> Distribution

- Biological Indicators (BIs)  
=> provide  $\geq 1.0 \times 10^6$  spores  
=> yes/no growth ( $\geq 6$  log-reduction)  
=> no/growth ( $< 6$  log-reduction)
- Chemical Indicators (CIs)  
measure log-reduction  
=> graduated result  
=> spectrophotometer reader
- Calculation of Distribution (Range):  
Max log-reduction – Min log-reduction



## Bioquell CI Bericht

Benutzername: Klein, Leonie  
Gerätekalibrierung gültig bis: 25.07.2025

Berichtsdatum: 16.01.2026  
Gerät S/N: 2010007223

Bioquell Isolator CI  
Chargennummer: 25261L1273

Verfallsdatum: 02.09.2026

| Name      | Datum      | Zeit  | dE    | Log Reduzierung | Ergebnis |
|-----------|------------|-------|-------|-----------------|----------|
| Sample 1  | 16.01.2026 | 11:06 | 16.11 | 4.95            | FAIL     |
| Sample 2  | 16.01.2026 | 11:06 | 13.01 | 5.52            | FAIL     |
| Sample 3  | 16.01.2026 | 11:06 | 17.02 | 4.78            | FAIL     |
| Sample 4  | 16.01.2026 | 11:06 | 16.85 | 4.81            | FAIL     |
| Sample 5  | 16.01.2026 | 11:07 | 13.99 | 5.34            | FAIL     |
| Sample 6  | 16.01.2026 | 11:07 | 18.54 | 4.50            | FAIL     |
| Sample 7  | 16.01.2026 | 11:07 | 16.78 | 4.82            | FAIL     |
| Sample 8  | 16.01.2026 | 11:07 | 13.97 | 5.34            | FAIL     |
| Sample 9  | 16.01.2026 | 11:07 | 18.13 | 4.57            | FAIL     |
| Sample 10 | 16.01.2026 | 11:07 | 15.56 | 5.05            | FAIL     |
| Sample 11 | 16.01.2026 | 11:07 | 16.46 | 4.88            | FAIL     |
| Sample 12 | 16.01.2026 | 11:08 | 16.29 | 4.91            | FAIL     |
| Sample 13 | 16.01.2026 | 11:08 | 17.73 | 4.65            | FAIL     |
| Sample 14 | 16.01.2026 | 11:08 | 19.01 | 4.41            | FAIL     |
| Sample 15 | 16.01.2026 | 11:08 | 16.18 | 4.93            | FAIL     |
| Sample 16 | 16.01.2026 | 11:08 | 15.46 | 5.06            | FAIL     |
| Sample 17 | 16.01.2026 | 11:08 | 18.51 | 4.50            | FAIL     |
| Sample 18 | 16.01.2026 | 11:09 | 16.95 | 4.79            | FAIL     |
| Sample 19 | 16.01.2026 | 11:09 | 15.84 | 4.99            | FAIL     |
| Sample 20 | 16.01.2026 | 11:09 | 19.42 | 4.33            | FAIL     |
| Sample 21 | 16.01.2026 | 11:10 | 16.51 | 4.87            | FAIL     |
| Sample 22 | 16.01.2026 | 11:10 | 17.85 | 4.62            | FAIL     |
| Sample 23 | 16.01.2026 | 11:10 | 15.89 | 4.99            | FAIL     |
| Sample 24 | 16.01.2026 | 11:11 | 15.82 | 5.00            | FAIL     |
| Sample 25 | 16.01.2026 | 11:11 | 14.38 | 5.26            | FAIL     |



# Test Design

Jan 2026: Testruns with Puro

Nov 2025: Testruns without Puro  
(N/A or Reference Runs)

## 2 loading patterns

empty cart / production load (D)

## Number of tested locations (CIs)

26 empty cart

33 production load (D)

## 2 decontamination phase durations

7 min / 13 min

## 2 Puro positions

floor / side wall

## Puro operating modes

continuous / iterative

## Puro programmes


contin. -4000 rpm

contin. +4000 rpm


iterative 3 s -8000 rpm / 10 s pause

iterative 10 s -8000 rpm / 10 s pause

iterative 10 s +8000 rpm / 10 s -8000 rpm



# Remarks



- Puro on side wall only tested in 3 runs (13 min) with production load.

(contin. -4000 rpm | iterative 3 s -8000 rpm / 10 s pause | contin. +4000 rpm)

For reasons of space, installation would not be feasible in the future.

- Duration of production deco phase: 24 min


Deco phase durations for testing (7 min / 13 min) were designed to not get a full kill to reveal differences.

- The key learnings on the combination of Puro, load and deco phase duration will be summarized on the next 2 slides for a clear understanding.
- Further results from the tests can be found in the back-up slides.

# Comparison of H<sub>2</sub>O<sub>2</sub> distribution w/wo Puro (on the floor)



- In the absence of Puro, the distribution is proven to be dependent on the load.
- Best load independence shown for iterative operating mode.
- Cont. operating mode is load dependent.
- General: H<sub>2</sub>O<sub>2</sub> distribution (range) is poorer with production load (D) than at minimum load (empty cart)

Range for different Puro settings  
7 min runs



# Comparison of H<sub>2</sub>O<sub>2</sub> distribution w/wo Puro (on the floor)

- Impact of time: Distribution (range) is poorer at minimum load (empty cart) than at maximum load (D).
- Volume of load items causes turbulences (Best Case). Reason for need to stick to the defined loading pattern.
- With Puro, the distribution at minimum load is better for both settings than at maximum load (opposite of the reference)
- For both Puro settings tested, the differences in distribution between minimum and maximum load are closer to each other.
- Alternating mode of the Puro has the “worst” distribution.



# Summary

## Reference Runs (N/A) without Puro

- As the duration increases, the difference in the distribution between minimum and maximum load is reversed.
- The minimum load represents the worst-case scenario due to the poor distribution.
- Conclusion: Distribution depends on both time and load

## Puro Runs

- Iterative mode (3s -8000 rpm / 10 s pause) produced the most uniform distribution for both maximum and minimum loading, as well as for the two deco times.
- Slight differences in distribution (maximum vs. minimum loading) can be eliminated through cycle duration.
- Conclusion: Iterative operating mode of Puro offers load-independent distribution by means of an appropriate deco phase duration.

# Conclusion

**How can a homogeneous H<sub>2</sub>O<sub>2</sub> Distribution enable time- and load-independent Decontamination in a Material Airlock?**

- By using the Puro, this can be achieved.
- Based on the available results, the validation of a new process can take place.
- This makes it as flexible as possible to adapt the loading to the needs of production.
- The homogeneous distribution through the Puro with an iterative program enables a load-independent and largely time-independent process that meets the regulatory requirements for spore log-reduction.

# Special Thanks to...

- the whole Levitronix team  
especially...
- Ivana Ferstera (online and site support)
- Adrian Ljusic (online and site support)
- Sabrina Müller (excellent organization)
- Stefan Girlich (executing tests in Nov 2025)
- Leonie Klein (executing tests in Jan 2026)

# Contact

---

- **Dr. Simone Bläsi**  
H<sub>2</sub>O<sub>2</sub> Validation Expert (MSAT Steriles)
- **EMAIL**  
simone.blaesi@roche.com
- **ORGANISATION**  
Roche - MSAT Steriles, Kaiseraugst
- [contact@levitronix.com](mailto:contact@levitronix.com)
- [www.levitronix.com](http://www.levitronix.com)



# Back-up Slides

## Digging deeper into data....

# Overview of Testruns


| Execution Date | Run No. / Deco Dur. | Load       | Position Puro | Programme Puro        |
|----------------|---------------------|------------|---------------|-----------------------|
| 29-Nov-2025    | Lauf 8 / 13 min     | D)         | N/A           | N/A                   |
| 15-Jan-2026    | Lauf 1 / 13 min     | D)         | Floor         | -4000                 |
| 15-Jan-2026    | Lauf 2 / 13 min     | D)         | Floor         | 3s -8000 / 10s Pause  |
| 15-Jan-2026    | Lauf 3 / 13 min     | D)         | Floor         | +4000                 |
| 16-Jan-2026    | Lauf 4 / 13 min     | D)         | Wall          | -4000                 |
| 16-Jan-2026    | Lauf 5 / 13 min     | D)         | Wall          | 3s -8000 / 10s Pause  |
| 16-Jan-2026    | Lauf 6 / 13 min     | D)         | Wall          | +4000                 |
|                |                     |            |               |                       |
| 20-Jan-2026    | Lauf 7 / 13 min     | D)         | Floor         | 10s -8000 / 10s Pause |
| 20-Jan-2026    | Lauf 8 / 13 min     | D)         | Floor         | -8000                 |
| 20-Jan-2026    | Lauf 9 / 13 min     | D)         | Floor         | 10s +8000 / 10s -8000 |
|                |                     |            |               |                       |
| 29-Nov-2025    | Lauf 7 / 7 min      | D)         | N/A           | N/A                   |
| 21-Jan-2026    | Lauf 10 / 7 min     | D)         | Floor         | 3s -8000 / 10s Pause  |
| 21-Jan-2026    | Lauf 11 / 7 min     | D)         | Floor         | -4000                 |
|                |                     |            |               |                       |
| 08-Nov-2025    | Lauf 1 / 7 min      | Empty Cart | N/A           | N/A                   |
| 22-Jan-2026    | Lauf 12 / 7 min     | Empty Cart | Floor         | -4000                 |
| 22-Jan-2026    | Lauf 13 / 7 min     | Empty Cart | Floor         | 3s -8000 / 10s Pause  |
|                |                     |            |               |                       |
| 08-Nov-2025    | Lauf 2 / 13 min     | Empty Cart | N/A           | N/A                   |
| 23-Jan-2026    | Lauf 14 / 13 min    | Empty Cart | Floor         | -4000                 |
| 23-Jan-2026    | Lauf 15 / 13 min    | Empty Cart | Floor         | 3s -8000 / 10s Pause  |

Overview  
of the runs performed (Jan 2026)  
and their corresponding reference runs  
(Nov 2025)

# Puro Position and Loading




Puro on floor  
empty cart




Puro on side wall  
Production load (D)




# Comparison Puro Position




# Comparision Puro Program



# Load + Deco Time: Impact on Distribution



# Conclusion for 7 min




## In general

- The distribution (range) is poorer at maximum load (D) than at minimum load (empty cart)

## With Puro vs without Puro

- When set to 3 s at -8000 rpm / 10 s pause, the difference is smaller compared to the reference run.
- at -4000 rpm, the difference is more significant compared to the reference run

# Conclusion for 13 min



## Reference Runs


- The distribution (range) is poorer at minimum load (empty cart) than at maximum load (D).

## With Puro vs without Puro

- With Puro, the distribution at minimum load is better for both settings than at maximum load (the opposite of the reference)
- For both Puro settings tested, the differences in distribution between minimum and maximum load are closer to each other.

# Summary

## Load and Deco Duration



### Reference Runs

As the duration increases, the difference in the distribution between minimum and maximum load is reversed.

The minimum load represents the worst-case scenario due to the poor distribution.

Conclusion: the distribution depends on both time and load

### Puro Programs

3s -8000 rpm / 10 s Pause produced the most uniform distribution for both maximum and minimum loading, as well as for the two deco times.

Conclusion: the best setting for reaching a distribution independent of time and loading