

Levitronix® Standard Article List and Order Information
LEVIFLOW® LFSC-D Series, Sets, Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note
LeviFlow Sets - Flowsensor with Converter			
LFSC-05D-001 + LFC-1C-PC (1/8" x 1/4" Tubing, 0-1 l/min, 4 Calibration Sets) Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-91113	LEVIFLOW LFSC Flowmeter	<i>Default activated calibration set is Silicone for 37°C water. Adaptor/extension cables to be ordered as separate article.</i>
LFSC-08D-001 + LFC-1C-PC (1/4" x 3/8" Tubing, 0-4 l/min, 4 Calibration Sets) Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-91076	LEVIFLOW LFSC Flowmeter	
LFSC-12D-007 + LFC-1C-PC (3/8" x 9/16" Tubing, 0-20 l/min, 4 Calibration Sets) Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-91012	LEVIFLOW LFSC Flowmeter	
LFSC-22D-005 + LFC-1C-PC (3/4" x 1" Tubing, 0-80 l/min, 4 Calibration Sets) Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-91013	LEVIFLOW LFSC Flowmeter	
LFSC-30D-001 + LFC-1C-PC (1" x 13/8" Tubing, 0-160 l/min, 4 Calibration Sets) Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-91183	LEVIFLOW LFSC Flowmeter	
LEVIFLOW® Articles for Configurations with 6-Channel Converters LFC-6C/LFC-6CIO			
LeviFlow Converter LFC-6CIO-PC (24VDC, 6-Chan., RS485, Analog/Dig. IO's) <i>Note 1: Connector and Termination Box to be ordered as separate article. Note 2: EMI ferrite for Termination wires to be ordered as separate article.</i>	100-30446	LEVIFLOW LFSC Flowmeter	<i>Article packages with flow sensors available on request.</i>
Accessories			
LeviFlow Interconnect Cable LFI-C.1-10 (1m, PVC) LeviFlow Interconnect Cable LFI-C.1-30 (3m, PVC) LeviFlow Interconnect Cable LFI-C.1-60 (6m, PVC)	190-10307 190-10308 190-10309	Accessories LEVIFLOW	
LeviFlow Adaptor Extension Cable LFE-C.2-10 (1m, PVC, Wall-Mountable Connector) LeviFlow Adaptor Extension Cable LFE-C.2-30 (3m, PVC, Wall-Mountable Connector) LeviFlow Adaptor Extension Cable LFE-C.2-60 (6m, PVC, Wall-Mountable Connector)	190-10310 190-10311 190-10312	Accessories LEVIFLOW	
LeviFlow Converter Starter-Kit LFC-1C-PC-SK (Converter, AC/DC Desktop Power Supply, Cabling, Config.-Software)	100-91072	Accessories LEVIFLOW	
LeviFlow Connector Box for LFC-6CIO-PC (for wiring)	100-30447	Accessories LEVIFLOW	<i>For LFC-6CIO-PC.</i>
LeviFlow Termination Box for LFC-6C (for termination of RS485)	100-30317	Accessories LEVIFLOW	<i>For LFC-6CIO-PC.</i>
LeviFlow Splitting Ferrite	100-30353	Accessories LEVIFLOW	<i>Is recommended for all sensor cables, when stacking multiple LFC-6CIO-PC in order to avoid crosstalk.</i>
USB-RS485 YN-485I-TR (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	100-30392	Accessories IPS	
Spare Parts			
LeviFlow Clamp-On Sensor LFSC-05D-001 (1/8"x1/4" Tubing, 0-1 l/min, 4 Calibr. Sets) Two Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-30407	Spare Parts LEVIFLOW	<i>Default activated calibration set is Silicone for 37°C water.</i>
LeviFlow Clamp-On Sensor LFSC-08D-001 (1/4"x3/8" Tubing, 0-4 l/min, 4 Calibr. Sets) Two Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-30396	Spare Parts LEVIFLOW	<i>Default activated calibration set is Silicone for 37°C water.</i>
LeviFlow Clamp-On Sensor LFSC-12D-007 (3/8"x9/16" Tubing, 0-20 l/min, 4 Calibr. Sets) Two Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-30390	Spare Parts LEVIFLOW	<i>Default activated calibration set is Silicone for 37°C water.</i>
LeviFlow Clamp-On Sensor LFSC-22D-005 (3/4"x1" Tubing, 0-80 l/min, 4 Calibr. Sets) Two Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-30391	Spare Parts LEVIFLOW	<i>Default activated calibration set is Silicone for 37°C water.</i>
LeviFlow Clamp-On Sensor LFSC-30D-001 (1" x 13/8" Tubing, 0-160 l/min, 4 Calibr. Sets) Two Calibration for Silicone and C-Flex Tubing with Water at 37°C and 20°C	100-30414	Spare Parts LEVIFLOW	<i>Default activated calibration set is Silicone for 37°C water.</i>
LeviFlow Converter LFC-1C-PC (24VDC, Single Channel, 2 Digital Outputs, 1 Digital Input, 1 Analog Output, RS485)	100-30374	Spare Parts LEVIFLOW	<i>Includes an EMI ferrite (article # 100-30353) for the flow sensor cable.</i>
LeviFlow Converter LFC-6CIO-PC (24VDC, 6-Channel, RS485, Analog/Digital IO's)	100-30446	LEVIFLOW LFSC Flowmeter	
Connector for LFC-6CIO / Analog Out	100-30351	Spare Parts LEVIFLOW	
Connector for LFC-6CIO / Digital I/O	100-30352	Spare Parts LEVIFLOW	