

**Levitronix® Standard Article List and Order Information**  
**Maglev Pumpsystem PuraLev® i30, Accessories and Spare Parts**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard Driver System Configurations</b>			
<b>Pump Driver System PLD-i30SU.1</b> (OEM Model) Integrated Pump Driver IPD-30.1-50-02 (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-90955	PuraLev SU	Note: Mounting Base Plate MBP-i30.1 has to be ordered separately.
<b>Pump Driver System PLD-i30SU.6</b> (OEM Model, ATEX/IECEX) Integrated Pump Driver IPD-30.12-50-02 (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91514	PuraLev SU	Note: Mounting Base Plate MBP-i30.1 has to be ordered separately.
<b>Pump Driver System PLD-i30SU.2</b> (EasyConnect Model) Integrated Pump Driver IPD-30.3-02 (24VDC, Epoxy Coating, PLC and 2x Fieldbus, incl. MBP-i30.1) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91025	PuraLev SU	
<b>Pump Driver System PLD-i30SU.3</b> (Stand-Alone Model) Integrated Pump Driver IPD-30.5-02 (24VDC, Epoxy Coating, Display with Buttons, PLC, incl. MBP-i30.1) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91026	PuraLev SU	
<b>Single-Use Pump Heads</b>			
<b>Single-Use Pumphead DCP-30.2</b> (Barb 3/8", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile) <b>Single-Use Pumphead DCP-30.1</b> (Triclamp 3/8", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile)	100-90968 100-90959	PuraLev SU	
<b>Single-Use Pump Head Gamma Irradiated DCP-30.2-G25</b> (Applied Gamma Dosage ≥ 25 kGy) <b>Single-Use Pump Head Gamma Irradiated DCP-30.1-G25</b> (Applied Gamma Dosage ≥ 25 kGy)	100-91071 100-91170	PuraLev SU	
<b>Single-Use Pump Head Sterile Fittings DCP-30.2-SF1-G25</b> (Sterile Fittings, Gamma Dosage ≥ 25 kGy)	100-91234	PuraLev SU	
<b>Accessories</b>			
<b>Mounting Base Plate MBP-i30.1</b> (PP+GF)	190-10313	Accessories IPS	
<b>Power Supply TPC 055-124</b> (Traco, 24V / 55W, DIN Rail Mounting)	100-40014	Accessories IPS	
<b>Desktop AC/DC Power Supply VEC50US24 HR30</b> (24VDC/50W, 90-264VAC, Hirose Connector)	100-40015	Accessories IPS	AC mains cable ACM-1.x for various countries to be ordered as separate article.
<b>AC Mains Cable AMC-1.1</b> (2m, US, Canada) <b>AC Mains Cable AMC-1.2</b> (2.5m, GE, DK, NO, FI, BE, NL, SE, AT) <b>AC Mains Cable AMC-1.3</b> (2.5m, Japan) <b>AC Mains Cable AMC-1.4</b> (2.5m, Switzerland) <b>AC Mains Cable AMC-1.5</b> (2.5m, United Kingdom)	190-10331 190-10332 190-10333 190-10334 190-10335	Accessories IPS	AC cable for desktop power supply VEC50US24.
<b>USB-RS485 YN-485I-TR</b> (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	100-30392	Accessories IPS	
<b>IPS Cable Power 3 Wires ICP-1.1-50</b> (5m PVC, Circular Driver Connector to Open Wires)	190-10342	Accessories IPS	For power supplies with open wire connection.
<b>IPS Cable Signal 6 Wires ICS-1.1-01</b> (0.1m PVC, Circular Driver Connector to Circular Connector) <b>IPS Cable Signal 6 Wires ICS-1.1-10</b> (1m) <b>IPS Cable Signal 6 Wires ICS-1.1-30</b> (3m)	190-10343 190-10344 190-10345	Accessories IPS	For connecting fieldbus of EasyConnect driver arrays.
<b>IPS Cable Signal 6 Wires ICS-1.2-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10346	Accessories IPS	For connecting devices to PLC of Stand-Alone driver or fieldbus of EasyConnect driver.
<b>IPS Cable Signal 12 Wires ICS-2.1-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connecting devices to PLC of EasyConnect driver.
<b>Fieldbus Termination Connector FTC-1.1</b> (PPS with FPM Sealing)	190-10348	Accessories IPS	Termination of fieldbus of EasyConnect drivers.
<b>ATEX Cable Sealing System ACS-A.1</b>	100-90292	Accessories IPS	
<b>Spare Parts</b>			
<b>Integrated Pump Driver IPD-30.1-50-02</b> (OEM Model) (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating)	100-10088	Spare Parts IPS	
<b>Integrated Pump Driver IPD-30.12-50-02</b> (OEM Model, ATEX/IECEX) (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating)	100-10203	Spare Parts IPS	
<b>Integrated Pump Driver IPD-30.3-02</b> (EasyConnect Model) (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus, includes Mounting Base Plate MBP-i30.1)	100-10097	Spare Parts IPS	
<b>Integrated Pump Driver IPD-30.5-02</b> (Stand-Alone Model) (24VDC, Epoxy Coating, Display with Buttons, PLC, includes Mounting Base Plate MBP-i30.1)	100-10098	Spare Parts IPS	
<b>Pump Head Socket PHS-i30.1</b> (Anodized Aluminum, Bayonet Type)	100-90947	Spare Parts IPS	