

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® 200MU (Multi-Use), Accessories and Spare Parts

| Article Name (and Description) | Part No. | Article Category | Note |
|--|------------------------------|------------------|---|
| Standard Systems | | | |
| Pump System PuraLev-200MU.26 (Standalone) Pumphead LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA Impeller, EPDM O-Ring, With Drain Port) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Controller LPC-200.1-02 (24V/200W, Integrated Display, PLC) Note: Direct cabling between motor and controller. | 100-90865 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.27 (PLC/USB) or PuraLev-200MU.28 (PLC/RS232) Pumphead LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA Impeller, EPDM O-Ring, With Drain Port) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Controller LPC-200.2-02 (24V/200W, PLC, USB) or LPC-200.3-02 (PLC, RS232) Note: Direct cabling between motor and controller. | 100-90866 or 100-90867 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.29 (Standalone) Pumphead LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA Impeller, EPDM, No Drain Port) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Controller LPC-200.1-02 (24V/200W, integrated Display, Enable Connector) Note: Direct cabling between motor and controller. | 100-90868 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.30 (PLC/USB) or PuraLev-200MU.31 (PLC/RS232) Pumphead LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA Impeller, EPDM O-Ring, No Drain Port) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Controller LPC-200.2-02 (24V/200W, PLC, USB) or LPC-200.3-02 (PLC, RS232) Note: Direct cabling between motor and controller. | 100-90869 or 100-90870 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.32 (ATEX, Standalone) Pumphead LPP-200.17 (PVDF Housing Tricl. 1/2", PFA Impeller, EPDM O-Ring, With Drain Port) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Connector, ETFE Housing Coating) Controller LPC-200.1-02 (24V/200W, integrated Display, Enable Connector) To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m | 100-90871 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.33 (ATEX, PLC/USB) or PuraLev-200MU.34 (ATEX, PLC/RS232) Pumphead LPP-200.17 (PVDF Housing Tricl. 1/2", PFA Impeller, EPDM O-Ring, With Drain Port) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Connector, ETFE Housing Coating) Controller LPC-200.2-02 (24V/200W, PLC, USB) or LPC-200.3-02 (PLC, RS232) To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m | 100-90872 or 100-90873 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.35 (ATEX, Standalone) Pumphead LPP-200.16 (PVDF Housing Tricl. 1/2", PFA Impeller, EPDM O-Ring, No Drain Port) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Connector, ETFE Housing Coating) Controller LPC-200.1-02 (24V/200W, integrated Display, Enable Connector) To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m | 100-90874 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |
| Pump System PuraLev-200MU.36 (ATEX, PLC/USB) or PuraLev-200MU.37 (ATEX, PLC/RS232) Pumphead LPP-200.16 (PVDF Housing Tricl. 1/2", PFA Impeller, Silicon O-Ring, No Drain Port) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Connector, ETFE Housing Coating) Controller LPC-200.2-02 (24V/200W, PLC, USB) or LPC-200.3-02 (PLC, RS232) To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m | 100-90875 or 100-90876 | PuraLev MU | <i>With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.</i> |

| Article Name (and Description) | Part No. | Article Category | Note |
|--|-----------|-----------------------------|--|
| Accessories | | | |
| Power Supply TSP 360-124 (Traco, 24V / 360W) | 100-40016 | Accessories Pump Systems | |
| Power Supply TSP 360-148 (Traco, 48V / 360W) | 100-40017 | Accessories Pump Systems | |
| Mounting Base Plate MBP-1.1 (PVDF) | 190-10004 | Accessories Pump Systems | |
| Air Cooling Module ACM-1.1 (PVDF, NPT1/4") | 190-10003 | Accessories Pump Systems | Discontinued by September 2023. Alternative is ACM-1.3. |
| Air Cooling Module ACM-1.3 (PP, NPT1/4") | 190-10193 | Accessories Pump | |
| ATEX Cable Sealing System ACS-A.1 | 100-90292 | Accessories Pump Systems | |
| Autoclaving Reinforcing Tool ART-200.1 (Anodized Alu., reinforcement of pump head during Autoclaving) | 190-10280 | Accessories Pump Systems | |
| User Panel LUI-A.1 (RS232 interface) | 100-30300 | Accessories Pump Systems | Used with LPC-200.3-02 |
| Adapter/Extension Cable MCA-1.2-05 (0.5m, FEP-Jacket, Circular/D-SUB) | 190-10147 | | |
| Adapter/Extension Cable MCA-1.2-30 (3m) | 190-10092 | | |
| Adapter/Extension Cable MCA-1.2-50 (5m) | 190-10128 | Accessories Pump Systems | |
| Adapter/Extension Cable MCA-1.2-70 (7m) | 190-10148 | | |
| Adapter/Extension Cable MCA-1.2-100 (10m) | 190-10149 | | |
| Adapter/Extension Cable MCA-1.3-05 (0.5m, FEP-Jacket, Circular Wallmountable/D-SUB) | 190-10234 | | |
| Adapter/Extension Cable MCA-1.3-30 (3m) | 190-10237 | | |
| Adapter/Extension Cable MCA-1.3-50 (5m) | 190-10200 | Accessories Pump Systems | |
| Adapter/Extension Cable MCA-1.3-70 (7m) | 190-10235 | | |
| Adapter/Extension Cable MCA-1.3-100 (10m) | 190-10236 | | |
| Adapter/Extension Cable MCA-1.4-05 (0.5m, PVC-Jacket, Circular/D-SUB) | 190-10228 | | |
| Adapter/Extension Cable MCA-1.4-30 (3m) | 190-10170 | | |
| Adapter/Extension Cable MCA-1.4-50 (5m) | 190-10187 | Accessories Pump Systems | |
| Adapter/Extension Cable MCA-1.4-70 (7m) | 190-10229 | | |
| Adapter/Extension Cable MCA-1.4-100 (10m) | 190-10230 | | |
| Adapter/Extension Cable MCA-1.5-05 (0.5m, PVC-Jacket, Circular Wallmountable/D-SUB) | 190-10225 | | |
| Adapter/Extension Cable MCA-1.5-30 (3m) | 190-10231 | | |
| Adapter/Extension Cable MCA-1.5-50 (5m) | 190-10199 | Accessories Pump Systems | |
| Adapter/Extension Cable MCA-1.5-70 (7m) | 190-10232 | | |
| Adapter/Extension Cable MCA-1.5-100 (10m) | 190-10233 | | |
| Spare Parts | | | |
| Motor BSM-1.4 (ETFE Coating, 5m FEP Cable with D-SUB Connector) | 100-10005 | Spare Parts Pump Systems | |
| Motor BSM-1.6 (ATEX) (ETFE Coating, 2m FEP Cable with Circular Connector, ATEX) | 100-10063 | Spare Parts Pump Systems | |
| Motor BSM-1.8 (ATEX) (Epoxy Coating, 2m PVC Cable with Circular Connector, ATEX) | 100-10090 | Spare Parts Pump Systems | |
| Pump Head LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA-Impeller LPI-200.5, EPDM O-Ring, No Drain Port) | 100-90863 | Spare Parts Pump Systems | With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP. |
| Pump Head LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA-Impeller LPI-200.5, EPDM O-Ring, With Drain Port) | 100-90864 | Spare Parts Pump Systems | With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP. |
| Controller LPC-200.1-02 (24V/200W, Display, Enable connector) (included: Enable connector) | 100-30030 | Spare Parts Pump Systems | |
| Controller LPC-200.2-02 (24V/200W, USB, PLC interface) (included: PLC connector) | 100-30031 | Spare Parts Pump Systems | |
| Controller LPC-200.3-02 (24V/200W, RS232, PLC interface) (included: PLC connector) | 100-30032 | Spare Parts Pump Systems | |