

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® 200SU (Single-Use), Accessories and Spare Parts

| Article Name (and Description) | Part No. | Article Category |
|---|--|-----------------------------|
| Driver System Packages (Motor with Pumphead Socket and Controller) | | |
| PuraLev Driver System PLD-200SU.1 (Standalone) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Type) Controller LPC-200.1-01 (24V/200W, Integrated Display, PLC) | 100-90647 | PuraLev SU |
| PuraLev Driver System PLD-200SU.2 (PLC/USB) or PLD-200SU.3 (PLC/RS232) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Type) Controller LPC-200.2-01 (24V/200W, PLC, USB) or LPC-200.3-01 (PLC, RS232) | 100-90648 or 100-90649 | PuraLev SU |
| PuraLev Driver System PLD-200SU.4 (ATEX, Standalone) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Hummel Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Type) Controller LPC-200.1-01 (24V/200W, Integrated Display, PLC) | 100-90650 | PuraLev SU |
| PuraLev Driver System PLD-200SU.5 (PLC/USB) or PLD-200SU.6 (PLC/RS232) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Hummel Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Type) Controller LPC-200.2-01 (24V/200W, PLC, USB) or LPC-200.3-01 (PLC, RS232) | 100-90651 or 100-90652 | PuraLev SU |
| Single-Use Pumpheads | | |
| Single-Use Pumphead DCP-200.2 Welded and Molded PP Housing, Triclamp 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile | 100-90734 | PuraLev SU |
| Single-Use Pump Head Gamma Irradiated DCP-200.2-G25 (Applied Gamma Dosage ≥ 25 kGy) | 100-91078 | PuraLev SU |
| Single-Use Pumphead DCP-200.3 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile | 100-90792 | PuraLev SU |
| Single-Use Pump Head Gamma Irradiated DCP-200.3-G25 (Applied Gamma Dosage ≥ 25 kGy) | 100-91122 | PuraLev SU |
| Single-Use Pumphead DCP-200.3-SF1-G25 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Sterile Fittings, Applied Gamma Dosage ≥ 25 kGy) | 100-91233 | PuraLev SU |
| Accessories | | |
| Power Supply TSP 360-124 (Traco, 24V / 360W) | 100-40016 | Accessories Pump Systems |
| Power Supply TSP 360-148 (Traco, 48V / 360W) | 100-40017 | Accessories Pump Systems |
| Mounting Base Plate MBP-1.1 (PVDF) | 190-10004 | Accessories Pump Systems |
| Air Cooling Module ACM-1.1 (PVDF, NPT1/4") | 190-10003 | Accessories Pump Systems |
| Air Cooling Module ACM-1.3 (PP, NPT1/4") | 190-10193 | Accessories Pump Systems |
| ATEX Cable Sealing System ACS-A.1 | 100-90292 | Accessories Pump Systems |
| User Panel LUI-A.1 (RS232 interface) | 100-30300 | Accessories Pump Systems |
| Adapter/Extension Cable MCA-1.2-05 (0.5m, FEP-Jacket, Circular/D-SUB) Adapter/Extension Cable MCA-1.2-30 (3m) Adapter/Extension Cable MCA-1.2-50 (5m) Adapter/Extension Cable MCA-1.2-70 (7m) Adapter/Extension Cable MCA-1.2-100 (10m) | 190-10147 190-10092 190-10128 190-10148 190-10149 | Accessories Pump Systems |
| Adapter/Extension Cable MCA-1.3-05 (0.5m, FEP-Jacket, Circular Wallmountable/D-SUB) Adapter/Extension Cable MCA-1.3-30 (3m) Adapter/Extension Cable MCA-1.3-50 (5m) Adapter/Extension Cable MCA-1.3-70 (7m) Adapter/Extension Cable MCA-1.3-100 (10m) | 190-10234 190-10237 190-10200 190-10235 190-10236 | Accessories Pump Systems |
| Adapter/Extension Cable MCA-1.4-05 (0.5m, PVC-Jacket, Circular/D-SUB) Adapter/Extension Cable MCA-1.4-30 (3m) Adapter/Extension Cable MCA-1.4-50 (5m) Adapter/Extension Cable MCA-1.4-70 (7m) Adapter/Extension Cable MCA-1.4-100 (10m) | 190-10228 190-10170 190-10187 190-10229 190-10230 | Accessories Pump Systems |
| Adapter/Extension Cable MCA-1.5-05 (0.5m, PVC-Jacket, Circular Wallmountable/D-SUB) Adapter/Extension Cable MCA-1.5-30 (3m) Adapter/Extension Cable MCA-1.5-50 (5m) Adapter/Extension Cable MCA-1.5-70 (7m) Adapter/Extension Cable MCA-1.5-100 (10m) | 190-10225 190-10231 190-10199 190-10232 190-10233 | Accessories Pump Systems |
| Spare Parts | | |
| Motor BSM-1.4 (ETFE Coating, 5m FEP Cable with D-SUB Connector) | 100-10005 | Spare Parts Pump Systems |
| Motor BSM-1.6 (ATEX) (ETFE Coating, 2m FEP Cable with Circular Connector, ATEX) | 100-10063 | Spare Parts Pump Systems |
| Motor BSM-1.8 (ATEX) (Epoxy Coating, 2m PVC Cable with Circular Connector, ATEX) | 100-10090 | Spare Parts Pump Systems |
| Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Socket Type) | 100-90695 | Spare Parts Pump Systems |
| Controller LPC-200.1-01 (24V/200W, Display, Enable connector) (included: Enable connector) | 100-30027 | Spare Parts Pump Systems |
| Controller LPC-200.2-01 (24V/200W, USB, PLC interface) (included: PLC connector) | 100-30028 | Spare Parts Pump Systems |
| Controller LPC-200.3-01 (24V/200W, RS232, PLC interface) (included: PLC connector) | 100-30029 | Spare Parts Pump Systems |

**on
and Spare Parts**

Note

Note: Motor has direct cabling to controller.

Note: Motor has direct cabling to controller.

*Note: To be ordered as separate article:Extension
Cables MCA-1.x-yy -> 0.5-10m.*

*Note: To be ordered as separate article:Extension
Cables MCA-1.x-yy -> 0.5-10m.*

| |
|--|
| |
| |
| |
| |
| |

| |
|--|
| |
| |
| |

*Discontinued by September 2023.
Alternative is ACM-1.3.*

| |
|--|
| |
| |

Used with LPC-200.3-01

| |
|--|
| |
| |
| |
| |
| |

| |
|--|
| |
| |
| |
| |
| |
| |