

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® 600MU (Multi-Use), Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note
Standard Systems			
Pump System PuraLev-600MU.1 (ATEX) Pump Head LPP-600.7 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, With Drain Port) Motor LPM-600.5 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.1-02 (48V/600W, Integrated Display, Enable Connector)	100-90590	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.2 (ATEX, USB/PLC) or PuraLev-600MU.3 (ATEX, RS232/PLC) Pump Head LPP-600.7 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, With Drain Port) Motor LPM-600.5 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.2-02 (48V/600W, USB/PLC) or LPC-600.3-02 (RS232/PLC)	100-90591 or 100-90592	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.4 (ATEX) Pump Head LPP-600.7 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, With Drain Port) Motor LPM-600.4 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.1-02 (48V/600W, Integrated Display, Enable Connector)	100-90593	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.5 (ATEX, USB/PLC) or PuraLev-600MU.6 (ATEX, RS232/PLC) Pump Head LPP-600.7 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, With Drain Port) Motor LPM-600.4 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.2-02 (48V/600W, USB/PLC) or LPC-600.3-02 (RS232/PLC)	100-90594 or 100-90595	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.7 (ATEX) Pump Head LPP-600.18 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.5 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.1-02 (48V/600W, Integrated Display, Enable Connector)	100-90596	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.8 (ATEX, USB/PLC) or PuraLev-600MU.9 (ATEX, RS232/PLC) Pump Head LPP-600.18 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.5 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.2-02 (48V/600W, USB/PLC) or LPC-600.3-02 (RS232/PLC)	100-90597 or 100-90598	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.10 (ATEX, Standalone) Pump Head LPP-600.18 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.4 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.1-02 (48V/600W, Integrated Display, Enable Connector)	100-90599	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.11 (ATEX, USB/PLC) or PuraLev-600MU.12 (ATEX, RS232/PLC) Pump Head LPP-600.18 (PVDF-Housing, Triclamp 1", PFA-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.4 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.2-02 (48V/600W, USB/PLC) or LPC-600.3-02 (RS232/PLC)	100-90632 or 100-90633	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.14 (ATEX, Stainless Steel Pump Head, Standalone) Pump Head LPP-600.27 (SS-Housing, Triclamp 1", SS-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.13 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.1-06 (48V/600W, Standalone Display, PLC)	100-91325	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.15 (ATEX, SS Pump Head, USB/PLC) or PuraLev-600MU.16 (ATEX, RS232/PLC) Pump Head LPP-600.27 (SS-Housing, Triclamp 1", SS-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.13 (3m PVC Cable, Circular Connectors, Epoxy Housing Coating) Controller LPC-600.2-06 (48V/600W, USB/PLC) or LPC-600.3-06 (RS232/PLC)	100-91326 or 100-91327	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.17 (ATEX, Stainless Steel Pump Head, Standalone) Pump Head LPP-600.27 (SS-Housing, Triclamp 1", SS-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.12 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.1-06 (48V/600W, Integrated Display, Enable Connector)	100-91328	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Pump System PuraLev-600MU.18 (ATEX, SS Pump Head, USB/PLC) or PuraLev-600MU.19 (ATEX, RS232/PLC) Pump Head LPP-600.27 (SS-Housing, Triclamp 1", SS-Impeller, EPDM O-Ring, No Drain Port) Motor LPM-600.12 (3m FEP Cable, Circular Connectors, ETFE Housing Coating) Controller LPC-600.2-06 (48V/600W, USB/PLC) or LPC-600.3-06 (RS232/PLC)	100-91329 or 100-91330	PuraLev MU	Extension cables MCAP-600.x.yy and MCAS-600.x.yy for motor to be ordered as separate article (see accessories).
Accessories			
Power Supply TSP 600-148-M (Traco, 48V DC, M = Modified Lev. Design)	100-40013	Accessories Pump Systems	
Air Cooling Module ACM-600.2 (PP, NPT 1/4" Fitting)	190-10140	Accessories Pump Systems	
Fan Cooling Module FCM-600.1 (PP+Talk Cover, 24V/3.4W IP68 Fan, 3m PP Cable/M12 Connector)	190-10401	Accessories Pump Systems	Accessories Pump Systems
Fan Cooling Module Cable, FCC-1.1-1.50 (5.0m, PP Cable, M12 Connector to Open Wire)	190-10407	Accessories Pump Systems	Accessories Pump Systems
Fan Cooling Module Cable, FCC-1.1-1.100 (10.0m)	190-10408	Accessories Pump Systems	Accessories Pump Systems
ATEX Cable Sealing System ACS-A.1	100-90292	Accessories Pump Systems	
Autoclaving Reinforcing Tool ART-600.1 (Anodized Aluminum, for Reinforcement of Pumphead during Autoclaving)	190-10281	Accessories Pump Systems	
User Panel LUH-A.1 (RS232 interface)	100-30300	Accessories Pump Systems	Used with LPC-600.3-02
Adaptor/Extension Power Cable, MCAP-600.2-05 (0.5m) (PVC-Jacket, Circular Hum. Wall-Mount. M23/COMBICON)	190-10227	Accessories Pump Systems	
Adaptor/Extension Power Cable, MCAP-600.2-30 (3.0m)	190-10240		
Adaptor/Extension Power Cable, MCAP-600.2-50 (5.0m)	190-10126		
Adaptor/Extension Power Cable, MCAP-600.2-70 (7.0m)	190-10106		
Adaptor/Extension Power Cable, MCAP-600.2-100 (10.0m)	190-10241		
Adaptor/Extension Sensor Cable, MCAS-600.2-05 (0.5m) (PVC-Jacket, Circular Hum. Wall-Mount. M23/D-SUB)	190-10226	Accessories Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.2-30 (3.0m)	190-10238		
Adaptor/Extension Sensor Cable, MCAS-600.2-50 (5.0m)	190-10127		
Adaptor/Extension Sensor Cable, MCAS-600.2-70 (7.0m)	190-10105		
Adaptor/Extension Sensor Cable, MCAS-600.2-100 (10.0m)	190-10239		
Adaptor/Extension Power Cable, MCAP-600.3-05 (0.5m) (PVC-Jacket, Circular Hummel/COMBICON)	190-10154	Accessories Pump Systems	
Adaptor/Extension Power Cable, MCAP-600.3-30 (3.0m)	190-10155		
Adaptor/Extension Power Cable, MCAP-600.3-50 (5.0m)	190-10129		
Adaptor/Extension Power Cable, MCAP-600.3-70 (7.0m)	190-10156		
Adaptor/Extension Power Cable, MCAP-600.3-100 (10m)	190-10157		
Adaptor/Extension Sensor Cable, MCAS-600.3-05 (0.5m) (PVC-Jacket, Circular Hummel/D-SUB)	190-10158	Accessories Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.3-30 (3.0m)	190-10159		
Adaptor/Extension Sensor Cable, MCAS-600.3-50 (5.0m)	190-10130		
Adaptor/Extension Sensor Cable, MCAS-600.3-70 (7.0m)	190-10160		
Adaptor/Extension Sensor Cable, MCAS-600.3-100 (10m)	190-10161		
Spare Parts			
Motor LPM-600.4 (ATEX, ETFE Coating, 3m FEP Cable with Hummel Connectors)	100-10038	Spare Parts Pump Systems	
Motor LPM-600.5 (ATEX, Epoxy Coating, 3m PVC Cable with Hummel Connectors)	100-10039		
Motor LPM-600.12 (ATEX, ETFE Coating, 3m FEP Cable with Hummel Connectors)	100-10149		
Motor LPM-600.13 (ATEX, Epoxy Coating, 3m PVC Cable with Hummel Connectors)	100-10150		To be used with pump head LPP-600.27 and controllers LPC-600.x-06
Pump Head LPP-600.7 (PVDF Housing, PFA-Impeller, 1" Triclamp Fitting, EPDM O-Ring, With Drain Port)	100-90287	Spare Parts Pump Systems	
Pump Head LPP-600.18 (PVDF Housing, PFA-Impeller, 1" Triclamp Fitting, EPDM O-Ring, No Drain Port)	100-90548		
Pump Head LPP-600.27 (Stainless Steel Housing, SS-Impeller, 1" Triclamp Fitting, EPDM O-Ring, No Drain Port)	100-91319		
Controller LPC-600.1-02 (48V/600W, Standalone Display, PLC, Power Supply Cable/Enable connector incl.)	100-30033	Spare Parts Pump Systems	
Controller LPC-600.2-02 (48V/600W, USB/PLC, Power Supply Cable and PLC Connector incl.)	100-30034		
Controller LPC-600.3-02 (48V/600W, RS232/PLC, Power Supply Cable and PLC Connector incl.)	100-30035		
Controller LPC-600.1-06 (48V/600W, Standalone Display, PLC, Power Supply Cable/Enable connector incl.)	100-30086	Spare Parts Pump Systems	
Controller LPC-600.2-06 (48V/600W, USB/PLC, Power Supply Cable and PLC Connector incl.)	100-30087		
Controller LPC-600.3-06 (48V/600W, RS232/PLC, Power Supply Cable and PLC Connector incl.)	100-30088		