

**Levitronix® Standard Article List and Order Information**  
**PuraLev® Consoles, System Packages, Components and Accessories**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard System Packages with PuraLev® Consoles</b>			
<b>PuraLev® Console System PCS-200SU.1</b> Console LCO-600.1-01 Motor BSM-1.6 (ETFE, 2m FEP, Circular Connector) and Pump Head Socket PHS-200.1 (Anodized Alu) <i>Note: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i>	100-91002	PuraLev Console	
<b>PuraLev® Console System PCS-600SU.2</b> Console LCO-600.1-01 Motor LPM-600.5 (Epoxy, 3m PVC, Circular Connector) and Pump Head Socket PHS-600.1 (Anodized Alu) <i>Note: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i>	100-91003	PuraLev Console	
<b>PuraLev® Console System PCS-200SU.6</b> Console LCO-600.1-01 Motor BSM-1.6 (ETFE, 2m FEP, Circular Connector) and Pump Head Socket PHS-200.1 (Anodized Alu) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min, 4 Set Calibration: Silicone/C-Flex, 20/37°C Water) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91138	PuraLev Console	Default activated calibration set is Silicone 37°C.
<b>PuraLev® Console System PCS-600SU.7</b> Console LCO-600.1-01 Motor LPM-600.5 (Epoxy, 3m PVC, Circular Connector) and Pump Head Socket PHS-600.1 (Anodized Alu) LeviFlow Clamp-On Sensor LFSC-22D-005 (80 l/min, 4 Set Calibration: Silicone/C-Flex, 20/37°C Water) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91142	PuraLev Console	Default activated calibration set is Silicone 37°C.
<b>PuraLev® Console System PCS-200MU.1 (No Drain Port) or PCS-200MU.2 (With drain port)</b> Console LCO-600.1-01 Motor BSM-1.6 (ETFE, 2m FEP, Circular Connector) Multi-Use Pump Head LPP-200.16 (No drain port) or LPP-200.17 (With drain port) <i>Note: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i>	100-91036 or 100-91006	PuraLev Console	
<b>PuraLev® Console System PCS-600MU.1 (No Drain Port) or PCS-600MU.2 (With drain port)</b> Console LCO-600.1-01 Motor LPM-600.5 (Epoxy, 3m PVC, Circular Connector) Multi-Use Pump Head LPP-600.18 (No drain port) or LPP-600.7 (With drain port) <i>Note: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i>	100-91139 or 100-91007	PuraLev Console	
<b>PuraLev® Console System PCS-200MU.11 (No drain port) PCS-200MU.12 (with drain port)</b> Console LCO-600.1-01 Motor BSM-1.6 (ETFE, 2m FEP, Circular Connector) Multi-Use Pump Head LPP-200.16 (No drain port) or LPP-200.17 (With drain port) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min, 4 Calibration Sets: Silicone/C-Flex, 20/37°C Water) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91136 or 100-91137	PuraLev Console	Default activated calibration set is Silicone 37°C.
<b>PuraLev® Console System PCS-600MU.11 (No drain port) or PCS-600MU.12 (With drain port)</b> Console LCO-600.1-01 Motor LPM-600.5 (Epoxy, 3m PVC, Circular Connector) Multi-Use Pump Head LPP-600.18 (No drain port) or LPP-600.7 (With drain port) LeviFlow Clamp-On Sensor LFSC-22D-005 (80 l/min, 4 Calibration Sets: Silicone/C-Flex, 20/37°C Water) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91140 or 100-91141	PuraLev Console	Default activated calibration set is Silicone 37°C.
<b>Note: PuraLev Console System Packages with LFSC-IX flowmeters on request. Adaptor cable ICS-3.1 is needed to connect LFSC-IX clamp-on flowmeters to console.</b>			
<b>PuraLev® Consoles</b>			
<b>PuraLev® Console LCO-600.1-01</b> (200-240 VAC, 600W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i>	100-30500	PuraLev Console	
<b>Standard Pump Motors and Pump Heads Compatible with Console</b>			
<b>Multi-Use Pump Head LPP-600.18 (Triclamp 1", PVDF/PFA, No drain port) or LPP-600.7 (With Drain Port)</b>	100-90548 or 100-90287	Spare Parts Pump System	PuraLev® 600SU/MU pump system components.
<b>Single-Use Pump Head DCP-600.2 (Triclamp 1", PP, Non-Sterile)</b> <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-90784	PuraLev SU	
<b>Single-Use Pump Head DCP-600.2-G25 (Gamma Dosage &gt; 25 kGy)</b> <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91079		
<b>Single-Use Pump Head DCP-600.2-SF1-G25 (Sterile Fittings)</b> <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91294		
<b>Irradiated Single-Use Pump head DCP-600.2-G25 (Dosage ≥ 25 kGy)</b> <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91079	PuraLev SU	
<b>Motor LPM-600.5 (ATEX, Epoxy Coating, 3m PVC Cable with Hummel Connectors)</b>	100-10039	Spare Parts Pump System	
<b>Pumphead Socket PHS-600.1 (Anodized Aluminium, Bayonet Socket Type)</b>	100-90696	Spare Parts Pump System	
<b>Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port)</b>	100-90863 or 100-90864	Spare Parts Pump System	PuraLev® 200SU/MU pump system components.
<b>Single-Use Pump head DCP-200.2 (Triclamp 1/2", PP, Non-Sterile)</b> <i>Note: Needs Pump Head Socket PHS-200.1</i>	100-90734	PuraLev SU	
<b>Single-Use Pump head DCP-200.2-G25 (Dosage ≥ 25 kGy)</b> <i>Note: Needs Pump Head Socket PHS-200.1</i>	100-91078		
<b>Single-Use Pump head DCP-200.3 (Barb 1/2", PP, Non-Sterile)</b> <i>Note: Needs Pump Head Socket PHS-200.1</i>	100-90792		
<b>Single-Use Pump head DCP-200.3-G25 (Dosage ≥ 25 kGy)</b> <i>Note: Needs Pump Head Socket PHS-200.1</i>	100-91122	PuraLev SU	
<b>Single-Use Pump head DCP-200.3-SF1-G25 (Sterile Fittings)</b> <i>Note: Needs Pump Head Socket PHS-200.1</i>	100-91233		
<b>Motor BSM-1.6 (ATEX) (ETFE Coating, 2m FEP Cable with Circular Connector, ATEX)</b>	100-10063	Spare Parts Pump System	
<b>Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Socket Type)</b>	100-90695	Spare Parts Pump System	
<b>Standard LEVIFLOW® Flow Sensors Compatible with Console</b>			
<b>Single-Use Sensor LFS-03SU-Z (Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)</b>	100-30375	LEVIFLOW LFS-SU Flowmeter	For other configurations like gamma irradiated flow sensors consult the LFS-SU product literature.
<b>Single-Use Sensor LFS-03SU-Z-SC1 (Extended Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)</b>	100-30418		
<b>Single-Use Sensor LFS-06SU-Z (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)</b>	100-30377		
<b>Single-Use Sensor LFS-06SU-Z-SC1 (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)</b>	100-30406		
<b>Single-Use Sensor LFS-10SU-Z (Water 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)</b>	100-30397		
<b>Single-Use Sensor LFS-10SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)</b>	100-30408		
<b>Single-Use Sensor LFS-15SU-Z (Water 20°C, 50 l/min, Triclamp 1/2" Fittings, Z-Shape)</b>	100-30412		
<b>Single-Use Sensor LFS-15SU-Z-SC1 (Extended Calibr at 20°C, 50 l/min, Triclamp 1" Fittings, Z-Shape)</b>	100-30431		
<b>Single-Use Sensor LFS-20SU.1-Z (Water 20°C, 80 l/min, Triclamp 1/2" Fittings, Z-Shape)</b>	100-30483		
<b>Single-Use Sensor LFS-20SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape)</b>	100-30464		
<b>LeviFlow Clamp-On Sensor LFSC-05D-001 (For Silicone/C-Flex 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min)</b>	100-30407	Spare Parts LEVIFLOW	For other configurations and sizes consult the LFSC-D product literature.
<b>LeviFlow Clamp-On Sensor LFSC-08D-001 (For Silicone/C-Flex 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min)</b>	100-30396		
<b>LeviFlow Clamp-On Sensor LFSC-12D-007 (For Silicone/C-Flex 3/8"x9/16", Water at 20°C/37°C, 0 - 20 l/min)</b>	100-30390		
<b>LeviFlow Clamp-On Sensor LFSC-22D-005 (For Silicone/C-Flex 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)</b>	100-30391		
<b>LeviFlow Clamp-On Sensor LFSC-30D-001 (For Silicone/C-Flex 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min)</b>	100-30414		
<b>LeviFlow Clamp-On Sensor LFSC-i06X-002 (For Silicone/C-Flex/TPE 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min)</b>	100-30485	LEVIFLOW LFSC Flowmeter	For other configurations and sizes consult the LFSC-IX product literature.
<b>LeviFlow Clamp-On Sensor LFSC-i10X-001 (For Silicone/C-Flex/TPE 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min)</b>	100-30477		
<b>LeviFlow Clamp-On Sensor LFSC-i16X-001 (For Silicone/C-Flex/TPE 3/8"x5/8", Water at 20°C/37°C, 0 - 20 l/min)</b>	100-30482		
<b>LeviFlow Clamp-On Sensor LFSC-i19X-001 (For Silicone/C-Flex/TPE 1/2"x3/4", Water at 20°C/37°C, 0 - 50 l/min)</b>	100-30479		
<b>LeviFlow Clamp-On Sensor LFSC-i25X-001 (For Silicone/C-Flex/TPE 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)</b>	100-30480		
<b>LeviFlow Clamp-On Sensor LFSC-i35X-001 (For Silicone/C-Flex/TPE 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min)</b>	100-30481		
<b>LeviFlow Clamp-On Converter LFSC-i05X.1-001 (For LFP-05SU/-05SU.1 1/8", 0-0.2/0-1 l/min, 4 Calibration Sets)</b>	100-30490	LEVIFLOW LFSC Flowmeter	Calibration for 20 and 37°C water.  For flow pipe spec./order info and for other configurations and sizes consult the LFSC-IX.1 product literature.
<b>LeviFlow Clamp-On Converter LFSC-i09X.1-001 (For LFP-09SU 1/4", 0-4 l/min, 2 Calibration Sets)</b>	100-30491		
<b>LeviFlow Clamp-On Converter LFSC-i13X.1-001 (For LFP-13SU 3/8", 0-20 l/min, 2 Calibration Sets)</b>	100-30492		
<b>LeviFlow Clamp-On Converter LFSC-i16X.1-001 (For LFP-16SU 1/2", 0-50 l/min, 2 Calibration Sets)</b>	100-30498		
<b>LeviFlow Clamp-On Converter LFSC-i20X.1-001 (For LFP-20SU 3/4", 0-80 l/min, 2 Calibration Sets)</b>	100-30499		
<b>LeviFlow Clamp-On Converter LFSC-i23X.1-001 (For LFP-23SU 7/8" and 1", 0-160 l/min, 4 Calibration Sets)</b>	100-30509		
<b>LeviFlow Clamp-On Converter LFSC-i36X.1-001 (For LFP-36SU 13/8", 0-400 l/min, 2 Calibration Sets)</b>	100-30487		
<b>LeviFlow Clamp-On Converter LFSC-i36X.1-001 (For LFP-36SU 13/8", 0-400 l/min, 2 Calibration Sets)</b>	100-30487		
<b>Standard Pressure and Temperature Sensors Compatible with Console</b>			
<b>Pendotech Single-Use Pressure Sensor PREPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)</b>	190-10375	Accessories Console	
<b>Pendotech Single-Use Pressure Sensor PREPS-N-5-5 (Polysulfone, 3/4" Sanitary Flange Fitting, Non-Sterile)</b>	190-10373		
<b>Pendotech Single-Use Pressure Sensor PREPS-N-1-1 (Polysulfone, 1" Sanitary Flange Fitting, Non-Sterile)</b>	190-10374		
<b>Pendotech Single-Use Pressure Sensor PREPS-N-15-15 (Polysulfone, 1.5" Sanitary Flange Fitting, Non-Sterile)</b>	190-10478		
<b>Pendotech Single-Use Temperature Sensor TEMPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)</b>	190-10377		
<b>Pendotech Single-Use Temperature Sensor TEMPS-N-050 (Polysulfone, 1/2" Sanitary Flange Fitting, Non-Sterile)</b>	190-10376		

Article Name (and Description)	Part No.	Article Category	Note
<b>Cables and other Accessories</b>			
LeviFlow Interconnect Cable LFI-C.1-30 (3m, PVC)	190-10308	Accessories LEVIFLOW	For connecting LFSC-D and LFS-SU flow sensors to the console.
IP Cable Signal 4 Wires ICS-3.1-30 (3m, PUR Jacket, Circular to Circular Connectors)	190-10447	Accessories Console	For connection of LFSC-iX sensors to console.
Pressure Sensor Adaptor Cable LPE-1.1-30 (3m, Circular Console to Pressure Sensor Connector)	190-10354	Accessories Console	For pressure sensors.
Temperature Sensor Adaptor Cable LTE-1.1-30 (3m, Circular Console to Temperature Sensor Connector)	190-10353		
IPS Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connection of devices to PLC of console.
Adaptor/Extension Cable MCI-1.3-40 (4m, PVC-Jacket, Circular to Circular)	190-10360	Accessories Pump Systems	For extension of BSM-1 motor cable.
Adaptor/Extension Cable Power / Sensor MCIP-600.1-30/MCIS-600.1-30 (3m, PVC, Circular to Circular)	190-10361 / 62	Accessories Pump Systems	For extension of LPM-600 motor cables.
AC Mains Cable AMC-2.1 (3 m, US, Canada)	190-10336	Accessories Console	Mains cables for connecting console to AC supply. Other countries on request.
AC Mains Cable AMC-2.2 (3 m, Germany, Norway, Finland, Belgium, Netherland, Sweden)	190-10337		
AC Mains Cable AMC-2.3 (3 m, Japan)	190-10338		
AC Mains Cable AMC-2.4 (3 m, Switzerland)	190-10339		
AC Mains Cable AMC-2.5 (3 m, United Kingdom)	190-10340		
Fan Cooling Module FCM-600.1 (PP+Talk Cover, 24V/3.4W IP68 Fan, 3m PP Cable/M12 Connector)	190-10401	Accessories Pump Systems	
FCM Interconnect Cable FCC-1.2-30 (3m, PP jacket, Circular Connectors)	190-10414	Accessories Pump Systems	Other standard length (5, 7m) on request.