

**Levitronix® Standard Article List and Order Information**  
**Maglev Pumpsystem PuraLev® i100MU, Accessories and Spare Parts**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard Pump Systems</b>			
<b>Pump System PuraLev-i100MU.1</b> (OEM Model, No Drain Port) or <b>PuraLev-i100MU.4</b> (Drain Port) Integrated Pump Driver IPD-100.1-50-03 (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating) Pump Head LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, No Drain Port) or Pump Head LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, With Drain Port)	<b>100-91227</b> or <b>100-91228</b>	PuraLev MU	
<b>Pump System PuraLev-i100MU.9</b> (OEM Model ATEX, No Drain Port) or <b>PuraLev-i100MU.8</b> (Drain Port) Integrated Pump Driver IPD-100.10-50-03 (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating) Pump Head LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, No Drain Port) or Pump Head LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, With Drain Port)	<b>100-91513</b> or <b>100-91512</b>	PuraLev MU	
<b>Pump System PuraLev-i100MU.2</b> (EasyConnect Model, No Drain Port) or <b>PuraLev-i100MU.5</b> (Drain Port) Integrated Pump Driver IPD-100.4-03 (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus) Pump Head LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, No Drain Port) or Pump Head LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, With Drain Port)	<b>100-91229</b> or <b>100-91230</b>	PuraLev MU	
<b>Pump System PuraLev-i100MU.3</b> (Stand-Alone Model, No Drain Port) or <b>PuraLev-i100MU.6</b> (Drain Port) Integrated Pump Driver IPD-100.6-03 (24VDC, Epoxy Housing Coating, Display with Buttons) Pump Head LPP-200.16 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, No Drain Port) or Pump Head LPP-200.17 (PVDF Housing, Triclamp 1/2", PFA-Impeller, EPDM O-Ring, With Drain Port)	<b>100-91231</b> or <b>100-91232</b>	PuraLev MU	
<b>Accessories</b>			
<b>Autoclaving Reinforcing Tool ART-200.1</b> (Anodized Alu., reinforcement of pump head during Autoclaving)	<b>190-10280</b>	Accessories Pump Systems	
<b>Power Supply TSP 180-124</b> (Traco, 24V / 180W, DIN Rail Mounting, 85-132 and 187-264VAC autoselec)	<b>100-40018</b>	Accessories Pump Systems	
<b>Desktop AC/DC Power Supply GST160A24-R7B IC915</b> (24VDC/160W, 85-264VAC, Intercontec Conn.)	<b>100-40020</b>	Accessories IPS	AC mains cable AMC-1.x for various countries to be ordered as separate article.
<b>AC Mains Cable AMC-1.1</b> (2m, US, Canada) <b>AC Mains Cable AMC-1.2</b> (2.5m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria) <b>AC Mains Cable AMC-1.3</b> (2.5m, Japan) <b>AC Mains Cable AMC-1.4</b> (2.5m, Switzerland) <b>AC Mains Cable AMC-1.5</b> (2.5m, United Kingdom)	<b>190-10331</b> <b>190-10332</b> <b>190-10333</b> <b>190-10334</b> <b>190-10335</b>	Accessories IPS	AC cable for desktop power supply.
<b>USB-RS485 YN-485I-TR</b> (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	<b>100-30392</b>	Accessories IPS	
<b>IPS Cable Power 2 Wires ICP-2.1-50</b> (5m PVC, Circular Driver Connector to Open Wires)	<b>190-10370</b>	Accessories IPS	For power supplies with open wire connection.
<b>IPS Cable Signal 6 Wires ICS-1.1-02</b> (0.2m PVC, Circular Driver Connector to Circular Connector) <b>IPS Cable Signal 6 Wires ICS-1.1-10</b> (1m) <b>IPS Cable Signal 6 Wires ICS-1.1-30</b> (3m)	<b>190-10381</b> <b>190-10344</b> <b>190-10345</b>	Accessories IPS	For connecting fieldbus of EasyConnect driver arrays.
<b>IPS Cable Signal 6 Wires ICS-1.2-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	<b>190-10346</b>	Accessories IPS	For connecting devices to PLC of Stand-Alone driver or fieldbus of Connectivity driver.
<b>IPS Cable Signal 12 Wires ICS-2.1-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	<b>190-10347</b>	Accessories IPS	For connecting devices to PLC of EasyConnect driver.
<b>Fieldbus Termination Connector FTC-1.1</b> (PPS with FPM Sealing)	<b>190-10348</b>	Accessories IPS	Termination of fieldbus of EasyConnect drivers.
<b>ATEX Cable Sealing System ACS-A.1</b>	<b>100-90292</b>	Accessories IPS	
<b>Spare Parts</b>			
<b>Pump Head LPP-200.16</b> (PVDF Housing, Triclamp 1/2", PFA-Impeller LPI-200.5, EPDM O-Ring, No Drain Port)	<b>100-90863</b>	Spare Parts Pump Systems	With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.
<b>Pump Head LPP-200.17</b> (PVDF Housing, Triclamp 1/2", PFA-Impeller LPI-200.5, EPDM O-Ring, With Drain Port)	<b>100-90864</b>	Spare Parts Pump Systems	With gas bubble robust and low shear force impeller. For liquid viscosities < 20cP.
<b>Integrated Pump Driver IPD-100.1-50-03</b> (OEM Model) (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating)	<b>100-10134</b>	Spare Parts IPS	
<b>Integrated Pump Driver IPD-100.10-50-03</b> (OEM Model, ATEX/IECEX) (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating)	<b>100-10202</b>	Spare Parts IPS	
<b>Integrated Pump Driver IPD-100.4-03</b> (EasyConnect Model) (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus)	<b>100-10135</b>	Spare Parts IPS	
<b>Integrated Pump Driver IPD-100.6-03</b> (Stand-Alone Model) (24VDC, Epoxy Coating, Display with Buttons, PLC)	<b>100-10136</b>	Spare Parts IPS	