

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® iF30SU, Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note
Standard Driver System Configurations			
Flow Control Driver System PLD-iF30SU.1 (OEM Model) Integrated Flow Control Driver IFD-30.3-02 (24VDC, Hybrid Connector, Epoxy Housing Coating) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91223	PuraLev SU	Note: Mounting Base Plate MBP-i30.1 has to be ordered separately.
Flow Control Driver System PLD-iF30SU.2 (EasyConnect Model) Integrated Flow Control Driver IFD-30.4-02 (24VDC, Epoxy Coating, PLC and 2x Fieldbus, incl. MBP-i30.1) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91224	PuraLev SU	
Flow Control Driver System PLD-iF30SU.3 (Stand-Alone Model) Integrated Flow Control Driver IFD-30.5-02 (24VDC, Epoxy Coating, Display/Buttons, PLC, incl. MBP-i30.1) Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-91225	PuraLev SU	
Single-Use Pump Heads			
Single-Use Pumphead DCP-30.2 (Barb 3/8", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile) Single-Use Pumphead DCP-30.1 (Triclamp 3/8", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile)	100-90968 100-90959	PuraLev SU	
Single-Use Pump Head Gamma Irradiated DCP-30.2-G25 (Applied Gamma Dosage ≥ 25 kGy) Single-Use Pump Head Gamma Irradiated DCP-30.1-G25 (Applied Gamma Dosage ≥ 25 kGy)	100-91071 100-91170	PuraLev SU	
Single-Use Pump Head Sterile Fittings DCP-30.2-SF1-G25 (Sterile Fittings, Gamma Dosage ≥ 25 kGy)	100-91234	PuraLev SU	
Standard LEVIFLOW® Flow Sensors Compatible with Flowcontrol Drivers			
Single-Use Sensor LFS-03-SU-Z-F1 (Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape) Single-Use Sensor LFS-03-SU-Z-F1-SC1 (Extended Calibr. at 20°C, 0.8 l/min, Triclamp 3/8", Z-Shape) Single-Use Sensor LFS-06-SU-Z-F1 (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape) Single-Use Sensor LFS-06-SU-Z-F1-SC1 (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8", Z-Shape)	100-30423 100-30470 100-30425 100-30445	LEVIFLOW LFS-SU Flowmeter	
Single-Use Sensor Gamma Irradiated LFS-03SU-Z-F1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor Gamma Irradiated LFS-03SU-Z-F1-SC1-G25 (Irrad. with Gamma Dosage ≥ 25 kGy) Single-Use Sensor Gamma Irradiated LFS-06SU-Z-F1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor Gamma Irradiated LFS-06SU-Z-F1-SC1-G25 (Irrad. with Gamma Dosage ≥ 25 kGy)	100-30424 100-30471 100-30426 100-30444	LEVIFLOW LFS-SU Flowmeter	
Accessories			
Mounting Base Plate MBP-i30.1 (PP+GF)	190-10313	Accessories IPS	
Power Supply TPC 055-124 (Traco, 24V / 55W, DIN Rail Mounting)	100-40014	Accessories IPS	
Desktop AC/DC Power Supply VEC50US24 HR30 (24VDC/50W, 90-264VAC, Hirose Connector)	100-40015	Accessories IPS	AC mains cable ACM-1.x for various countries to be ordered as separate article.
AC Mains Cable AMC-1.1 (2m, US, Canada) AC Mains Cable AMC-1.2 (2.5m, DE, DK, NO, FI, BE, NL, SE, AT) AC Mains Cable AMC-1.3 (2.5m, Japan) AC Mains Cable AMC-1.4 (2.5m, Switzerland) AC Mains Cable AMC-1.5 (2.5m, United Kingdom)	190-10331 190-10332 190-10333 190-10334 190-10335	Accessories IPS	AC cable for desktop power supply VEC50US24.
USB-RS485 YN-485I-TR (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	100-30392	Accessories IPS	
User Panel LUI-B.1-01 (User panel with display and menu buttons)	100-30448	Accessories Pump Systems	Connection to EasyConnect drivers with ICS-1.1 cables.
IPS Cable Power 3 Wires ICP-1.1-50 (5m PVC, Circular Driver Connector to Open Wires)	190-10342	Accessories IPS	For power supplies with open wire connection.
IPS Cable Signal 6 Wires ICS-1.1-01 (0.1m PVC, Circular Driver Connector to Circular Connector) IPS Cable Signal 6 Wires ICS-1.1-10 (1m) IPS Cable Signal 6 Wires ICS-1.1-30 (3m)	190-10343 190-10344 190-10345	Accessories IPS	For connecting fieldbus of EasyConnect driver arrays. For connecting LUI-B.1 to EasyConnect driver.
IPS Cable Signal 6 Wires ICS-1.2-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10346	Accessories IPS	For connecting devices to fieldbus of EasyConnect driver
IPS Cable Signal 6 Wires ICS-1.3-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10389	Accessories IPS	For connecting devices to PLC of Stand-Alone driver.
IPS Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connecting devices to PLC of EasyConnect driver.
IPS Cable Hybrid 15 Wires ICH-1.1-30 (3m PVC, Circular Hybrid Driver Connector to Open Wire) IPS Cable Hybrid 15 Wires ICH-1.1-50 (5m)	190-10386 190-10341	Accessories IPS	For connecting supply and signals to OEM driver.
Fieldbus Termination Connector FTC-1.1 (PPS with FPM Sealing)	190-10348	Accessories IPS	Termination of fieldbus of EasyConnect drivers.
LeviFlow Mounting Kit LMK-1.1 (Anodized Alu., for LFS-03SU and LFS-06SU)	100-91179	Accessories LEVIFLOW	For mounting of SU flow sensors.
Spare Parts			
Integrated Pump Driver IFD-30.3-02 (OEM Model) (24VDC, Circular Hybrid Connector, Epoxy Housing Coating)	100-10130	Spare Parts IPS	
Integrated Pump Driver IFD-30.4-02 (EasyConnect Model) (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus, includes Mounting Base Plate MBP-i30.1)	100-10131	Spare Parts IPS	
Integrated Pump Driver IFD-30.5-02 (Stand-Alone Model) (24VDC, Epoxy Coating, Display with Buttons, PLC, includes Mounting Base Plate MBP-i30.1)	100-10132	Spare Parts IPS	
Pump Head Socket PHS-i30.1 (Anodized Aluminum, Bayonet Type)	100-90947	Spare Parts IPS	