

Article Name (and Description)	Part No.	Article Category	Note
Accessories			
Air Cooling Module ACM-4.2 (PP, NPT1/4", Molded)	190-10139	Accessories Pump Systems	Not allowed in HazLoc or ATEX/IECEx applications.
Air Cooling Module ACM-4.3 (PP-EL-S, NPT1/4", Machined, for ATEX Motors)	190-10243	Accessories Pump Systems	
ATEX Cable Sealing System ACS-A.1	100-90292	Accessories Pump Systems	
Fan Cooling Module FCM-2000.1 (PP+Talc Cover, 24V/33.5W IP68 Fan, 3m PP Cable/M12 Connector)	190-10390	Accessories Pump Systems	Not allowed in HazLoc or ATEX/IECEx applications.
Fan Cooling Module Cable, FCC-1.1-50 (5.0m, PP Cable, M12 Connector to Open Wire)	190-10407	Accessories Pump Systems	
Fan Cooling Module Cable, FCC-1.1-100 (10.0m)	190-10408	Accessories Pump Systems	
Autoclaving Reinforcing Tool ART-2000.1 (Anodized Aluminum, Reinforcement of Pumphead during Autoclaving)	190-10282	Accessories Pump Systems	
User Panel LUI-A.1 (RS232 interface)	100-30300	Accessories Pump Systems	Used with LPC-2000.3-05 or LPC-2000.3-13.
Adaptor/Extension Power Cable, MCAP-2000.2-05 (0.5m) (PVC-Jacket, Circ. Hum. Wall-Mount. M23/COMBICON)	190-10213	Accessories Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.2-30 (3.0m)	190-10215		
Adaptor/Extension Power Cable, MCAP-2000.2-50 (5.0m)	190-10216		
Adaptor/Extension Power Cable, MCAP-2000.2-70 (7.0m)	190-10217		
Adaptor/Extension Power Cable, MCAP-2000.2-100 (10.0m)	190-10218		
Adaptor/Extension Sensor Cable, MCAS-600.2-05 (0.5m) (PVC-Jacket, Circ. Hum. Wall-Mountable M23/D-SUB)	190-10226	Accessories Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.2-30 (3.0m)	190-10238		
Adaptor/Extension Sensor Cable, MCAS-600.2-50 (5.0m)	190-10127		
Adaptor/Extension Sensor Cable, MCAS-600.2-70 (7.0m)	190-10105		
Adaptor/Extension Sensor Cable, MCAS-600.2-100 (10.0m)	190-10239		
Adaptor/Extension Power Cable, MCAP-2000.3-05 (0.5m) (PVC-Jacket, Circ. Hummel M23/COMBICON)	190-10219	Accessories Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.3-30 (3.0m)	190-10221		
Adaptor/Extension Power Cable, MCAP-2000.3-50 (5.0m)	190-10222		
Adaptor/Extension Power Cable, MCAP-2000.3-70 (7.0m)	190-10223		
Adaptor/Extension Power Cable, MCAP-2000.3-100 (10.0m)	190-10224		
Adaptor/Extension Sensor Cable, MCAS-600.3-05 (0.5m) (PVC-Jacket, Circular Hummel M23/D-SUB)	190-10158	Accessories Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.3-30 (3.0m)	190-10159		
Adaptor/Extension Sensor Cable, MCAS-600.3-50 (5.0m)	190-10130		
Adaptor/Extension Sensor Cable, MCAS-600.3-70 (7.0m)	190-10160		
Adaptor/Extension Sensor Cable, MCAS-600.3-100 (10.0m)	190-10161		
Interconnect Power Cable, MCIP-2000.1-05 (0.5m) (PVC-Jacket, Circular Hummel/Circular Hummel)	190-10396	Accessories Pump Systems	This cable is required to connect the HazLoc motors in addition to MCAP-2000.2/3-YY.
Interconnect Power Cable, MCIP-2000.1-30 (3.0m)	190-10397		
Interconnect Power Cable, MCIP-2000.1-50 (5.0m)	190-10398		
Interconnect Power Cable, MCAP-2000.1-70 (7.0m)	190-10399		
Interconnect Power Cable, MCAP-2000.1-100 (10.0m)	190-10400		
Interconnect Sensor Cable, MCIS-2000.1-05 (0.5m) (PVC-Jacket, Circular Hummel/Circular Hummel)	190-10391	Accessories Pump Systems	This cable is required to connect the HazLoc motors in addition to MCAS-2000.2/3-YY.
Interconnect Sensor Cable, MCIS-2000.1-30 (3.0m)	190-10392		
Interconnect Sensor Cable, MCIS-2000.1-50 (5.0m)	190-10393		
Interconnect Sensor Cable, MCIS-2000.1-70 (7.0m)	190-10394		
Interconnect Sensor Cable, MCIS-2000.1-100 (10.0m)	190-10395		
Spare Parts			
Motor LPM-2000.7 (ATEX, Epoxy Coating, 3m PVC Cable with Circular Connectors)	100-10059	Spare Parts Pump System	These motors can be used with plastics pump heads only.
Motor LPM-2000.8 (ATEX, ETFE Coating, 3m FEP Cable with Circular Connectors)	100-10060		
Motor LPM-2000.9 (HazLoc Class1 Div2, ETFE Coating, Circular Connectors)	100-10112		
Motor LPM-2000.14 (ATEX, Epoxy Coating, 3m PVC Cable with Circular Connectors)	100-10163	Spare Parts Pump Systems	These motors can be used with stainless steel pump heads only.
Motor LPM-2000.15 (ATEX, ETFE Coating, 3m FEP Cable with Circular Connectors)	100-10164		
Motor LPM-2000.17 (HazLoc) (HazLoc Class1 Div2, ETFE Coating, Circular Connectors)	100-10166		
Pump Head CP-4.11 (PVDF Housing, Triclamp 1" Fitting, PFA-Impeller, EPDM Seal, SS Screws, With Drain Port)	100-90245	Spare Parts Pump Systems	Needs specific controllers and motors.
Pump Head CP-4.19 (PVDF Housing, Triclamp 1" Fitting, PFA-Impeller, EPDM Seal, SS Screws, No Drain Port)	100-90458	Spare Parts Pump Systems	Needs specific controllers and motors.
Pump Head LPP-2000.25 (Stainl. Steel Hous./Imp., Triclamp 1" Fitting, EPDM O-Ring, SS Screws, With Drain Port)	100-91545	Spare Parts Pump Systems	Needs specific controllers and motors.
Pump Head LPP-2000.19 (Stainl. Steel Hous./Imp., Triclamp 1" Fitting, EPDM O-Ring, SS Screws, No Drain Port)	100-91446	Spare Parts Pump Systems	Needs specific controllers and motors.
Controller LPC-2000.1-05 (Stand-alone, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30036	Spare Parts Pump Systems	These controllers can be used with plastics pump heads only.
Controller LPC-2000.2-05 (USB and PLC Interface, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30037		
Controller LPC-2000.3-05 (RS232+PLC Interface, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30038		
Controller LPC-2000.1-13 (Stand-alone, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30096	Spare Parts Pump Systems	These controllers can be used with stainless steel pump heads only.
Controller LPC-2000.2-13 (USB and PLC Interface, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30095		
Controller LPC-2000.3-13 (RS232+ PLC Interface, 2kW, High Flow, Power Supply and PLC Connectors included)	100-30097		