

**Levitronix® Standard Article List and Order Information**  
**Maglev Pumpsystem PuraLev® iF100SU, Accessories and Spare Parts**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard Driver System Configurations</b>			
<b>Flow Control Driver System PLD-iF100SU.1 (OEM Model)</b> Integrated Flow Control Driver IFD-100.4-02 (24VDC, Hybrid Connector, Epoxy Housing Coating) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91482	PuraLev SU	
<b>Flow Control Driver System PLD-iF100SU.2 (EasyConnect Model)</b> Integrated Flow Control Driver IFD-100.5-02 (24VDC, Epoxy Coating, PLC and 2x Fieldbus) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91483	PuraLev SU	
<b>Flow Control Driver System PLD-iF100SU.3 (Stand-Alone Model)</b> Integrated Flow Control Driver IFD-100.6-02 (24VDC, Epoxy Coating, Display/Buttons, PLC) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91484	PuraLev SU	
<b>Single-Use Pump Heads</b>			
<b>Single-Use Pumphead DCP-200.3</b> (Barb 1/2", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile) <b>Single-Use Pumphead DCP-200.2</b> (Triclamp 1/2", PP FDA, USP-VI, BSE/TSE, ADI Free, Non-Sterile)	100-90792 100-90734	PuraLev SU	
<b>Single-Use Pump Head Gamma Irradiated DCP-200.3-G25</b> (Applied Gamma Dosage ≥ 25 kGy) <b>Single-Use Pump Head Gamma Irradiated DCP-200.2-G25</b> (Applied Gamma Dosage ≥ 25 kGy)	100-91122 100-91078	PuraLev SU	
<b>Single-Use Pump Head Sterile Fittings DCP-200.3-SF1-G25</b> (Sterile Fittings, Gamma Dosage ≥ 25 kGy)	100-91233	PuraLev SU	
<b>Standard LEVIFLOW® Flow Sensors Compatible with Flowcontrol Drivers</b>			
<b>Single-Use Sensor LFS-06SU-Z</b> (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape) <b>Single-Use Sensor LFS-06SU-Z-SC1</b> (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8", Z-Shape) <b>Single-Use Sensor LFS-10SU-Z</b> (Calibration at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape) <b>Single-Use Sensor LFS-10SU-Z-SC1</b> (Calibr. at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30377 100-30394 100-30397 100-30408	LEVIFLOW LFS-SU Flowmeter	
<b>Single-Use Sensor Gamma Irradiated LFS-06SU-Z-G25</b> (Irradiated with Gamma Dosage ≥ 25 kGy) <b>Single-Use Sensor Gamma Irradiated LFS-06SU-Z-SC1-G25</b> (Irrad. with Gamma Dosage ≥ 25 kGy) <b>Single-Use Sensor Gamma Irradiated LFS-10SU-Z-G25</b> (Irradiated with Gamma Dosage ≥ 25 kGy) <b>Single-Use Sensor Gamma Irradiated LFS-10SU-Z-SC1-G25</b> (Irrad. with Gamma Dosage ≥ 25 kGy)	100-30400 100-30406 100-30405 100-30416	LEVIFLOW LFS-SU Flowmeter	
<b>Accessories</b>			
<b>Power Supply TSP 180-124</b> (Traco, 24V / 180W, DIN Rail Mounting, 85-132 and 187-264VAC autoselec)	100-40018	Accessories Pump Systems	
<b>Desktop AC/DC Power Supply GST160A24-R7B IC915</b> (24VDC/160W, 85-264VAC, Intercontec Conn.)	100-40020	Accessories IPS	AC mains cable AMC-1.x for various countries to be ordered as separate article.
<b>AC Mains Cable AMC-1.1</b> (2m, US, Canada) <b>AC Mains Cable AMC-1.2</b> (2.5m, DE, DK, NO, FI, BE, NL, SE, AT) <b>AC Mains Cable AMC-1.3</b> (2.5m, Japan) <b>AC Mains Cable AMC-1.4</b> (2.5m, Switzerland) <b>AC Mains Cable AMC-1.5</b> (2.5m, United Kingdom)	190-10331 190-10332 190-10333 190-10334 190-10335	Accessories IPS	AC cable for desktop power supply VEC50US24.
<b>USB-RS485 YN-485I-TR</b> (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	100-30392	Accessories IPS	
<b>User Panel LUI-B.1-01</b> (User panel with display and menu buttons)	100-30448	Accessories Pump Systems	Connection to EasyConnect drivers with ICS-1.1 cables.
<b>IPS Cable Power 3 Wires ICP-1.1-50</b> (5m PVC, Circular Driver Connector to Open Wires)	190-10342	Accessories IPS	For power supplies with open wire connection.
<b>IPS Cable Signal 6 Wires ICS-1.1-01</b> (0.1m PVC, Circular Driver Connector to Circular Connector) <b>IPS Cable Signal 6 Wires ICS-1.1-10</b> (1m) <b>IPS Cable Signal 6 Wires ICS-1.1-30</b> (3m)	190-10343 190-10344 190-10345	Accessories IPS	For connecting fieldbus of EasyConnect driver arrays. For connecting LUI-B.1 to EasyConnect driver.
<b>IPS Cable Signal 6 Wires ICS-1.2-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10346	Accessories IPS	For connecting devices to fieldbus of EasyConnect driver
<b>IPS Cable Signal 6 Wires ICS-1.3-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10389	Accessories IPS	For connecting devices to PLC of Stand-Alone driver.
<b>IPS Cable Signal 12 Wires ICS-2.1-50</b> (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connecting devices to PLC of EasyConnect driver.
<b>IPS Cable Hybrid 15 Wires ICH-1.1-30</b> (3m PVC, Circular Hybrid Driver Connector to Open Wire) <b>IPS Cable Hybrid 15 Wires ICH-1.1-50</b> (5m)	190-10386 190-10341	Accessories IPS	For connecting supply and signals to OEM driver.
<b>Fieldbus Termination Connector FTC-1.1 (PPS with FPM Sealing)</b>	190-10348	Accessories IPS	Termination of fieldbus of EasyConnect drivers.
<b>LeviFlow Mounting Kit LMK-1.2</b> (Anodized Alu., for LFS-03SU and LFS-06SU) <b>LeviFlow Mounting Kit LMK-2.2</b> (Anodized Alu., for LFS-10SU)	100-91478 100-91479	Accessories LEVIFLOW	For mounting of SU flow sensors.
<b>Spare Parts</b>			
<b>Integrated Pump Driver IFD-100.4-02 (OEM Model)</b> (24VDC, Circular Hybrid Connector, Epoxy Housing Coating)	100-10138	Spare Parts IPS	
<b>Integrated Pump Driver IFD-100.5-02 (EasyConnect Model)</b> (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus, includes Mounting Base Plate MBP-i30.1)	100-10139	Spare Parts IPS	
<b>Integrated Pump Driver IFD-100.6-02 (Stand-Alone Model)</b> (24VDC, Epoxy Coating, Display with Buttons, PLC, includes Mounting Base Plate MBP-i30.1)	100-10140	Spare Parts IPS	
<b>Pump Head Socket PHS-i100.1</b> (Anodized Aluminum, Bayonet Type)	100-91053	Spare Parts IPS	