

**Levitronix® Standard Article List and Order Information**  
**PuraLev® Console LCO-i100, System Packages, Components and Accessories**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard System Packages with PuraLev® Consoles</b>			
<b>PuraLev® Console System PCS-i30SU.1</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i30.1 (Anod. Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91098	PuraLev Console	
<b>PuraLev® Console System PCS-i100SU.1</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector) and Pump Head Socket PHS-i100.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91097	PuraLev Console	
<b>PuraLev® Console System PCS-i600SU.1</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC) Driver IPD-600.3-30-02 (Epoxy Coat., 3m PVC Cable, Connector) and Pump Head Socket PHS-600.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Split cable ICY-4.1-02 and ICS-1.5-x to connect power supply and driver to console must be ordered as separate article.</i>	100-91494	PuraLev Console	
<b>PuraLev® Console System PCS-i30SU.5 or PCS-i30SU.3 or PCS-i30SU.4</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i30.1 (Anod. Alu) LexiFlow Clamp-On Sensor LFSC-05D-001 (1 l/min) or LFSC-08D-001 (4 l/min) or LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91349 or 100-91099 or 100-91100	PuraLev Console	
<b>PuraLev® Console System PCS-i100SU.2 or PCS-i100SU.3</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i100.1 (Anod. Alu) LexiFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) or LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91096 or 100-91092	PuraLev Console	
<b>PuraLev® Console System PCS-i30MU.1</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-06 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) Multi-Use Pump Head LPP-30.7 (Triclamp 3/8", PEEK/PFA/EPDM, No drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91562	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.1 or PCS-i100MU.2</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91268 or 100-91269	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.3 or PCS-i100MU.4</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LexiFlow Clamp-On Sensor LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91270 or 100-91271	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.5 or PCS-i100MU.6</b> Console LCO-1100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LexiFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91272 or 100-91273	PuraLev Console	
<i>Note: PuraLev Console System Packages with LFSC-iX flowmeters on request. Adaptor cable ICS-3.1 is needed to connect LFSC-iX clamp-on flowmeters to console.</i>			
<b>PuraLev® Consoles</b>			
<b>PuraLev® Console LCO-1100.1-01</b> (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i>	100-30501	PuraLev Console	
<b>PuraLev® Console LCO-1100.1-03 - FDA CFR 21 Part 11 Ready</b> (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i>	100-30539	PuraLev Console	
<b>Standard Pump Drivers and Pump Heads Compatible with Console</b>			
<b>Single-Use Pump Head DCP-30.1</b> (Triclamp 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-90959	PuraLev SU	PuraLev® i30SU/MU pump system components.  For other configurations consult the PuraLev® i30SU/MU product literature.
<b>Single-Use Pump Head DCP-30.1-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91170		
<b>Single-Use Pump Head DCP-30.2</b> (Triclamp 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-90968	PuraLev SU	
<b>Single-Use Pump Head DCP-30.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91071		
<b>Single-Use Pump Head DCP-30.2-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91234	PuraLev SU	
<b>Single-Use Pump Head DCP-30.5</b> (Barb 1/4", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91691		
<b>Single-Use Pump Head DCP-30.5-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91700	PuraLev MU	
<b>Multi-Use Pump Head LPP-30.7</b> (Triclamp 3/8", PEEK/PFA/EPDM, No drain port)	100-91557		
<b>Driver IPD-30.8-03-02</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, SU Firmware)	100-10104	Spare Parts Pump System	
<b>Driver IPD-30.8-03-06</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, MU Firmware)	100-10212		
<b>Pumphead Socket PHS-i30.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-90947	Spare Parts Pump System	
<b>Single-Use Pump head DCP-200.2</b> (Triclamp 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-90734	PuraLev SU	
<b>Single-Use Pump Head DCP-200.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91078		
<b>Single-Use Pump head DCP-200.3</b> (Barb 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-90792	PuraLev SU	
<b>Single-Use Pump Head DCP-200.3-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91122		
<b>Single-Use Pump Head DCP-200.3-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91233	PuraLev MU	
<b>Multi-Use Pump Head LPP-200.16</b> (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or <b>LPP-200.17</b> (With drain port)	100-90863 or 100-90864		
<b>Driver IPD-100.3-03-02</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connector, SU Firmware)	100-10105	Spare Parts Pump System	
<b>Driver IPD-100.3-03-03</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connector, MU Firmware)	100-10137		
<b>Pumphead Socket PHS-i100.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-91053	Spare Parts Pump System	
<b>Single-Use Pump Head DCP-600.2</b> (Triclamp 1", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-90784	PuraLev SU	
<b>Single-Use Pump Head DCP-600.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91079		
<b>Single-Use Pump Head DCP-600.2-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91294	PuraLev® i600SU pump system components.	
<b>Driver IPD-600.3-30-02</b> (Epoxy Coating, 3m PVC Cable with Circular Connector, SU Firmware)	100-10197		
<b>Pumphead Socket PHS-600.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-90696	Spare Parts Pump Systems	

Article Name (and Description)	Part No.	Article Category	Note		
<b>Standard LEVIFLOW® Flow Sensors Compatible with Console</b>					
Single-Use Sensor LFS-03SU-Z (Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30375	LEVIFLOW LFS-SU Flowmeter	For other configurations like gamma irradiated flow sensors consult the LFS-SU product literature.		
Single-Use Sensor LFS-03SU-Z-SC1 (Extended Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30418				
Single-Use Sensor LFS-06SU-Z (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30377				
Single-Use Sensor LFS-06SU-Z-SC1 (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30406				
Single-Use Sensor LFS-10SU-Z (Water 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30397				
Single-Use Sensor LFS-10SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30408				
Single-Use Sensor LFS-15SU-Z (Water 20°C, 50 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30412				
Single-Use Sensor LFS-15SU-Z-SC1 (Extended Calibr at 20°C, 50 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30431				
Single-Use Sensor LFS-20SU-1-Z (Water 20°C, 80 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30483				
Single-Use Sensor LFS-20SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape)	100-30464				
LeviFlow Clamp-On Sensor LFSC-05D-001 (For Silicone/C-Flex 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min)	100-30407	Spare Parts LEVIFLOW	For other configurations and sizes consult the LFSC-D product literature.		
LeviFlow Clamp-On Sensor LFSC-08D-001 (For Silicone/C-Flex 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min)	100-30396				
LeviFlow Clamp-On Sensor LFSC-12D-007 (For Silicone/C-Flex 3/8"x9/16", Water at 20°C/37°C, 0 - 20 l/min)	100-30390				
LeviFlow Clamp-On Sensor LFSC-22D-005 (For Silicone/C-Flex 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)	100-30391				
LeviFlow Clamp-On Sensor LFSC-30D-001 (For Silicone/C-Flex 1"x1.318", Water at 20°C/37°C, 0 - 160 l/min)	100-30414				
LeviFlow Clamp-On Sensor LFSC-06X-002 (For Silicone/C-FlexTPE 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min)	100-30485				
LeviFlow Clamp-On Sensor LFSC-10X-001 (For Silicone/C-FlexTPE 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min)	100-30477				
LeviFlow Clamp-On Sensor LFSC-16X-001 (For Silicone/C-FlexTPE 3/8"x5/8", Water at 20°C/37°C, 0 - 20 l/min)	100-30482				
LeviFlow Clamp-On Sensor LFSC-19X-001 (For Silicone/C-FlexTPE 1/2"x3/4", Water at 20°C/37°C, 0 - 50 l/min)	100-30479				
LeviFlow Clamp-On Sensor LFSC-125X-001 (For Silicone/C-FlexTPE 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)	100-30480				
LeviFlow Clamp-On Sensor LFSC-135X-001 (For Silicone/C-FlexTPE 1"x1.318", Water at 20°C/37°C, 0 - 160 l/min)	100-30481	LEVIFLOW LFSC Flowmeter	For other configurations and sizes consult the LFSC-iX product literature.		
LeviFlow Clamp-On Converter LFSC-05X-1-001 (For LFP-05SU-05SU.1 1/8", 0-0.2/0-1 l/min, 2 Calibration Sets)	100-30490				
LeviFlow Clamp-On Converter LFSC-09X-1-001 (For LFP-09SU 1/4", 0-4 l/min, 2 Calibration Sets)	100-30491				
LeviFlow Clamp-On Converter LFSC-13X-1-001 (For LFP-13SU 3/8", 0-20 l/min, 2 Calibration Sets)	100-30492				
LeviFlow Clamp-On Converter LFSC-16X-1-001 (For LFP-16SU 1/2", 0-50 l/min, 2 Calibration Sets)	100-30498				
LeviFlow Clamp-On Converter LFSC-20X-1-001 (For LFP-20SU 3/4", 0-80 l/min, 2 Calibration Sets)	100-30499				
LeviFlow Clamp-On Converter LFSC-23X-1-001 (For LFP-23SU 7/8" and 1", 0-160 l/min, 4 Calibration Sets)	100-30509				
LeviFlow Clamp-On Converter LFSC-36X-1-001 (For LFP-36SU 1.318", 0-400 l/min, 2 Calibration Sets)	100-30487				
<b>Standard Pressure and Temperature Sensors Compatible with Console</b>					
Pendotech Single-Use Pressure Sensor PREPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)	190-10375			Accessories Console	
Pendotech Single-Use Pressure Sensor PREPS-N-5-5 (Polysulfone, 3/4" Sanitary Flange Fitting, Non-Sterile)	190-10373				
Pendotech Single-Use Pressure Sensor PREPS-N-1-1 (Polysulfone, 1" Sanitary Flange Fitting, Non-Sterile)	190-10374				
Pendotech Single-Use Pressure Sensor PREPS-N-15-15 (Polysulfone, 1.5" Sanitary Flange Fitting, Non-Sterile)	190-10478				
Pendotech Single-Use Temperature Sensor TEMPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)	190-10377				
Pendotech Single-Use Temperature Sensor TEMPS-N-050 (Polysulfone, 1/2" Sanitary Flange Fitting, Non-Sterile)	190-10376				
<b>Cables and Desktop AC/DC Power Supplies</b>					
LeviFlow Interconnect Cable LFC-1-1-30 (3m, PVC)	190-10308	Accessories LEVIFLOW	For connecting LEVIFLOW® Sensors to the console.		
IP Cable Signal 4 Wires ICS-3.1-30 (3m, PUR Jacket, Circular to Circular Connectors)	190-10447	Accessories Console	For connection of LFSC-iX sensors to console.		
Pressure Sensor Adaptor Cable LPE-1.1-30 (3m, Circular Console to Pressure Sensor Connector)	190-10354	Accessories Console	For pressure sensors.		
Temperature Sensor Adaptor Cable LTE-1.1-30 (3m, Circular Console to Temperature Sensor Connector)	190-10353	Accessories Console	For temperature sensors.		
IPS Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connection of devices to PLC of console.		
IP Adaptor Cable Hybrid ICH-1.2-30 (3m, PVC-Jacket, Circular to Circular)	190-10329	Accessories IPS	For extension of driver cable.		
IPS Split Cable ICY-4.1-02 (0.2m, PVC cables, PC Split Box, PUR Sealing, circ. Connectors, Integrated 48 to 24VDC converter)	190-10476	Accessories IPS	Splits cable to connect power supply and user panel to driver.		
IPS Cable Signal 4 Wires ICS-1.5-30 (3m PVC, Circular Hirose to Console Connector)	190-10468	Accessories Console	For connecting split box to console power connector.		
Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC/50W, 90-264VAC, Circular Intercontec Connector)	100-40019	Accessories IPS	AC mains cable ACM-1.x for various countries to be ordered as separate article.		
Desktop AC/DC Power Supply GST160A24-R7B IC915 (24VDC/160W, 85-264VAC, Circular Intercontec Connector)	100-40020	Accessories IPS			
AC Mains Cable AMC-1.1 (2m, US, Canada)	190-10331	Accessories Console	Mains cables for connecting console to AC supply. Cables for other countries on request.		
AC Mains Cable AMC-1.2 (2.5m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria)	190-10332				
AC Mains Cable AMC-1.3 (2.5m, Japan)	190-10333				
AC Mains Cable AMC-1.4 (2.5m, Switzerland)	190-10334				
AC Mains Cable AMC-1.5 (2.5m, United Kingdom)	190-10335				
AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC, Circ. Connectors)	100-40026	Accessories IPS	AC mains cable AMC-2.x for various countries to be ordered as separate article.		
AC Mains Cable AMC-2.1 (3m, US, Canada)	190-10336	Accessories IPS	Mains cable for AC/DC power supply DDP-600-US48-SC		
AC Mains Cable AMC-2.2 (3m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria)	190-10337				
AC Mains Cable AMC-2.3 (3m, Japan)	190-10338				
AC Mains Cable AMC-2.4 (3m, Switzerland)	190-10339				
AC Mains Cable AMC-2.5 (3m, United Kingdom)	190-10340				