

| Levitronix® Standard Article List and Order Information PuraLev® Console LCO-i100, System Packages, Components and Accessories | | | |
|---|--|--------------------------|---|
| Article Name (and Description) | Part No. | Article Category | Note |
| Standard System Packages with PuraLev® Consoles | | | |
| PuraLev® Console System PCS-i30SU.1 Console LCO-i100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i301.) and Pump Head Socket PHS-i30.1 (Anod. Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> | 100-91098 | PuraLev Console | |
| PuraLev® Console System PCS-i100SU.1 Console LCO-i100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector) and Pump Head Socket PHS-i100.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> | 100-91097 | PuraLev Console | |
| PuraLev® Console System PCS-i600SU.1 Console LCO-i100.1-01 with AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC) Driver IPD-600.3-30-02 (Epoxy Coat., 3m PVC Cable, Connector) and Pump Head Socket PHS-600.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Split cable ICY-4.1-02 and ICS-1.5-xx to connect power supply and driver to console must be ordered as separate article.</i> | 100-91494 | PuraLev Console | |
| PuraLev® Console System PCS-i30SU.5 or PCS-i30SU.3 or PCS-i30SU.4 Console LCO-i100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i301.) and Pump Head Socket PHS-i30.1 (Anod. Alu) LeviFlow Clamp-On Sensor LFSC-05D-001 (1 l/min) or LFSC-08D-001 (4 l/min) or LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as seperate article.</i> | 100-91349 or 100-91099 or 100-91100 | PuraLev Console | |
| PuraLev® Console System PCS-i100SU.2 or PCS-i100SU.3 Console LCO-i100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i301.) and Pump Head Socket PHS-i100.1 (Anod. Alu) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) or LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as seperate article.</i> | 100-91096 or 100-91092 | PuraLev Console | |
| PuraLev® Console System PCS-i30MU.1 Console LCO-i100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-06 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) Multi-Use Pump Head LPP-30.7 (Triclamp 3/8", PEEK/PTFE/EPDM, No drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> | 100-91562 | PuraLev Console | |
| PuraLev® Console System PCS-i100MU.1 or PCS-i100MU.2 Console LCO-i100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> | 100-91268 or 100-91269 | PuraLev Console | |
| PuraLev® Console System PCS-i100MU.3 or PCS-i100MU.4 Console LCO-i100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LeviFlow Clamp-On Sensor LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as seperate article.</i> | 100-91270 or 100-91271 | PuraLev Console | |
| PuraLev® Console System PCS-i100MU.5 or PCS-i100MU.6 Console LCO-i100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C.1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i> | 100-91272 or 100-91273 | PuraLev Console | |
| Note: PuraLev Console System Packages with LFSC-iX flowmeters on request. Adaptor cable ICS-3.1 is needed to connect LFSC-iX clamp-on flowmeters to console. | | | |
| PuraLev® Consoles | | | |
| PuraLev® Console LCO-i100.1-01 (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i> | 100-30501 | PuraLev Console | |
| PuraLev® Console LCO-i100.1-03 - FDA CFR 21 Part 11 Ready (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i> | 100-30539 | PuraLev Console | |
| Standard Pump Drivers and Pump Heads Compatible with Console | | | |
| Single-Use Pump Head DCP-30.1 (Triclamp 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i> Single-Use Pump Head DCP-30.1-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i> | 100-90959 100-91170 | PuraLev SU | PuraLev® i30SU/MU pump system components. For other configurations consult the PuraLev® i30SU/MU product literature. |
| Single-Use Pump Head DCP-30.2 (Barb 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i> Single-Use Pump Head DCP-30.2-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i> Single-Use Pump Head DCP-30.2-SF1-G25 (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i30.1</i> | 100-90968 100-91071 100-91234 | PuraLev SU | |
| Single-Use Pump Head DCP-30.5 (Barb 1/4", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i> Single-Use Pump Head DCP-30.5-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i> | 100-91691 100-91700 | PuraLev SU | |
| Multi-Use Pump Head LPP-30.7 (Triclamp 3/8", PEEK/PTFE/EPDM, No drain port) | 100-91557 | PuraLev MU | |
| Driver IPD-30.8-03-02 (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, SU Firmware) Driver IPD-30.8-03-06 (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, MU Firmware) | 100-10104 100-10212 | Spare Parts Pump System | |
| Pumphead Socket PHS-i30.1 (Anodized Aluminium, Bayonet Socket Type) | 100-90947 | Spare Parts Pump System | |
| Single-Use Pump head DCP-200.2 (Triclamp 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i> Single-Use Pump Head DCP-200.2-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i> Single-Use Pump head DCP-200.3 (Barb 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i> Single-Use Pump Head DCP-200.3-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i> Single-Use Pump Head DCP-200.3-SF1-G25 (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i100.1</i> | 100-90734 100-91078 100-90792 100-91122 100-91233 | PuraLev SU | PuraLev® i100SU/MU pump system components. For other configurations consult the PuraLev® i100SU/MU product literature. |
| Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) | 100-90863 or 100-90864 | PuraLev MU | |
| Driver IPD-100.3-03-02 (Epoxy Coating, 0.3m PVC Cable with Circular Connector, SU Firmware) Driver IPD-100.3-03-03 (Epoxy Coating, 0.3m PVC Cable with Circular Connector, MU Firmware) | 100-10105 100-10137 | Spare Parts Pump System | |
| Pumphead Socket PHS-i100.1 (Anodized Aluminium, Bayonet Socket Type) | 100-91053 | Spare Parts Pump System | |
| Single-Use Pump head DCP-600.2 (Triclamp 1", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-600.1</i> Single-Use Pump Head DCP-600.2-G25 (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-600.1</i> Single-Use Pump Head DCP-600.2-SF1-G25 (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-600.1</i> | 100-90784 100-91079 100-91294 | PuraLev SU | PuraLev® i600SU pump system components. |
| Driver IPD-600.3-30-02 (Epoxy Coating, 3m PVC Cable with Circular Connector, SU Firmware) | 100-10197 | Spare Parts IPS | |
| Pumphead Socket PHS-600.1 (Anodized Aluminium, Bayonet Socket Type) | 100-90696 | Spare Parts Pump Systems | |

| Article Name (and Description) | Part No. | Article Category | Note |
|---|-----------|---------------------------|--|
| Standard LEVIFLOW® Flow Sensors Compatible with Console | | | |
| Single-Use Sensor LFS-03SU-Z (Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape) | 100-30375 | LEVIFLOW LFS-SU Flowmeter | For other configurations like gamma irradiated flow sensors consult the LFS-SU product literature. |
| Single-Use Sensor LFS-03SU-Z-SC1 (Extended Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape) | 100-30418 | | |
| Single-Use Sensor LFS-06SU-Z (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape) | 100-30377 | | |
| Single-Use Sensor LFS-06SU-Z-SC1 (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape) | 100-30406 | | |
| Single-Use Sensor LFS-10SU-Z (Water 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape) | 100-30397 | | |
| Single-Use Sensor LFS-10SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape) | 100-30408 | | |
| Single-Use Sensor LFS-15SU-Z (Water 20°C, 50 l/min, Triclamp 1/2" Fittings, Z-Shape) | 100-30412 | | |
| Single-Use Sensor LFS-15SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape) | 100-30431 | | |
| Single-Use Sensor LFS-20SU.1-Z (Water 20°C, 80 l/min, Triclamp 1/2" Fittings, Z-Shape) | 100-30483 | | |
| Single-Use Sensor LFS-20SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape) | 100-30464 | | |
| LeviFlow Clamp-On Sensor LFSC-05D-001 (For Silicone/C-Flex 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min) | 100-30407 | | |
| LeviFlow Clamp-On Sensor LFSC-08D-001 (For Silicone/C-Flex 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min) | 100-30396 | | |
| LeviFlow Clamp-On Sensor LFSC-12D-007 (For Silicone/C-Flex 3/8"x9/16", Water at 20°C/37°C, 0 - 20 l/min) | 100-30390 | | |
| LeviFlow Clamp-On Sensor LFSC-22D-005 (For Silicone/C-Flex 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min) | 100-30391 | LEVIFLOW LFSC Flowmeter | For other configurations and sizes consult the LFSC-iX product literature. |
| LeviFlow Clamp-On Sensor LFSC-30D-001 (For Silicone/C-Flex 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min) | 100-30414 | | |
| LeviFlow Clamp-On Sensor LFSC-i06X-002 (For Silicone/C-Flex/TPE 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min) | 100-30485 | | |
| LeviFlow Clamp-On Sensor LFSC-i10X-001 (For Silicone/C-Flex/TPE 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min) | 100-30477 | | |
| LeviFlow Clamp-On Sensor LFSC-i16X-001 (For Silicone/C-Flex/TPE 3/8"x5/8", Water at 20°C/37°C, 0 - 20 l/min) | 100-30482 | | |
| LeviFlow Clamp-On Sensor LFSC-i19X-001 (For Silicone/C-Flex/TPE 1/2"x3/4", Water at 20°C/37°C, 0 - 50 l/min) | 100-30479 | | |
| LeviFlow Clamp-On Sensor LFSC-i25X-001 (For Silicone/C-Flex/TPE 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min) | 100-30480 | | |
| LeviFlow Clamp-On Sensor LFSC-i35X-001 (For Silicone/C-Flex/TPE 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min) | 100-30481 | LEVIFLOW LFSC Flowmeter | Calibration for 20 and 37°C water. |
| LeviFlow Clamp-On Converter LFSC-i05X.1-001 (For LFP-05SU/-05SU.1 1/8", 0-0.2/0-1 l/min, 4 Calibration Sets) | 100-30490 | | |
| LeviFlow Clamp-On Converter LFSC-i09X.1-001 (For LFP-09SU 1/4", 0-4 l/min, 2 Calibration Sets) | 100-30491 | | |
| LeviFlow Clamp-On Converter LFSC-i13X.1-001 (For LFP-13SU 3/8", 0-20 l/min, 2 Calibration Sets) | 100-30492 | | |
| LeviFlow Clamp-On Converter LFSC-i16X.1-001 (For LFP-16SU 1/2", 0-50 l/min, 2 Calibration Sets) | 100-30498 | | |
| LeviFlow Clamp-On Converter LFSC-i20X.1-001 (For LFP-20SU 3/4", 0-80 l/min, 2 Calibration Sets) | 100-30499 | | |
| LeviFlow Clamp-On Converter LFSC-i23X.1-001 (For LFP-23SU 7/8" and 1", 0-160 l/min, 4 Calibration Sets) | 100-30509 | | |
| LeviFlow Clamp-On Converter LFSC-i36X.1-001 (For LFP-36SU 13/8", 0-400 l/min, 2 Calibration Sets) | 100-30487 | | |
| Standard Pinch Valves of the LPVS-i Series Compatible with Console | | | |
| Levitronix® Pinch Valve System LPVS-i100.1 (Up to 0.375"=3/8" Tubing) | 100-30566 | LPVS Pinch Valves | Tubing inserts have to be ordered as separate articles. |
| Levitronix® Pinch Valve System LPVS-i200.1 (Up to 0.625"=5/8" Tubing) | 100-30495 | | |
| Levitronix® Pinch Valve System LPVS-i400.1 (Up to 1" Tubing) | 100-30496 | | |
| Levitronix® Pinch Valve System LPVS-i400P.1 (Up to 1.405" Tubing) | 100-30562 | | |
| Levitronix® Pinch Valve System LPVS-i800.1 (Up to 1.5" Tubing) | 100-30497 | | |
| Standard Pressure and Temperature Sensors Compatible with Console | | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-000 (Polysulfone, 1/8" Luer Fitting, Non-Sterile) | 190-10572 | Accessories Consoles | |
| Pendotech Single-Use Pressure Sensor PREPS-N-012 (Polysulfone, 1/8" Hose-Barb Fitting, Non-Sterile) | 190-10570 | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-025 (Polysulfone, 1/4" Hose-Barb Fitting, Non-Sterile) | 190-10571 | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile) | 190-10375 | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-5-5 (Polysulfone, 3/4" Sanitary Flange Fitting, Non-Sterile) | 190-10373 | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-1-1 (Polysulfone, 1" Sanitary Flange Fitting, Non-Sterile) | 190-10374 | | |
| Pendotech Single-Use Pressure Sensor PREPS-N-15-15 (Polysulfone, 1.5" Sanitary Flange Fitting, Non-Sterile) | 190-10478 | | |
| Pendotech Single-Use Temperature Sensor TEMPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile) | 190-10377 | | |
| Pendotech Single-Use Temperature Sensor TEMPS-N-050 (Polysulfone, 1/2" Sanitary Flange Fitting, Non-Sterile) | 190-10376 | | |
| Cables and Desktop AC/DC Power Supplies | | | |
| LeviFlow Interconnect Cable LFI-C.1-30 (3m, PVC) | 190-10308 | Accessories LEVIFLOW | For connecting LEVIFLOW® Sensors to the console. |
| IP Cable Signal 4 Wires ICS-3.1-30 (3m, PUR Jacket, Circular to Circular Connectors) | 190-10447 | Accessories Consoles | For connection of LFSC-iX sensors to console. |
| Pressure Sensor Adaptor Cable LPE-1.1-30 (3m, Circular Console to Pressure Sensor Connector) | 190-10354 | Accessories Consoles | For pressure sensors. |
| Temperature Sensor Adaptor Cable LTE-1.1-30 (3m, Circular Console to Temperature Sensor Connector) | 190-10353 | Accessories Consoles | For temperature sensors. |
| IP Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector) | 190-10347 | Accessories IPS | For connection of devices to PLC of console. |
| IP Adaptor Cable Hybrid ICH-1.2-30 (3m, PVC-Jacket, Circular to Circular) | 190-10329 | Accessories IPS | For extension of driver cable. |
| IPS Split Cable ICY-4.1-02 (0.2m, PVC cables, PC Split Box, PUR Sealing, circ. Connectors, Integrated 48 to 24VDC converter) | 190-10476 | Accessories IPS | Splits cable to connect power supply and user panel to IPD-600 driver. |
| IP Cable Signal 4 Wires ICS-1.5-30 (3m PVC, Circular Hirose to Console Connector) | 190-10468 | Accessories Consoles | For connecting split box to console power connector. |
| IP Split Cable ICY-2.1-02 (PVC cable jacket with PA+GF split block, TPE sealing, circular type connectors) | 190/10466 | Accessories IPS | Split pinch valves to power supply and console |
| IP Cable Signal 12 Wires ICS-6.3-30 (3m PVC, Circular to Circular Connector) | 190-10551 | Accessories LPVS | Connection of pinch valves to split cable ICP-Y-2.1. |
| IP Cable Signal 12 Wires ICS-6.3-50 (3m) | 190-10542 | | |
| IP Cable Signal 5 Wires ICS-3.3-30 (3m PVC, Circular to Circular Connector) | 190-10583 | Accessories Consoles | Connection of pinch valves or other RS485 devices to Remote port of console. |
| IP Cable Signal 12 Wires ICS-2.8-30 (3m PVC, Circular to Circular Connector) | 190-10580 | Accessories Consoles | For connection of PLC of two consoles. |
| IP Cable Signal 12 Wires ICS-2.5-30 (3m PVC, Circular to Circular Connector) | 190-10554 | Accessories Consoles | For connection of EasyConnect drivers of i30/i100 pump systems to console PLC. |
| Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC/50W, 90-264VAC, Circular Intercontec Connector) | 100-40019 | Accessories IPS | AC mains cable ACM-1.x for various countries to be ordered as seperate article. |
| Desktop AC/DC Power Supply GST160A24-R7B IC915 (24VDC/160W, 85-264VAC, Circular Intercontec Connector) | 100-40020 | Accessories IPS | |
| AC Mains Cable AMC-1.1 (2m, US, Canada) | 190-10331 | Accessories Consoles | Mains cables for connecting console to AC supply. Cables for other countries on request. |
| AC Mains Cable AMC-1.2 (2.5m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria) | 190-10332 | | |
| AC Mains Cable AMC-1.3 (2.5m, Japan) | 190-10333 | | |
| AC Mains Cable AMC-1.4 (2.5m, Switzerland) | 190-10334 | | |
| AC Mains Cable AMC-1.5 (2.5m, United Kingdom) | 190-10335 | | |
| AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC , Circ. Connectors) | 100-40026 | Accessories IPS | AC mains cable AMC-2.x for various countries to be ordered as seperate article. |
| AC Mains Cable AMC-2.1 (3m, US, Canada) | 190-10336 | Accessories IPS | Mains cable for AC/DC power supply DDP-600-US48-SC |
| AC Mains Cable AMC-2.2 (3m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria) | 190-10337 | | |
| AC Mains Cable AMC-2.3 (3m, Japan) | 190-10338 | | |
| AC Mains Cable AMC-2.4 (3m, Switzerland) | 190-10339 | | |
| AC Mains Cable AMC-2.5 (3m, United Kingdom) | 190-10340 | | |