

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® 200SU (Single-Use), Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note
Driver System Packages (Motor with Pumphead Socket and Controller)			
PuraLev Driver System PLD-200SU.1 (Standalone) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminum, Bayonet Type) Controller LPC-200.1-01 (24V/200W, Integrated Display, PLC)	100-90647	PuraLev SU	Note: Motor has direct cabling to controller.
PuraLev Driver System PLD-200SU.2 (PLC/USB) or PLD-200SU.3 (PLC/RS232) Motor BSM-1.4 (5m FEP Cable, D-Sub Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminum, Bayonet Type) Controller LPC-200.2-01 (24V/200W, PLC, USB) or LPC-200.3-01 (PLC, RS232)	100-90648 or 100-90649	PuraLev SU	Note: Motor has direct cabling to controller.
PuraLev Driver System PLD-200SU.4 (ATEX, Standalone) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Hummel Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminum, Bayonet Type) Controller LPC-200.1-01 (24V/200W, Integrated Display, PLC)	100-90650	PuraLev SU	Note: To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m.
PuraLev Driver System PLD-200SU.5 (PLC/USB) or PLD-200SU.6 (PLC/RS232) Motor BSM-1.6 (ATEX, 2m FEP Cable, Circular Hummel Connector, ETFE Housing Coating) Pumphead Socket PHS-200.1 (Anodized Aluminum, Bayonet Type) Controller LPC-200.2-01 (24V/200W, PLC, USB) or LPC-200.3-01 (PLC, RS232)	100-90651 or 100-90652	PuraLev SU	Note: To be ordered as separate article: Extension Cables MCA-1.x-yy -> 0.5-10m.
Single-Use Pumpheads			
Single-Use Pumphead DCP-200.2 Welded and Molded PP Housing, Triclamp 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile	100-90734	PuraLev SU	
Single-Use Pump Head Gamma Irradiated DCP-200.2-G25 (Applied Gamma Dosage ≥ 25 kGy)	100-91078	PuraLev SU	
Single-Use Pumphead DCP-200.3 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile	100-90792	PuraLev SU	
Single-Use Pump Head Gamma Irradiated DCP-200.3-G25 (Applied Gamma Dosage ≥ 25 kGy)	100-91122	PuraLev SU	
Single-Use Pumphead DCP-200.3-SF1-G25 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Sterile Fittings, Applied Gamma Dosage ≥ 25 kGy)	100-91233	PuraLev SU	
Accessories			
Power Supply TSP 360-124 (Traco, 24V / 360W)	100-40016	Accessories Pump Systems	
Power Supply TSP 360-148 (Traco, 48V / 360W)	100-40017	Accessories Pump Systems	
Mounting Base Plate MBP-1.1 (PVDF)	190-10004	Accessories Pump Systems	
Air Cooling Module ACM-1.1 (PVDF, NPT1/4")	190-10003	Accessories Pump Systems	Discontinued by September 2023. Alternative is ACM-1.3.
Air Cooling Module ACM-1.3 (PP, NPT1/4")	190-10193	Accessories Pump Systems	
ATEX Cable Sealing System ACS-A.1	100-90292	Accessories Pump Systems	
User Panel LUI-A.1 (RS232 interface)	100-30300	Accessories Pump Systems	Used with LPC-200.3-01
Adapter/Extension Cable MCA-1.2-05 (0.5m, FEP-Jacket, Circular/D-SUB) Adapter/Extension Cable MCA-1.2-30 (3m) Adapter/Extension Cable MCA-1.2-50 (5m) Adapter/Extension Cable MCA-1.2-70 (7m) Adapter/Extension Cable MCA-1.2-100 (10m)	190-10147 190-10092 190-10128 190-10148 190-10149	Accessories Pump Systems	
Adapter/Extension Cable MCA-1.3-05 (0.5m, FEP-Jacket, Circular Wallmountable/D-SUB) Adapter/Extension Cable MCA-1.3-30 (3m) Adapter/Extension Cable MCA-1.3-50 (5m) Adapter/Extension Cable MCA-1.3-70 (7m) Adapter/Extension Cable MCA-1.3-100 (10m)	190-10234 190-10237 190-10200 190-10235 190-10236	Accessories Pump Systems	
Adapter/Extension Cable MCA-1.4-05 (0.5m, PVC-Jacket, Circular/D-SUB) Adapter/Extension Cable MCA-1.4-30 (3m) Adapter/Extension Cable MCA-1.4-50 (5m) Adapter/Extension Cable MCA-1.4-70 (7m) Adapter/Extension Cable MCA-1.4-100 (10m)	190-10228 190-10170 190-10187 190-10229 190-10230	Accessories Pump Systems	
Adapter/Extension Cable MCA-1.5-05 (0.5m, PVC-Jacket, Circular Wallmountable/D-SUB) Adapter/Extension Cable MCA-1.5-30 (3m) Adapter/Extension Cable MCA-1.5-50 (5m) Adapter/Extension Cable MCA-1.5-70 (7m) Adapter/Extension Cable MCA-1.5-100 (10m)	190-10225 190-10231 190-10199 190-10232 190-10233	Accessories Pump Systems	
Spare Parts			
Motor BSM-1.4 (ETFE Coating, 5m FEP Cable with D-SUB Connector)	100-10005	BPS-200	
Motor BSM-1.6 (ATEX) (ETFE Coating, 2m FEP Cable with Circular Connector, ATEX)	100-10063	PuraLev MU	
Motor BSM-1.8 (ATEX) (Epoxy Coating, 2m PVC Cable with Circular Connector, ATEX)	100-10090	PuraLev SU	
Pumphead Socket PHS-200.1 (Anodized Aluminium, Bayonet Socket Type)	100-90695	PuraLev SU	
Controller LPC-200.1-01 (24V/200W, Display, Enable connector) (included: Enable connector)	100-30027	PuraLev SU	
Controller LPC-200.2-01 (24V/200W, USB, PLC interface) (included: PLC connector)	100-30028	PuraLev SU	
Controller LPC-200.3-01 (24V/200W, RS232, PLC interface) (included: PLC connector)	100-30029	PuraLev SU	