

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem DuraLev® 2000, Accessories and Spare Parts

*Controllers LPC-2000.x are discontinued by 30-Jun-2026 (product notification PCN #PR-2505-79) and replaced by the LPC-2000S.x models with newest design.
 Note that all new controller models LPC-2000S.x can run with the same firmware as the old LPC-2000.x models.
 The pump systems DuraLev-2000.x are hence replaced by the new pump systems DuraLev-2000.xS (see separate article list) with the same index "x" but new article number.*

Article Name (and Description)	Part No.	Article Category	Note
Systems			
Pump System DuraLev-2000.1 Motor LPM-2000.1 (3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.1-01 (Standalone Display, Enable Connector) Pump Head LPP-2000.4 (PVDF Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.1 and MCAP-2000.1 (1-10m) to be ordered as separate article	100-90491	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.1S (100-91911).</i>
Pump System DuraLev-2000.2 or DuraLev-2000.3 Motor LPM-2000.1 (3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.2-01 (USB, Extended PLC Interface) or LPC-2000.3-01 (RS232, Ext. PLC) Pump Head LPP-2000.4 (PVDF Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.1 and MCAP-2000.1 (1-10m) to be ordered as separate article	100-90492 or 100-90493	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.2S (100-91912) and DuraLev-2000.3S (100-91913) respectively.</i>
Pump System DuraLev-2000.4 Motor LPM-2000.1 (3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.1-01 (Standalone Display, Enable Connector) Pump Head LPP-2000.2 (PP+GF30 Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.1 and MCAP-2000.1 (1-10m) to be ordered as separate article	100-90494	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.4S (100-91914).</i>
Pump System DuraLev-2000.5 or DuraLev-2000.6 Motor LPM-2000.1 (3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.2-01 (USB, Extended PLC Interface) or LPC-2000.3-01 (RS232, Ext. PLC) Pump Head LPP-2000.2 (PP+GF30 Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.1 and MCAP-2000.1 (1-10m) to be ordered as separate article	100-90495 or 100-90496	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.5S (100-91915) and DuraLev-2000.6S (100-91916) respectively.</i>
Pump System DuraLev-2000.7 (ATEX) Motor LPM-2000.7 (ATEX, 3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.1-01 (Standalone Display, Enable Connector) Pump Head LPP-2000.4 (PVDF Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.3 and MCAP-2000.3 (1-10m) to be ordered as separate article	100-90497	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.7S (100-91917).</i>
Pump System DuraLev-2000.8 or DuraLev-2000.9 (ATEX) Motor LPM-2000.7 (ATEX, 3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.2-01 (USB, Extended PLC Interface) or LPC-2000.3-01 (RS232, Ext. PLC) Pump Head LPP-2000.4 (PVDF Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.3 and MCAP-2000.3 (1-10m) to be ordered as separate article	100-90498 or 100-90499	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.8S (100-91918) and DuraLev-2000.9S (100-91919) respectively.</i>
Pump System DuraLev-2000.10 (ATEX) Motor LPM-2000.7 (ATEX, 3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.1-01 (Standalone Display, Enable Connector) Pump Head LPP-2000.2 (PP+GF30 Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.3 and MCAP-2000.3 (1-10m) to be ordered as separate article	100-90532	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.10S (100-91920).</i>
Pump System DuraLev-2000.11 or DuraLev-2000.12 (ATEX) Motor LPM-2000.7 (ATEX, 3m PVC Cable, Epoxy Coating, AMP Connector) Controller LPC-2000.2-01 (USB, Extended PLC Interface) or LPC-2000.3-01 (RS232, Ext. PLC) Pump Head LPP-2000.2 (PP+GF30 Housing, NPT 1", PFA-Impeller, FPM O-Ring) Extension Cables MCAS-600.3 and MCAP-2000.3 (1-10m) to be ordered as separate article	100-90533 or 100-90534	DuraLev	<i>Discontinued. Replacement system is DuraLev-2000.11S (100-91921) and DuraLev-2000.12S (100-91922) respectively.</i>

Article Name (and Description)	Part No.	Article Category	Note
Accessories			
Air Cooling Module ACM-4.2 (PP, NPT1/4", Molded)	190-10139	Accessories Pump Systems	
Air Cooling Module ACM-4.3 (PP-EL-S, NPT1/4", Machined, for ATEX Motors)	190-10243	Accessories Pump Systems	
Fan Cooling Module FCM-2000.1 (PP+Talk Cover, 24V/33.5W IP68 Fan, 3m PP Cable/M12 Connector)	190-10390	Accessories Pump Systems	
Fan Cooling Module Cable, FCC-1.1-50 (5.0m, PP Cable, M12 Connector to Open Wire)	190-10407	Accessories	
Fan Cooling Module Cable, FCC-1.1-100 (10.0m)	190-10408	Pump Systems	
ATEX Cable Sealing System ACS-A.1	100-90292	Accessories Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.1-05 (0.5m, PVC-Jacket, Circular AMP/COMBICON)	190-10208		
Adaptor/Extension Power Cable, MCAP-2000.1-30 (3.0m)	190-10210		
Adaptor/Extension Power Cable, MCAP-2000.1-50 (5.0m)	190-10211	Accessories	
Adaptor/Extension Power Cable, MCAP-2000.1-70 (7.0m)	190-10205	Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.1-100 (10.0m)	190-10212		
Adaptor/Extension Sensor Cable, MCAS-600.1-05 (0.5m, PVC-Jacket, Circular AMP/D-SUB)	190-10122		
Adaptor/Extension Sensor Cable, MCAS-600.1-30 (3.0m)	190-10123		
Adaptor/Extension Sensor Cable, MCAS-600.1-50 (5.0m)	190-10124	Accessories	
Adaptor/Extension Sensor Cable, MCAS-600.1-70 (7.0m)	190-10101	Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.1-100 (10.0m)	190-10125		
Adaptor/Extension Power Cable, MCAP-2000.2-05 (0.5m) (PVC-Jacket, Circ. Wall-Mountable M23/COMBICON)	190-10213		
Adaptor/Extension Power Cable, MCAP-2000.2-30 (3.0m)	190-10215		
Adaptor/Extension Power Cable, MCAP-2000.2-50 (5.0m)	190-10216	Accessories	
Adaptor/Extension Power Cable, MCAP-2000.2-70 (7.0m)	190-10217	Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.2-100 (10.0m)	190-10218		
Adaptor/Extension Sensor Cable, MCAS-600.2-05 (0.5m) (PVC-Jacket, Circular Wall-Mountable M23/D-SUB)	190-10226		
Adaptor/Extension Sensor Cable, MCAS-600.2-30 (3.0m)	190-10238		
Adaptor/Extension Sensor Cable, MCAS-600.2-50 (5.0m)	190-10127	Accessories	
Adaptor/Extension Sensor Cable, MCAS-600.2-70 (7.0m)	190-10105	Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.2-100 (10.0m)	190-10239		
Adaptor/Extension Power Cable, MCAP-2000.3-05 (0.5m) (PVC-Jacket, Circular M23/COMBICON)	190-10219		
Adaptor/Extension Power Cable, MCAP-2000.3-30 (3.0m)	190-10221		
Adaptor/Extension Power Cable, MCAP-2000.3-50 (5.0m)	190-10222	Accessories	
Adaptor/Extension Power Cable, MCAP-2000.3-70 (7.0m)	190-10223	Pump Systems	
Adaptor/Extension Power Cable, MCAP-2000.3-100 (10.0m)	190-10224		
Adaptor/Extension Sensor Cable, MCAS-600.3-05 (0.5m) (PVC-Jacket, Circular Hummel M23/D-SUB)	190-10158		
Adaptor/Extension Sensor Cable, MCAS-600.3-30 (3.0m)	190-10159		
Adaptor/Extension Sensor Cable, MCAS-600.3-50 (5.0m)	190-10130	Accessories	
Adaptor/Extension Sensor Cable, MCAS-600.3-70 (7.0m)	190-10160	Pump Systems	
Adaptor/Extension Sensor Cable, MCAS-600.3-100 (10.0m)	190-10161		
Spare Parts			
Impeller Exchange Kit IEK-2000.3 (for DuraLev-2000) (LPI-2000.2, 8x SS Pump screws, 4 SS Motor Screws, Kalrez O-Ring, Impeller Exchange Tool IET-3.1)	100-90531	Spare Parts Pump Systems	
Motor LPM-2000.1 (Epoxy Coating, 3m FEP Cable with AMP Connectors)	100-10049	BPS-2000	
Motor LPM-2000.7 (ATEX) (ATEX, Epoxy Coating, 3m FEP Cable with AMP Connectors)	100-10059	PuraLev MU	
Controller LPC-2000.1-01 (Stand-alone, 2kW)	100-30018	BPS-2000	Includes power supply and PLC connector. Discontinued: Replacement article is LPC-2000S.1-01 (100-30164)
Controller LPC-2000.2-01 (USB and PLC Interface, 2kW)	100-30021	BPS-2000	Includes power supply and PLC connector. Discontinued: Replacement article is LPC-2000S.2-01 (100-30136)
Controller LPC-2000.3-01 (RS232 and PLC Interface, 2kW)	100-30024	BPS-2000	Includes power supply and PLC connector. Discontinued: Replacement article is LPC-2000S.3-01 (100-30127)
Pump Head LPP-2000.2 (PP+GF30) (PP+GF30 Housing, NPT 1" Fitting, PFA-Impeller LPI-2000.1, FPM Sealing O-Ring, SS Screws)	100-90414	Spare Parts Pump Systems	
Pump Head LPP-2000.4 (PVDF) (PVDF Housing, NPT 1" Fitting, PFA-Impeller LPI-2000.1, FPM Sealing O-Ring, SS Screws)	100-90416	Spare Parts Pump Systems	