

**Levitronix® Standard Article List and Order Information**  
**PuraLev® Console LCO-i100, System Packages, Components and Accessories**

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard System Packages with PuraLev® Consoles</b>			
<b>PuraLev® Console System PCS-i30SU.1</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i30.1 (Anod. Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91098	PuraLev Console	
<b>PuraLev® Console System PCS-i100SU.1</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector) and Pump Head Socket PHS-i100.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91097	PuraLev Console	
<b>PuraLev® Console System PCS-i600SU.1</b> Console LCO-100.1-01 with AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC) Driver IPD-600.3-30-02 (Epoxy Coat., 3m PVC Cable, Connector) and Pump Head Socket PHS-600.1 (Anodized Alu) <i>Note 1: AC mains cable AMC-2.x for various countries have to be ordered as separate article.</i> <i>Note 2: Split cable ICY-4.1-02 and ICS-1.5-xx to connect power supply and driver to console must be ordered as separate article.</i>	100-91494	PuraLev Console	
<b>PuraLev® Console System PCS-i30SU.5 or PCS-i30SU.3 or PCS-i30SU.4</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i30.1 (Anod. Alu) LeviFlow Clamp-On Sensor LFSC-05D-001 (1 l/min) or LFSC-08D-001 (4 l/min) or LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C-1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91349 or 100-91099 or 100-91100	PuraLev Console	
<b>PuraLev® Console System PCS-i100SU.2 or PCS-i100SU.3</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-02 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) and Pump Head Socket PHS-i100.1 (Anod. Alu) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) or LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C-1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91096 or 100-91092	PuraLev Console	
<b>PuraLev® Console System PCS-i30MU.1</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC, 50W, 90-264VAC) Driver IPD-30.8-03-06 (Epoxy Coat., 0.3m PVC Cable, Connector, MBP-i30.1) Multi-Use Pump Head LPP-30.7 (Triclamp 3/8", PEEK/PFA/EPDM, No drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91562	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.1 or PCS-i100MU.2</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i>	100-91268 or 100-91269	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.3 or PCS-i100MU.4</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LeviFlow Clamp-On Sensor LFSC-08D-001 (4 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C-1-xx for connecting LFSC-D Clamp-On Sensor to console to be ordered as separate article.</i>	100-91270 or 100-91271	PuraLev Console	
<b>PuraLev® Console System PCS-i100MU.5 or PCS-i100MU.6</b> Console LCO-100.1-01 with Desktop AC/DC Power Supply GST160A24 IC915 (24VDC, 160W, 85-264 VAC) Driver IPD-100.3-03-03 (Epoxy Coat., 0.3m PVC Cable, Connector) Multi-Use Pump Head LPP-200.16 (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or LPP-200.17 (With drain port) LeviFlow Clamp-On Sensor LFSC-12D-007 (20 l/min) <i>Note 1: AC mains cable AMC-1.x for various countries have to be ordered as separate article.</i> <i>Note 2: Cable ICH-1.2-x for extension of driver cable to be ordered as separate article.</i> <i>Note 3: Adaptor cables LFI-C-1-xx for connecting Clamp-On Sensor to console to be ordered as separate article.</i>	100-91272 or 100-91273	PuraLev Console	
<b>Note: PuraLev Console System Packages with LFSC-IX flowmeters on request. Adaptor cable ICS-3.1 is needed to connect LFSC-IX clamp-on flowmeters to console.</b>			
<b>PuraLev® Consoles</b>			
<b>PuraLev® Console LCO-i100.1-01</b> (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i>	100-30501	PuraLev Console	
<b>PuraLev® Console LCO-i100.1-03 - FDA CFR 21 Part 11 Ready</b> (24 VDC, 100W, Interfaces: 3 Pressure Sensors, 2 Flow Sensor, 1 Temperature Sensor, PLC, RS485, LAN, USB) <i>Note: AC/DC desktop supply to be ordered as separate article.</i>	100-30539	PuraLev Console	
<b>Standard Pump Drivers and Pump Heads Compatible with Console</b>			
<b>Single-Use Pump Head DCP-30.1</b> (Triclamp 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-90959	PuraLev SU	PuraLev® i30SU/MU pump system components.
<b>Single-Use Pump Head DCP-30.1-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91170	PuraLev SU	
<b>Single-Use Pump Head DCP-30.2</b> (Barb 3/8", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-90968	PuraLev SU	
<b>Single-Use Pump Head DCP-30.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91071	PuraLev SU	
<b>Single-Use Pump Head DCP-30.2-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91234	PuraLev SU	
<b>Single-Use Pump Head DCP-30.5</b> (Barb 1/4", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91691	PuraLev SU	
<b>Single-Use Pump Head DCP-30.5-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i30.1</i>	100-91700	PuraLev SU	
<b>Multi-Use Pump Head LPP-30.7</b> (Triclamp 3/8", PEEK/PFA/EPDM, No drain port)	100-91557	PuraLev MU	
<b>Driver IPD-30.8-03-02</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, SU Firmware)	100-10104	Spare Parts Pump System	
<b>Driver IPD-30.8-03-06</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connectors, MU Firmware)	100-10212	Spare Parts Pump System	
<b>Pumphead Socket PHS-i30.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-90947	Spare Parts Pump System	
<b>Single-Use Pump head DCP-200.2</b> (Triclamp 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-90734	PuraLev SU	PuraLev® i100SU/MU pump system components.
<b>Single-Use Pump Head DCP-200.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91078	PuraLev SU	
<b>Single-Use Pump Head DCP-200.3</b> (Barb 1/2", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-90792	PuraLev SU	
<b>Single-Use Pump Head DCP-200.3-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91122	PuraLev SU	
<b>Single-Use Pump Head DCP-200.3-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-i100.1</i>	100-91233	PuraLev SU	
<b>Multi-Use Pump Head LPP-200.16</b> (Triclamp 1/2", PVDF/PFA/EPDM, No drain port) or <b>LPP-200.17</b> (With drain port)	100-90863 or 100-90864	PuraLev MU	For other configurations consult the PuraLev® i100SU/MU product literature.
<b>Driver IPD-100.3-03-02</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connector, SU Firmware)	100-10105	Spare Parts Pump System	
<b>Driver IPD-100.3-03-03</b> (Epoxy Coating, 0.3m PVC Cable with Circular Connector, MU Firmware)	100-10137	Spare Parts Pump System	
<b>Pumphead Socket PHS-i100.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-91053	Spare Parts Pump System	
<b>Single-Use Pump Head DCP-600.2</b> (Triclamp 1", PP, Non-Sterile) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-90784	PuraLev SU	PuraLev® i600SU pump system components.
<b>Single-Use Pump Head DCP-600.2-G25</b> (Gamma Dosage ≥ 25 kGy) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91079	PuraLev SU	
<b>Single-Use Pump Head DCP-600.2-SF1-G25</b> (Sterile Fittings) <i>Note: Needs Pump Head Socket PHS-600.1</i>	100-91294	PuraLev SU	
<b>Driver IPD-600.3-30-02</b> (Epoxy Coating, 3m PVC Cable with Circular Connector, SU Firmware)	100-10197	Spare Parts IPS	
<b>Pumphead Socket PHS-600.1</b> (Anodized Aluminium, Bayonet Socket Type)	100-90696	Spare Parts Pump Systems	

Article Name (and Description)	Part No.	Article Category	Note
<b>Standard LEVIFLOW® Flow Sensors Compatible with Console</b>			
Single-Use Sensor LFS-03SU-Z (Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30375	LEVIFLOW LFS-SU Flowmeter	For other configurations like gamma irradiated flow sensors consult the LFS-SU product literature.
Single-Use Sensor LFS-03SU-Z-SC1 (Extended Calibration at 20°C, 0.8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30418		
Single-Use Sensor LFS-06SU-Z (Calibration at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30377		
Single-Use Sensor LFS-06SU-Z-SC1 (Extended Calibr. at 20°C, 8 l/min, Triclamp 3/8" Fittings, Z-Shape)	100-30406		
Single-Use Sensor LFS-10SU-Z (Water 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30397		
Single-Use Sensor LFS-10SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30408		
Single-Use Sensor LFS-15SU-Z (Water 20°C, 50 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30412		
Single-Use Sensor LFS-15SU-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape)	100-30431		
Single-Use Sensor LFS-20SU-1-Z (Water 20°C, 80 l/min, Triclamp 1/2" Fittings, Z-Shape)	100-30483		
Single-Use Sensor LFS-20SU-1-Z-SC1 (Extended Calibr at 20°C, 20 l/min, Triclamp 1" Fittings, Z-Shape)	100-30464		
LeviFlow Clamp-On Sensor LFSC-05D-001 (For Silicone/C-Flex 1/4"x3/8", Water at 20°C/37°C, 0 - 1 l/min)	100-30407	Spare Parts LEVIFLOW	For other configurations and sizes consult the LFSC-D product literature.
LeviFlow Clamp-On Sensor LFSC-05D-001 (For Silicone/C-Flex 1/4"x3/8", Water at 20°C/37°C, 0 - 1 l/min)	100-30396		
LeviFlow Clamp-On Sensor LFSC-12D-007 (For Silicone/C-Flex 3/8"x9/16", Water at 20°C/37°C, 0 - 20 l/min)	100-30390		
LeviFlow Clamp-On Sensor LFSC-12D-007 (For Silicone/C-Flex 3/8"x9/16", Water at 20°C/37°C, 0 - 20 l/min)	100-30391		
LeviFlow Clamp-On Sensor LFSC-22D-005 (For Silicone/C-Flex 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)	100-30414		
LeviFlow Clamp-On Sensor LFSC-30D-001 (For Silicone/C-Flex 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min)	100-30414		
LeviFlow Clamp-On Sensor LFSC-106X-002 (For Silicone/C-Flex/TPE 1/8"x1/4", Water at 20°C/37°C, 0 - 1 l/min)	100-30485		
LeviFlow Clamp-On Sensor LFSC-110X-004 (For Silicone/C-Flex/TPE 1/4"x3/8", Water at 20°C/37°C, 0 - 4 l/min)	100-30477		
LeviFlow Clamp-On Sensor LFSC-116X-001 (For Silicone/C-Flex/TPE 3/8"x5/8", Water at 20°C/37°C, 0 - 20 l/min)	100-30482		
LeviFlow Clamp-On Sensor LFSC-119X-001 (For Silicone/C-Flex/TPE 1/2"x3/4", Water at 20°C/37°C, 0 - 50 l/min)	100-30479		
LeviFlow Clamp-On Sensor LFSC-125X-001 (For Silicone/C-Flex/TPE 3/4"x1", Water at 20°C/37°C, 0 - 80 l/min)	100-30480	LEVIFLOW LFSC Flowmeter	For other configurations and sizes consult the LFSC-IX product literature.
LeviFlow Clamp-On Sensor LFSC-135X-001 (For Silicone/C-Flex/TPE 1"x13/8", Water at 20°C/37°C, 0 - 160 l/min)	100-30481		
LeviFlow Clamp-On Converter LFSC-105X-1-001 (For LFP-05SU-05SU.1 1/8", 0-0.2/0-1 l/min, 4 Calibration Sets)	100-30490		
LeviFlow Clamp-On Converter LFSC-109X-1-001 (For LFP-09SU 1/4", 0-4 l/min, 2 Calibration Sets)	100-30491		
LeviFlow Clamp-On Converter LFSC-113X-1-001 (For LFP-13SU 3/8", 0-20 l/min, 2 Calibration Sets)	100-30492		
LeviFlow Clamp-On Converter LFSC-116X-1-001 (For LFP-16SU 1/2", 0-50 l/min, 2 Calibration Sets)	100-30498		
LeviFlow Clamp-On Converter LFSC-120X-1-001 (For LFP-20SU 3/4", 0-80 l/min, 2 Calibration Sets)	100-30499		
LeviFlow Clamp-On Converter LFSC-123X-1-001 (For LFP-23SU 7/8" and 1", 0-160 l/min, 4 Calibration Sets)	100-30509		
LeviFlow Clamp-On Converter LFSC-136X-1-001 (For LFP-36SU 13/8", 0-400 l/min, 2 Calibration Sets)	100-30487		
<b>Standard Pinch Valves of the LPVS-i Series Compatible with Console</b>			
Levitronix® Pinch Valve System LPVS-i100.1 (Up to 0.375"-3/8" Tubing)	100-30566	LPVS Pinch Valves	Tubing inserts have to be ordered as separate articles.  For accessories and more detailed specifications consult the LPVS-i product literature.
Levitronix® Pinch Valve System LPVS-i100.1 (Up to 0.625"-5/8" Tubing)	100-30495		
Levitronix® Pinch Valve System LPVS-i400.1 (Up to 1" Tubing)	100-30496		
Levitronix® Pinch Valve System LPVS-i400P.1 (Up to 1.405" Tubing)	100-30562		
Levitronix® Pinch Valve System LPVS-i800.1 (Up to 1.5" Tubing)	100-30497		
<b>Standard Pressure and Temperature Sensors Compatible with Console</b>			
Pendotech Single-Use Pressure Sensor PREPS-N-000 (Polysulfone, 1/8" Luer Fitting, Non-Sterile)	190-10572	Accessories Consoles	Legacy product from Pendotech. Replaced by 190-10645.
Pendotech Single-Use Pressure Sensor 30961415 (Polysulfone, 1/8" Luer Fitting, Non-Sterile)	190-10645		
Pendotech Single-Use Pressure Sensor PREPS-N-012 (Polysulfone, 1/8" Hose-Barb Fitting, Non-Sterile)	190-10570		
Pendotech Single-Use Pressure Sensor PREPS-N-025 (Polysulfone, 1/4" Hose-Barb Fitting, Non-Sterile)	190-10571		
Pendotech Single-Use Pressure Sensor PREPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)	190-10375		
Pendotech Single-Use Pressure Sensor PREPS-N-5-5 (Polysulfone, 3/4" Sanitary Flange Fitting, Non-Sterile)	190-10373		
Pendotech Single-Use Pressure Sensor PREPS-N-1-1 (Polysulfone, 1" Sanitary Flange Fitting, Non-Sterile)	190-10374		
Pendotech Single-Use Pressure Sensor PREPS-N-15-15 (Polysulfone, 1.5" Sanitary Flange Fitting, Non-Sterile)	190-10478		
Pendotech Single-Use Temperature Sensor TEMPS-N-038 (Polysulfone, 3/8" Hose-Barb Fitting, Non-Sterile)	190-10377		
Pendotech Single-Use Temperature Sensor TEMPS-N-050 (Polysulfone, 1/2" Sanitary Flange Fitting, Non-Sterile)	190-10376		
<b>Cables and Desktop AC/DC Power Supplies</b>			
LeviFlow Interconnect Cable LFI-C-1-30 (3m, PVC)	190-10308	Accessories LEVIFLOW	For connecting LEVIFLOW® Sensors to the console.
IP Cable Signal 4 Wires ICS-3.1-30 (3m, PUR Jacket, Circular to Circular Connectors)	190-10447	Accessories Consoles	For connection of LFSC-IX sensors to console.
Pressure Sensor Adaptor Cable LPE-1.1-30 (3m, Circular Console to Pressure Sensor Connector)	190-10354	Accessories Consoles	For pressure sensors.
Temperature Sensor Adaptor Cable LTE-1.1-30 (3m, Circular Console to Temperature Sensor Connector)	190-10353	Accessories Consoles	For temperature sensors.
IP Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connection of devices to PLC of console.
IP Adaptor Cable Hybrid ICH-1.2-30 (3m, PVC-Jacket, Circular to Circular)	190-10329	Accessories IPS	For extension of driver cable.
IPS Split Cable ICY-4.1-02 (0.2m, PVC cables, PC Split Box, PUR Sealing, circ. Connectors, Integrated 48 to 24VDC converter)	190-10476	Accessories IPS	Splits cable to connect power supply and user panel to IPD-600 driver.
IP Cable Signal 4 Wires ICS-1.5-30 (3m PVC, Circular Hirose to Console Connector)	190-10468	Accessories Consoles	For connecting split box to console power connector.
IP Split Cable ICY-2.1-02 (PVC cable jacket with PA+GF split block, TPE sealing, circular type connectors)	190/10466	Accessories IPS	Split pinch valves to power supply and console
IP Cable Signal 12 Wires ICS-6.3-30 (3m PVC, Circular to Circular Connector)	190-10551	Accessories LPVS	Connection of pinch valves to split cable ICP-Y-2.1.
IP Cable Signal 12 Wires ICS-6.3-50 (3m)	190-10542		
IP Cable Signal 5 Wires ICS-3.3-30 (3m PVC, Circular to Circular Connector)	190-10583	Accessories Consoles	Connection of pinch valves or other RS485 devices to Remote port of console.
IP Cable Signal 12 Wires ICS-2.8-30 (3m PVC, Circular to Circular Connector)	190-10580	Accessories Consoles	For connection of PLC of two consoles.
IP Cable Signal 12 Wires ICS-2.5-30 (3m PVC, Circular to Circular Connector)	190-10554	Accessories Consoles	For connection of EasyConnect drivers of i30/i100 pump systems to console PLC.
Desktop AC/DC Power Supply VEC50US24 IC915 (24VDC/50W, 90-264VAC, Circular Intercontec Connector)	100-40019	Accessories IPS	AC mains cable ACM-1.x for various countries to be ordered as separate article.
Desktop AC/DC Power Supply GST160A24-R7B IC915 (24VDC/160W, 85-264VAC, Circular Intercontec Connector)	100-40020	Accessories IPS	
AC Mains Cable AMC-1.1 (2m, US, Canada)	190-10331	Accessories Consoles	Mains cables for connecting console to AC supply. Cables for other countries on request.
AC Mains Cable AMC-1.2 (2.5m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria)	190-10332		
AC Mains Cable AMC-1.3 (2.5m, Japan)	190-10333		
AC Mains Cable AMC-1.4 (2.5m, Switzerland)	190-10334		
AC Mains Cable AMC-1.5 (2.5m, United Kingdom)	190-10335		
AC/DC Power Supply DDP-600-US48-SC (48VDC/600W, 100-277VAC, Circ. Connectors)	100-40026	Accessories IPS	AC mains cable AMC-2.x for various countries to be ordered as separate article.
AC Mains Cable AMC-2.1 (3m, US, Canada)	190-10336	Accessories IPS	Mains cable for AC/DC power supply DDP-600-US48-SC
AC Mains Cable AMC-2.2 (3m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria)	190-10337		
AC Mains Cable AMC-2.3 (3m, Japan)	190-10338		
AC Mains Cable AMC-2.4 (3m, Switzerland)	190-10339		
AC Mains Cable AMC-2.5 (3m, United Kingdom)	190-10340		