

Declaration of Conformity for CE/UKCA Marking



Machinery Directive 2014/42/EU
 Pressure Equipment Directive 2014/68/EU
 EMC Directive 2014/30/EU

UK Regulation SI 2008 Nr. 1597
 UK Regulation SI 2016 Nr. 1105
 UK Regulation SI 2016 Nr. 1091

Issuer: Levitronix GmbH
 Bändliweg 30
 CH-8048 Zürich, Switzerland
 Tel.: +41-44-974 40 00

We herewith declare that the **LEVIFLOW®** single-use flowmeters, in its various configurations as listed below is in conformity with the mentioned *European Directives* and *UK Statutory Instruments*. The flowmeter systems are thought to be used in hydraulic circuits of the production equipment of Bioprocessing and Pharmaceutical industry.

Part Name	Description
LFS-xx-SU	Flowsensor Series (rated for 5bar): LFS-03SU (0 – 0.8 lpm), LFS-06SU (0 – 8 lpm), LFS-10SU (0 – 20 lpm), LFS-15SU (0 – 50 lpm), LFS-20SU (0 – 80 lpm)
LFC-1C-PC LFC-6CIO-PC	Single-Channel Converter: 24V DC supply input 6-Channel Converter: 24V DC supply input

EMC Directive 2014/30/EU – EMC Regulation SI 2016 No. 1091:

The following particular harmonized standards of the *EMC Directive* are tested and confirmed at a certified laboratory.

EN 61000-6-4	Generic Emission standard for industrial environments
EN 61000-6-2	Generic Immunity standard for industrial environments

Machinery Directive 2006/42/EU – Machinery Regulation SI 2008 Nr. 1597

The machinery directive essentially has been followed by a risk analysis, according mitigation actions and a user manual for safe operation. For the design and testing the following standards are used as a guideline:


EN ISO 15494	Plastic piping system for industrial applications
ISO12100	Safety for machinery – principles for risk assessments: used for system risk analysis.

Pressure Equipment Directive 2014/68/EU – Pressure Equipment Regulation SI 2016 No. 1105:

The pressure equipment directive is used as guide line according to *Category I Modul A* (internal quality control). For maximum pressure validation the following standard is followed:

EN ISO 15494	Plastic piping system for industrial applications
--------------	---

It has to be noted that the maximum pressure for the clamp-on sensor has been tested without tubing limitations. For the pressure rating of the clamp-on sensor with the mounted tubing, the specification of the used flexible tubing material has to be consulted in order to avoid bursting effects not covered by the clamp-on sensor specifications.

Start of CE-Marking:	2015-Jul-22	Place of issue:	Zurich, Switzerland
Last Update:	2022-Sep-27		
Levitronix GmbH  Dr. N. Barletta / CTO		Levitronix GmbH  Th. Eberle / Quality Manager	