



Polymun Scientific Immunbiologische Forschung GmbH

LNP Process Development and Analytics for mRNA Vaccines and Targeted LNPs

Andreas Wagner

**Levitronix Bioprocessing Conference
June 18th 2026**



Polymun Scientific Immunbiologische Forschung GmbH



■ A PRIVATE COMPANY

Developing and Manufacturing Biopharmaceuticals
and Liposomal Formulations for Human Application

- ▶ CEO: Dr. Dietmar Katinger
 - ▶ Founded: 1992
 - ▶ Employees: 98
-
- regularly inspected by the Austrian regulatory authority AGES on behalf of EMA, last inspection in December 2025
 - inspected by FDA in October 2013 / July 2023
 - numerous audits by clients (~10 per year)

Core Activities

- **Contract Development & Manufacturing of Biopharmaceuticals**
for human application with focus on mammalian cell culture, process development & GMP production
- **Contract Development & Manufacturing of LNPs and Liposomal Formulations**
LNP & liposomal formulation development for APIs and vaccine antigens & GMP production
- **Formulation of mRNA and oligonucleotides in liposomes/LNPs**
siRNA, saRNA, miRNA and mRNA formulated up to 300 g API input per batch
- **Liposomal adjuvants, liposomal vaccines**
liposomal formulation of MPLA as well as other TLR4 agonists in combination with other adjuvants like saponins, CpG,..
- **Covid-19 mRNA vaccine collaborations with:**
 - BioNTech/Pfizer
 - CureVac
 - Imperial College London
 - Arcturus Therapeutics
- **Research Reagents**
manufacturing and distribution of HIV antibodies and antigens
- **Own R&D Projects**
funded by revenues from contract development and contract manufacturing

Core Activities

- **Contract Development & Manufacturing of Biopharmaceuticals**
for human application with focus on mammalian cell culture, process development & GMP production
- **Contract Development & Manufacturing of LNPs and Liposomal Formulations**
LNP & liposomal formulation development for APIs and vaccine antigens & GMP production
- **Formulation of mRNA and oligonucleotides in liposomes/LNPs**
siRNA, saRNA, miRNA and mRNA formulated up to 300 g API input per batch
- **Liposomal adjuvants, liposomal vaccines**
liposomal formulation of MPLA as well as other TLR4 agonists in combination with other adjuvants like saponins, CpG,..
- **Covid-19 mRNA vaccine collaborations with:**
 - BioNTech/Pfizer
 - CureVac
 - Imperial College London
 - Arcturus Therapeutics
- **Research Reagents**
manufacturing and distribution of HIV antibodies and antigens
- **Own R&D Projects**

funded by revenues from contract development and contract manufacturing




If One Leading Coronavirus Vaccine Works, Thank This Tiny Firm in Rural Austria

Pfizer and BioNTech reliance on Polymun's product shows the fragility of the vaccine supply chain

At a Polymun laboratory in Klosterneuburg, Austria, the size distribution of lipid nanoparticles is measured. MARYLISE VIGNEAU FOR THE WALL STREET JOURNAL

mRNA Vaccines – Achievements within < 1 year

- Process set-up, optimization and scale-up
- Production of 10 different vaccines for tox studies
- Production of 5 different vaccines to initiate clinical trials
- Production of 2 different vaccines for phase 3 (> 40 000 subjects)
- Tech transfer to BioNTech/Pfizer network
- Analytical method validations
- (multi-center) process validation
- Regulatory support
- Production of 15 million doses to be used in EUA program in late 2020 / early 2021
- Continued bulk DP production in 2021




Production History




- 200+ GMP productions runs for oligo formulations
- 100+ GMP productions for
 - mRNA vaccines (incl. market product)
 - Gene editing formulations (mRNA/gRNA)
 - mRNA therapeutics (IV and pulmonary application)
- Intermediate batch volume > 500 L
- mRNA input amounts > 50 g / batch
- High reproducibility, High yields



Components of mRNA LNPs




Production of mRNA LNP Vaccines








Production of mRNA LNP Vaccines



- Process optimization
- Equipment optimization
- Scale-up; Scale-down



Mixing



**Development and optimization of a mixing process for the
preparation of lipid nanoparticle drug delivery systems for nucleic
acids using magnetic levitation pumps**

Master's thesis for obtaining the academic degree

Master of Science

in the study programme: **Clinical Research**

submitted by

Daniel Kumar, B.Sc.

Department for Health Sciences, Medicine and Research

at University for Continuing Education Krems


Advisor: Dr. Andreas Wagner

Vienna, July 2025




LNP processing – ongoing activities / optimization

- **Comparison:** piston pump batch with MagLev pump batches





- **Variations tested:**
 - Module size (0.5 mm vs. 1.0 mm)
 - Flow rates (low, standard, high)
 - Mixing angle (180° vs. 90°)




Analytics

- **CQAs**
 - Particle size & Size distribution (Pdl):
Dynamic Light Scattering (DLS)
 - Encapsulation efficiency:
RiboGreen assay



Baseline Comparison Piston vs. MagLev

- Flowrates of: 40/120 ml/min
- 180° mixing / 0.5 mm standard module






Piston pump (reference):

- Particle size: 66–71 nm
- Pdl: 0.11–0.15
- Encapsulation efficiency: 97-98%


MagLev pump:

- Particle size: 69-70 nm
- Pdl: 0.16–0.17
- EE: 95–96% (≈2% lower)





Larger module

- Flowrates of: 40/120 ml/min
- 180° mixing / 1 mm module




- Larger module (1.0 mm): stable, less backpressure
- Particle size : 70-71 nm
- EE: 96-97%
- Pdl: 0,14-0,15





Larger module – increased flow rates

- Flowrates of: 80/240 ml/min
- 180° mixing / 1 mm module




- Doubling flow rates possible only with 1.0 mm module - enables successful scale-up
- Particle size: 68-70 nm
- Pdl: 0,15-0,16
- EE: 96-97%






90° mixing

- Flowrates of: 40/120 ml/min
- 90° mixing / 0.5 mm standard module




- Particle size: 68–70 nm
- Pdl slightly improved 0.13–0.14
- EE stable at 96–97%



Summary


Key Observations:

- All batches yielded particle sizes around 70 nm
→ high reproducibility
- MagLev batches: comparable sizes (~68–72 nm), Pdl < 0.2
- Optimized MagLev settings (low flow / larger module / 90° angle): Pdl similar to piston pumps
- Only weak point: small module + high flow = back pressure → process not feasible
- Takeaway: MagLev delivers robust and reproducible sizes and Pdl values, with only minor differences compared to piston references.



Summary


ENCAPSULATION EFFICIENCY



Key Observations:

- Piston references: consistently > 97% EE across all runs
- MagLev batches: Direct comparison (MagLev 1–2): slightly lower, ~95–96% final product
- Optimized conditions (larger module, reduced flow, 90° mixing): improved back to 96–98%
- Takeaway: MagLev pumps achieve high EE comparable to piston pumps when operated under optimized conditions.

Emerging Trends in the LNP space – in-vivo CAR-Therapy



Polymun Scientific's Core Competences

Next Generation LNPs

**Modification of lipid
nanoparticle's surface by
conjugation of targeting
moieties**

*e.g.: aptamers, Fab-fragments,
peptides, antibodies*

leads to active targeting moiety
that binds to specific ligand

Develop
manu
biophar
using ma
tec

Fermentati Developme

- Batch
- Fed-bat
- Perfuse
- density
- systems

Fermentati GMP Produ

- Stirred
 - Bat
 - up
 - Per
 - up
- Dispos
 - 25

mun's innovative
osome / LNP
technology and
expertise


mun Technology

characteristics of our
egy

scalability
ptic process
ditions
le step process
nogeneous, uniform
cles
ellent batch to batch
sistency
procedure - stability


Why is targeting important?

Passive targeting (or non-targeting) LNP uptake mechanism



<https://doi.org/10.1016/j.omtm.2025.101436>




Active targeting LNP uptake mechanism “extrahepatic delivery”



https://ars.els-cdn.com/content/image/1-s2.0-S0168365924002438-gal_lrg.jpg

Lipid nanoparticles – Targeted delivery

LNP preparation





Targeting moieties:

- Antibodies, antibody fragment, peptides, etc.



Conjugation techniques


- Maleimide-thiol binding
- EDC/NHS coupling
- Biotin-avidin



The merge of Polymun's analytical expertise

How can we prove successful conjugation?

- Characterization using standard analytical methods for LNP release testing
- Analytical profiling of tLNPs needed to prove successful conjugation and efficiency (and stability!)
- Use of orthogonal methods




Standard Quality Control

- *mRNA content and EE% by Ribogreen Assay*
- *LNP size / size distribution (DLS)*
- *mRNA identity and integrity (Capillary electrophoreses)*
- *lipid identity and quantity (HPLC-CAD)*
- *mRNA content and ratio by IPRP HPLC*
- *pH*
- *Osmolality*
- *Bioburden testing*
- *Sterility testing*
- *Endotoxin testing*
- *Subvisible particles*
- *residual ethanol (GC)*




Last step after conjugation: Tangential Flow Filtration of tLNPs




Why?

- ▶ To remove unbound targeting moieties
- ▶ To remove organic residues from LNP preparation
- ▶ Exchange to desired buffer including cryoprotectants for frozen storage conditions



**Comparison of permeates during ultra-
/diafiltration of different conjugation strategies**



Lipid nanoparticles with targeting antibody moiety (2G12)



Legend:


1 – 2G12 control
2 – 2G12 thiol modification
(high excess)
3 – 2G12 thiol modification
(normal excess)

4 – IV coupled
5 – UDR low conc.
6 – 0.2 μ m low conc.
7 – UDR high conc.
8 – 0.2 μ m high conc.




- ✓ **Coupling of LNP and targeting moiety successful**
- Particle sizes around 65 nm, PDI ≤ 0.07
- %EE ≥ 95 %
- **Clear upwards shift** (antibody coupled to PEG-MAL) in SDS-PAGE and Western blot

SEC-HPLC with UV detection




- **From top to bottom:**
 - 1) not conjugated intermediate volume with antibody spike to show resolution between LNP and antibody (not-conjugated control)
 - 2) not conjugated intermediate volume (LNP control)
 - 3) not conjugated antibody (2G12) (Antibody control)

SEC-HPLC with UV detection



- Blue: conjugated
- Orange: not conjugated

SEC-HPLC with UV detection





Thank you

www.polymun.com

