

Levitronix® Standard Article List and Order Information

LEVIFLOW® LFS-SU Series - Single-Use Sensor, Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note		
LeviFlow Single-Use Sensors and Converter					
Single-Use Sensor LFS-03SU-Z (Calibration at 20°C, 0.8 l/min, Triclamp 3/8", Z-Shape) Single-Use Sensor LFS-03SU-Z-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor LFS-03SU-Z-SC1 (Calibration at 20°C, 0.8 l/min, Triclamp 3/8", Z-Shape) Single-Use Sensor LFS-03SU-Z-SC1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy)	100-30375 100-30399 100-30418 100-30419	LEVIFLOW LFS-SU Flowmeter	<i>Calibration for High Turn Down Ratio</i>		
Single-Use Sensor LFS-06SU-Z (Calibration at 20°C, 8 l/min, Triclamp 3/8", Z-Shape) Single-Use Sensor LFS-06SU-Z-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor LFS-06SU-Z-SC1 (Special Calibr. at 20°C, 8 l/min, Triclamp 3/8", Z-Shape) Single-Use Sensor LFS-06SU-Z-SC1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy)	100-30377 100-30400 100-30394 100-30406				
Single-Use Sensor LFS-10SU-Z (Water 20°C, 20 l/min, Triclamp 1/2", Z-Shape) Single-Use Sensor LFS-10SU-Z-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor LFS-10SU-Z-SC1 (Special Calibr. at 20°C, 20 l/min, Triclamp 1/2", Z-Shape) Single-Use Sensor LFS-10SU-Z-SC1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy)	100-30397 100-30405 100-30408 100-30416				
Single-Use Sensor LFS-15SU-Z (Water 20°C, 50 l/min, Triclamp 1", Z-Shape) Single-Use Sensor LFS-15SU-Z-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor LFS-15SU-Z-SC1 (Special Calibr. at 20°C, 50 l/min, Triclamp 1", Z-Shape) Single-Use Sensor LFS-15SU-Z-SC1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy)	100-30412 100-30311 100-30431 100-30432				
Single-Use Sensor LFS-20SU-Z (Water 20°C, 80 l/min, Triclamp 1", Z-Shape) Single-Use Sensor LFS-20SU-Z-G25 (Irradiated with Gamma Dosage ≥ 25 kGy) Single-Use Sensor LFS-20SU-Z-SC1 (Water 20°C, 80 l/min, Triclamp 1", Z-Shape) Single-Use Sensor LFS-20SU-Z-SC1-G25 (Irradiated with Gamma Dosage ≥ 25 kGy)	100-30379 100-30404 Coming Soon Coming Soon				
<i>Note: Adaptor/extension cables to be ordered as separate article.</i>					
LeviFlow Converter LFC-1C-PC (24VDC, Single Channel, 2 Digital Outputs, 1 Digital Input, 1 Analog Output, RS485)	100-30374			Spare Parts LEVIFLOW	Includes an EMI ferrite (article # 100-30353) for the flow sensor cable.
LeviFlow Converter Starter-Kit LFC-1C-PC-SK (Converter, AC/DC Desktop Power Supply, Cabling, Config.-Software)	100-91072			Spare Parts LEVIFLOW	
Accessories					
LeviFlow Interconnect Cable LFI-C.1-10 (1m, PVC) LeviFlow Interconnect Cable LFI-C.1-30 (3m, PVC) LeviFlow Interconnect Cable LFI-C.1-60 (6m, PVC)	190-10307 190-10308 190-10309			Accessories LEVIFLOW	
LeviFlow Adaptor Extension Cable LFE-C.2-10 (1m, PVC, Wall-Mountable Connector) LeviFlow Adaptor Extension Cable LFE-C.2-30 (3m, PVC, Wall-Mountable Connector) LeviFlow Adaptor Extension Cable LFE-C.2-60 (6m, PVC, Wall-Mountable Connector)	190-10310 190-10311 190-10312			Accessories LEVIFLOW	
LeviFlow Mounting Kit LMK-1.1 (Anodized Alu., for LFS-03SU and LFS-06SU) LeviFlow Mounting Kit LMK-2.1 (Anodized Alu., for LFS-10SU) LeviFlow Mounting Kit LMK-3.1 (Anodized Alu., for LFS-15SU) LeviFlow Mounting Kit LMK-4.1 (Anodized Alu., for LFS-20SU)	100-91179 100-91180 100-91181 100-91182			Accessories LEVIFLOW	
LeviFlow YN-485I-TR (RS485/USB Adaptor, Cable with wiring connector, Termination, Isolated)	100-30392			Accessories IPS	
Spare Parts					
NA	NA	NA			
NA	NA	NA			