

Levitronix® Standard Article List and Order Information
Maglev Pumpsystem PuraLev® i100, Accessories and Spare Parts

Article Name (and Description)	Part No.	Article Category	Note
Standard Driver System Configurations			
Pump Driver System PLD-i100SU.1 (OEM Model) Integrated Pump Driver IPD-100.1-50-02 (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91054	PuraLev SU	
Pump Driver System PLD-i100SU.5 (OEM Model, ATEX/IECEX) Integrated Pump Driver IPD-100.10-50-02 (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91511	PuraLev SU	
Pump Driver System PLD-i100SU.2 (EasyConnect Model) Integrated Pump Driver IPD-100.4-02 (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91095	PuraLev SU	
Pump Driver System PLD-i100SU.3 (Stand-Alone Model) Integrated Pump Driver IPD-100.6-02 (24VDC, Epoxy Housing Coating, Display with Buttons) Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91055	PuraLev SU	
Single-Use Pump Heads			
Single-Use Pumphead DCP-200.2 Welded and Molded PP Housing, Triclip 1/2", Overmolded PP Impeller FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile	100-90734	PuraLev SU	
Single-Use Pump Head Gamma Irradiated DCP-200.2-G25 (Applied Gamma Dosage \geq 25 kGy)	100-91078	PuraLev SU	
Single-Use Pumphead DCP-200.3 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller FDA, USP-VI, BSE/TSE and Animal Free, Gamma Sterilizable, Non-Sterile	100-90792	PuraLev SU	
Single-Use Pump Head Gamma Irradiated DCP-200.3-G25 (Applied Gamma Dosage \geq 25 kGy)	100-91122	PuraLev SU	
Single-Use Pumphead DCP-200.3-SF1-G25 Welded and Molded PP Housing, Barb 1/2", Overmolded PP Impeller, FDA, USP-VI, BSE/TSE and Animal Free, Sterile Fittings, Applied Gamma Dosage \geq 25 kGy)	100-91233	PuraLev SU	
Accessories			
Power Supply TSP 180-124 (Traco, 24V / 180W, DIN Rail Mounting, 85-132 and 187-264VAC autoselec)	100-40018	Accessories Pump Systems	
Desktop AC/DC Power Supply GST160A24-R7B IC915 (24VDC/160W, 85-264VAC, Intercontec Conn.)	100-40020	Accessories IPS	AC mains cable AMC-1.x for various countries to be ordered as separate article.
AC Mains Cable AMC-1.1 (2m, US, Canada) AC Mains Cable AMC-1.2 (2.5m, Germany, Denmark, Norway, Finland, Belgium, Netherl., Sweden, Austria) AC Mains Cable AMC-1.3 (2.5m, Japan) AC Mains Cable AMC-1.4 (2.5m, Switzerland) AC Mains Cable AMC-1.5 (2.5m, United Kingdom)	190-10331 190-10332 190-10333 190-10334 190-10335	Accessories IPS	AC cable for desktop power supply.
USB-RS485 YN-485I-TR (USB to RS485 Adaptor, Isolated, Cable with Termination Resistance)	100-30392	Accessories IPS	
IPS Cable Power 2 Wires ICP-2.1-50 (5m PVC, Circular Driver Connector to Open Wires)	190-10370	Accessories IPS	For power supplies with open wire connection.
IPS Cable Signal 6 Wires ICS-1.1-02 (0.2m PVC, Circular Driver Connector to Circular Connector) IPS Cable Signal 6 Wires ICS-1.1-10 (1m) IPS Cable Signal 6 Wires ICS-1.1-30 (3m)	190-10381 190-10344 190-10345	Accessories IPS	For connecting fieldbus of EasyConnect driver arrays.
IPS Cable Signal 6 Wires ICS-1.2-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10346	Accessories IPS	For connecting devices to PLC of Stand-Alone driver or fieldbus of Connectivity driver.
IPS Cable Signal 12 Wires ICS-2.1-50 (5m PVC, Circular Driver Connector to Screw Type Connector)	190-10347	Accessories IPS	For connecting devices to PLC of EasyConnect driver.
Fieldbus Termination Connector FTC-1.1 (PPS with FPM Sealing)	190-10348	Accessories IPS	Termination of fieldbus of EasyConnect drivers.
ATEX Cable Sealing System ACS-A.1	100-90292	Accessories IPS	
Spare Parts			
Pump Head Socket PHS-i100.1 (Anodized Aluminum, Bayonet Type)	100-91053	Spare Parts IPS	
Integrated Pump Driver IPD-100.1-50-02 (OEM Model) (24VDC, 5m PVC Cable Open Wires, Epoxy Housing Coating)	100-10070	Spare Parts IPS	
Integrated Pump Driver IPD-100.10-50-02 (OEM Model, ATEX/IECEX) (24VDC, 5m FEP Cable Open Wires, ETFE Housing Coating)	100-10201	Spare Parts IPS	
Integrated Pump Driver IPD-100.4-02 (EasyConnect Model) (24VDC, Epoxy Housing Coating, PLC and 2x Fieldbus)	100-10111	Spare Parts IPS	
Integrated Pump Driver IPD-100.6-02 (Stand-Alone Model) (24VDC, Epoxy Coating, Display with Buttons, PLC)	100-10106	Spare Parts IPS	